



How-To-Manual

Psycho-acoustical and physiological test procedures for the assessment of hearing impairment
Scientific background, methodology, practical use

Senti & Sentiero

Manufacturer

PATH MEDICAL GmbH

Landsberger Str. 65 D-82110 Germering Germany

Email info@pathme.de
Telephone +49 89 800 765 02
Fax +49 89 800 765 03

Manual Information

Article number: 100893-EN Release date: 2017-09-08

Revision: 05

Valid from: Firmware Rev. 2.2¹, Mira PC Software Rev. 1.7

All mentioned items, products, brands and trademarks are registered or owned by the mentioned companies.

All information, illustrations, and specifications provided within this manual are based on the latest product information available at the time of publication. PATH MEDICAL reserves the right to make changes at any time without notice.

The latest revision of the How-To-Manual is available online at www.pathme.de/support/download.

Errors and omissions excepted.

Copyright Notice

No part of this manual may be reproduced, translated, stored, or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise, without the prior written consent of PATH MEDICAL GmbH.

The copyright of mentioned literature remains with the authors. Please refer to the literature reference list.

Copyright © 2017 PATH MEDICAL GmbH

0124

_

¹ Please note that some planned firmware features are already described in this document.

Table of Contents

1	Introduct	tion	5
2	Overview	v	7
3	Methodo	ology, Practical Use, Scientific Background	13
3	.1 Psyd	cho-Acoustical Test Procedures	14
	3.1.1	Pure-Tone Audiometry (PTA)	15
	3.1.2	Multiple-Choice Auditory Graphical Interactive Check (MAGIC)	20
	3.1.3	Overview: Speech Intelligibility Tests	25
	3.1.4	Speech Understanding in Noise (SUN)	26
	3.1.5	Mainzer Audiometric Test for Children (MATCH)	30
	3.1.6	Universal Speech Test (UST)	35
	3.1.7	Speech CD Player	40
	3.1.8	Live Speech	42
	3.1.9	Munich Auditory Screening Test for Processing Disorders (MAUS)	44
	3.1.10	Bochum Auditory Speech Discrimination Test (BASD)	48
	3.1.11	Hearing Handicap Inventory for the Elderly (HHIE-S)	53
3	.2 Phy	siological Test Procedures	55
	3.2.1	Overview: Otoacoustic Emissions	56
	3.2.2	Transient Evoked Otoacoustic Emissions (TEOAEs)	63
	3.2.3	Distortion Product Otoacoustic Emissions (DPOAEs)	69
	3.2.4	Overview: Auditory Evoked Potentials	84
	3.2.5	Auditory Brainstem Reponses (ABRs)	88
	3.2.6	Auditory Steady-State Responses (ASSR)	106
	3.2.7	Electrocochleography (ECochG)	114
	3.2.8	Overview: Immittance Audiometry	117
	3.2.9	Tympanometry	121
	3.2.10	Acoustic Reflex Test	128
	2 2 11	Fustachian Tuhe Function Tests	12/

1 Introduction

Senti and Sentiero devices are PC-independent, portable medical instruments, designed for easy handling and detailed assessment of hearing impairment at different stages of the auditory pathway in infants, children, and adults. Various transducers for eliciting stimuli and electrodes for measuring auditory potentials from the scalp can be connected. A touch-screen display allows for patient data management, parameter setting, control of measurement procedures, and visualization of results. Data can be transferred to a PC for further analysis.

All devices are modular, so that subsets of available test methods can be combined. Some of the methods mentioned in this manual may not be included on your device. The availability of test modules depends on your device type (Senti, Senti Desktop, Sentiero, Sentiero Advanced, and Sentiero Desktop) and on your device license. Please contact your distributor if you would like to upgrade your license to include more test modules. Screenshots are provided as examples and may deviate from the representation on your device.

The purpose of this How-to Manual is to get familiar with all the methods available on Senti and Sentiero devices. In section 2: Overview all methods are quickly described together with typical clinical applications and the device types for which the respective module is available. Information about general aspects of the method including the scientific background and its practical use are explained for all test methods in section 3: Methodology, Practical Use, Scientific Background.

The *How-To-Manual* provides information on Senti and Sentiero devices seen from a daily use perspective. It does not replace the *Senti and Sentiero User Manual*, which is provided together with the device and which describes the main device function in a condensed version. For general information about device handling and functions, intended use, service and maintenance, cleaning, accessories, warranty, notes on safety, technical specifications, and electromagnetic compatibility, please refer to the *Senti and Sentiero User Manual*. For information on technical data and available test module parameters please refer to the *Senti and Sentiero Technical Specifications*. For additional information on the device and its test modules please refer to the **PATH MEDICAL website's** *Learning* section, which includes scientific papers, device tutorials, quick guides, and device videos.

The *How-To-Manual* is provided "as-is". PATH MEDICAL has taken care that the contents and descriptions are correct. The content was gathered and proof-read under scientific contribution of professionals and long term experts in the field. Selected literature is provided to enable the reader to read more on the scientific background of the methods implemented on Senti and Sentiero. PATH MEDICAL cannot take any liability arising from (mis-)interpretation or (mis-)application of this *How-To-Manual* and the provided references therein.

Please note that this *How-To-Manual* is meant to assist audiologically qualified personnel and does not replace adequate training and qualification.

2 Overview

Table 1 describes briefly all test modules available on Senti and Sentiero devices (including Senti [model SIH100097], Senti Desktop [models SID100419, SID100433], Sentiero [models SOH100098, SOH100360], and Sentiero Desktop [model SOD100497]) together with their typical clinical application. Please refer to the Senti and Sentiero User Manual for the device's intended use and contraindications.

Test module	Short description / Clinical application	Devices			
	PURE-TONE AUDIOMETRY				
PTA Pure-tone audiometry	Pure-tone audiometry according to IEC 60645-1 including screening and threshold mode, air and bone conduction stimulation, contralateral masking noise, high-frequency audiometry up to 16 kHz and automatic threshold detection procedures.				
	Frequency-specific determination of hearing status in cooperative subjects for audiological screening and diagnostics and occupational medicine.				
MAGIC Multiple-choice graphic interactive check	Image-based and self-controlled pure-tone audiometry including air and bone conduction stimulation. - MAGIC Screen: screening mode with fixed user-defined stimulus levels - MAGIC Audio: threshold mode with automatically controlled stimulus levels.	All			
	Frequency-specific determination of hearing status in cooperative subjects (especially in pre-school and school children) for audiological screening and diagnostics.				
SUN Speech Understanding In Noise	Self-controlled speech audiometry with logatoms and ipsilateral noise. - SUN Predefined: screening mode with fixed user-defined speech level - SUN Adaptive: threshold mode with automatically controlled speech levels	All			
	Determination of speech intelligibility in noise in cooperative subjects for audiological screening (especially in the elderly) (SUN Fixed) and diagnostics (SUN Adaptive).				
MATCH Mainzer Audiometric Test for Children	Image-based and self-controlled speech audiometry with realistic words and optional ipsilateral noise. - MATCH Fixed: fixed user-defined speech level - MATCH Adaptive: threshold mode with automatically controlled speech levels	All			
	Determination of speech intelligibility in quiet and in noise in cooperative subjects for audiological diagnostics. Words and images are optimized for two-year-old children.				

Test module	Short description / Clinical application	Devices
UST Universal Speech Test e.g. Freiburger (German), Dr. Tato (Spanish), NU-6 (English), CID W-22	Speech audiometry with words (e.g. monosyllabic words, numbers) and optional ipsilateral or contralateral noise. The test module is available for different word lists. Note: there are specific workflow adaptations to fit the typical workflow of a speech test.	All ¹
(English), PBK-50 (English), Maryland CNC (English), Mots français CAD (French)	Determination of speech intelligibility in quiet and in noise in cooperative subjects for audiological diagnostics. Intended age group depends on the speech material of the particular test.	
Speech CD Player e.g. NVA (Dutch)	Speech audiometry with words (e.g. monosyllabic words, numbers) played continuously as on a CD player and with optional ipsilateral or contralateral noise. Note: there are specific workflow adaptations to fit the typical workflow of a speech test.	All ¹
	Determination of speech intelligibility in quiet and in noise in cooperative subjects for audiological diagnostics. Intended age group depends on the speech material of the particular test.	
Live Speech	Speech audiometry with live speech and multiple test modes including speech detection threshold, speech recognition threshold, and word recognition.	All
	Determination of speech detection thresholds and speech/word recognition thresholds in quiet.	
	OTHER SUBJECTIVE AUDIOMETRY TESTS	
MAUS Munich Auditory Screening Test for Processing Disorders	Speech-based test with subtests regarding sequence of syllables, words in noise, phoneme differentiation and identification.	All ¹
DISOTUETS	Screening for auditory processing disorders (ADP) in cooperative subjects (especially in children from 6 to 11 years).	
BASD Bochum Auditory Speech Discrimination Test	Auditory discrimination test battery with subtests regarding discrimination of consonants and high resolution changes in frequency, level, time, and amplitude modulation and with monaural, binaural, dichotic, and interaural stimulus presentation.	All
	Diagnostics of auditory processing and speech development disorders in cooperative subjects.	
HHIE-S Hearing Handicap Inventory for the Elderly-	Questionnaire concerning activity limitations, participation restrictions, associated feelings and emotions of daily life.	All
Screening	Determination of subjective hearing status in cooperative adult subjects for audiological screening.	

Test module	Short description / Clinical application	Devices			
	OTOACOUSTIC EMISSIONS				
TEOAE Transient Evoked Otoacoustic Emissions	Objective test for checking outer hair cell functionality within 0.7 to 4 kHz (either in entire frequency range or in five frequency bands) using a short-term broadband stimulus and statistical response detection. - TEOAE Quick: fixed stop criterion - TEOAE Diagnostic: user-defined stop criterion in different frequency bands	Sentiero [SOH100098, SOH100360] Sentiero Desktop [SOD100497]			
	Determination of outer hair cell status in subjects of all age groups (especially in newborns and infants) for screening and audiological diagnostics.				
DPOAE Distortion Product Otoacoustic Emissions	Objective test for checking outer hair cell functionality at specific frequencies using two sine stimuli overlapping on the basilar membrane and statistical response detection. - DPOAE Quick: one fixed user-defined stimulus level - DPOAE Diagnostic: multiple fixed user-defined stimulus levels (Sentiero Desktop: optionally with static pressure offset) - DPOAE Threshold: threshold estimation with automatically controlled stimulus levels Additional optional features: - DPHIRES: user-defined frequency selection (start, stop, step size): allows for assessing DPOAE fine structure; multichannel measurement: allows for simultaneous test at multiple frequencies - FMDPOAE: frequency-modulated DPOAE, allows for reducing the impact of the 2 nd DPOAE source	Sentiero [SOH100098, SOH100360] Sentiero Desktop [SOD100497]			
	Determination of outer hair cell status in subjects of all age groups for audiological purposes: e.g., follow-up after newborn hearing screening, confirmation of cochlear hearing loss, topological hearing diagnostics, detection of early stage cochlear impairment during noise overexposure or ototoxic drug administration; frequency-specific determination of hearing thresholds (<i>DPOAE Threshold</i>).				
	AUDITORY EVOKED POTENTIALS				
ABR Auditory Brainstem Responses	Objective test for checking cochlear and neural sound processing using various short-term stimuli. - Quick ABR: one fixed user-defined stimulus level, statistical wave V evaluation - ABR: multiple fixed user-defined stimulus levels with multiple stimulus types (click, chirp), air and bone conduction stimulation, contralateral masking, statistical wave V detection Additional optional features: - ABR-FS: low/mid/high chirp, tone bursts - ABR-BIN: binaural measurement	Sentiero [SOH100360]			

Test module	Short description / Clinical application	Devices
	Determination of hearing status along the auditory pathway up to the brainstem in subjects of all age groups for audiological purposes: e.g. follow-up diagnostics after newborn hearing screening, neurological diagnostics, topological hearing diagnostics, determination of hearing thresholds.	
E-ABR (in preparation) Electrically Evoked Auditory Brainstem Responses	Objective test for checking cochlear and neural sound processing with statistical response detection by means of recording electrically evoked potentials externally triggered by a cochlear implant test system (provided by the cochlear implant manufacturer).	
	Determination of hearing status along the auditory pathway up to the brainstem in subjects of all age groups with a cochlear implant for audiological purposes including cochlear implant fitting.	
ECochG (in preparation) Electrocochleography	Objective test for checking cochlear and neural sound processing using short-term click stimuli and statistical response detection.	Sentiero [SOH100360]
	Determination of hearing status along the auditory pathway including summating potential (SP), action potential (AP) and cochlear microphonic (CM) in subjects of all age groups for audiological purposes.	
ASSR Auditory Steady-State Responses	Objective test for checking cochlear and neural sound processing using frequency-specific chirp stimuli (cochlear travelling wave compensated), 40/80 Hz and mixed stimulus rates, multichannel testing, contralateral masking, and statistical response detection. - ASSR Fixed: fixed user-defined stimulus levels - ASSR Threshold: threshold estimation with automatically controlled stimulus levels	Sentiero [SOH100360]
	Determination of hearing status along the auditory pathway up to the brainstem (80 Hz) and subcortical regions (40 Hz) in subjects of all age groups for audiological purposes: e.g. follow-up diagnostics after newborn hearing screening, neurological diagnostics, topological hearing diagnostics, frequency-specific determination of hearing thresholds.	
	IMMITTANCE AUDIOMETRY	
Tympanometry	Objective test for measuring middle ear impedance dependent on static pressure (+/-) in the outer ear canal.	Sentiero Desktop
	Determination of middle ear status in subjects of all age groups for audiological purposes: e.g. detection of middle ear pathologies (e.g. otitis media, middle ear effusion, perforation of eardrum, dysfunction of auditory ossicles, Eustachian tube dysfunction, otosclerosis, tympanosclerosis) or blockage of the ear canal.	[SOD100497]

Test module	Short description / Clinical application	Devices	
Acoustic Reflex Threshold	acoustic reflex elicited by tones with different levels and frequencies or noise. - Automatic mode (after performing tympanometry) - Manual mode		
	Determination of acoustic reflex thresholds in subjects of all age groups for audiological purposes: e.g. topological diagnostics and detection of middle ear pathologies (e.g. otosclerosis) or neural pathologies (e.g. facial nerve dysfunction).		
Acoustic Reflex Decay	Objective test for measuring the decay of ipsilateral and contralateral acoustic reflexes elicited by tones with different levels and frequencies or noise.	Sentiero Desktop [SOD100497]	
	Determination of acoustic reflex decay in subjects of all age groups for audiological purposes, e.g. detection of retrocochlear pathologies.		
ETF Eustachian Tube Function Tests	Objective impedance audiometry tests with multiple subtests: non-perforated eardrum, perforated eardrum, patulous Eustachian tube.	Sentiero Desktop [SOD100497]	
	Determination of Eustachian tube function at different test conditions (non-perforated eardrum, perforated eardrum, patulous Eustachian tube)		
only available for Senti [model SIH100097] and Sentiero [model SIH100098] devices with extended memory (i.e. PCB rev. ≥ 67)			

Table 1: Overview of test methods available on Senti and Sentiero devices

3 Methodology, Practical Use, Scientific Background

A description of the methodology and a practical guide of the methods used in Senti and Sentiero devices are presented in the following. This includes conventional pure-tone audiometry (PTA or Audio), image-based pure-tone audiometry (MAGIC), speech tests (e.g. SUN, MATCH), and other psycho-acoustical test methods (e.g. MAUS, BASD). Physiological test methods comprise otoacoustic emissions (OAE), auditory brainstem responses (ABR), electrocochleography (ECochG, in preparation), auditory steady-state responses (ASSR), and middle ear tests (e.g. tympanometry, acoustic reflex, Eustachian tube function tests). The highly sophisticated mechanisms of OAE, ABR, and ASSR generation, recording, reproducibility, stimulus parameter settings, and clinical applications are described in detail. Especially, the relationship between these measures and behavioral puretone thresholds is discussed. Additionally, clinical case examples are given. A list of references to further reading is also included for each module.

Please note that detailed technical information as e.g. available transducers, maximum levels, and other technical details are presented for each module in a separate *Technical Specification* document available for download from www.pathme.de/support/. For more information about available test module parameters please refer to the device online help.

In general, it is recommended to conduct all tests in an acoustically shielded booth or in a separate quiet room (especially when presenting low-level stimuli). When using ear probes or insert earphones, the transducer must be placed with a tight fit without any leakage between ear tip and outer ear canal. Proper occlusion of the ear canal by the ear tip diminishes the influence of external sounds.

For reliable AEP measurements an electromagnetically shielded booth should be used. A shielded electrode cable also reduces the impact of electromagnetic noise. Moreover, it is recommended to keep the electrode leads close to the subject's body and to carefully separate the electrode cable from the sound transducer cable.

Ear probe or insert earphone channels must not be blocked e.g. with cerumen. Ear probes or insert earphones should not be used in cases of external otitis (outer ear canal infection) or in cases where the outer ear canal is occluded with cerumen. Prior to any test, it is recommended to perform an otoscopy in order to determine if there is any blockage in the ear canal or any other visible middle ear disorder (e.g. perforated ear drum). Also a deformation of the ear canal may deteriorate sound transmission and may therefore deteriorate test performance.

It is recommended to perform as many audiological tests as possible in order to get a complete overview on the hearing status along the auditory pathway, especially if there are doubts about the true hearing status or if there are any inconsistencies or discrepancies. Qualified personnel needs to decide which tests are appropriate given the age, cooperativeness, and clinical history/anamnesis of the tested subject. Please note that interpretative hints and recommended parameter settings are provided in this *How-To-Manual* on an informative basis only and do not give any hint towards diagnosis. It is up to the qualified examiner to analyze results and to set parameters according to individual circumstances.

3.1 Psycho-Acoustical Test Procedures

Psycho-acoustical test procedures include pure-tone audiometry, speech audiometry and other behavioral tests which require patient feedback.

For all psycho-acoustical tests, the ability, willingness, and concentration of the subject to be tested are essential. Before performing any psycho-acoustical test procedure, the subject must be instructed about the task. The examiner must ensure that the subject has understood the task and is capable of performing the task. Always place the transducer on the subject's ears *after* instruction so that the subject's hearing ability is not reduced due to the transducer's sound attenuation.

Low ambient noise conditions are also essential for an appropriate test performance since hearing performance may be artificially reduced by high ambient noise.

In general, proper transducer placement and calibration is important for a reliable outcome. Incorrect placement of a transducer or usage of an incorrectly calibrated transducer may adulterate the result. A suitable transducer should be connected and properly fitted. Variation in earphone position up-down or front-back and leakage between earphone and ear may yield inappropriate results as e.g. falsely poor hearing thresholds, particularly at high frequencies.

If not specified otherwise all psycho-acoustical tests can be performed in subjects of all ages as long as appropriate participation and cooperation can be expected. This means that these tests are not expected to be suitable e.g. for very small children, mentally handicapped subjects or other non-cooperative subjects. For these subjects it is recommended to apply physiological test methods.

3.1.1 Pure-Tone Audiometry (PTA)

METHODOLOGY

PTA is the most commonly used psycho-acoustic method for testing a subject's ability to hear various sound frequencies. PTA is the key hearing test for assessing air conduction (AC) and bone conduction (BC) hearing thresholds enabling the determination of degree and type (conductive or sensorineural) of the hearing loss. Hearing thresholds typically indicate the softest sound audible to an individual. AC stimuli can be delivered via headphone, insert earphones or loudspeakers. BC stimuli can be delivered via an electro-mechanical bone conductor which is positioned on the forehead or mastoid of the subject to be tested. If there is better hearing for bone conduction than for air conduction (airbone gap) a conductive hearing loss is likely. If both air and bone conduction thresholds are out of the normal range a sensorineural (cochlear or retro-cochlear impairment) is likely. PTA can be used to assess the hearing status for tonal stimuli in cooperative subjects who are able to perform the task, e.g., in clinics (conventional audiology, bed-side use, ambulances) and occupational medicine. In order to differentiate a detected sensorineural hearing loss, supra-threshold measurements (e.g. Fowler, Carhart, categorical loudness scaling) may be performed. Otoacoustic emissions (OAE) and auditory evoked potentials (AEP) may be used for further topological diagnostics. Please note that for subjects who are not expected to be able to perform the test or are not able to react to the sound in a proper way to get frequency-specific information of the hearing loss, DPOAE or ASSR Threshold measurement is advised.

The test procedure demands an examiner who controls the test and who evaluates the response of the tested subject. Hearing is commonly tested at octave frequencies varying from low (250 Hz) to high frequencies (8 kHz). A typical test procedure (butterfly) starts at 1 kHz, moves up to higher frequencies, goes back to 1 kHz (re-check) and moves down to lower frequencies. It should be emphasized that the audiometric frequency range is just a part of the entire human auditory range, which extends between 20 Hz and 20 kHz. The core method of PTA is to present a pure-tone close to threshold, i.e. keep dropping the sound pressure level until the person stops responding (raising a hand or pressing a button) and then increasing the sound pressure level until the person starts responding again (see ISO 8253-1).

During some conditions subjects must be tested with masking noise on the non-test ear in order to prevent that sound going from the stimulated side over to the non-stimulated side can be heard on the non-stimulated side (cross hearing). The threshold of the test ear is measured at the same time as presenting the masking noise to the non-test ear. Thresholds obtained with masking provide an accurate representation of the true hearing threshold of the test ear. The interaural attenuation can go down to 0 dB for bone conductors because the bones of the skull are very efficient at transmitting sound. Therefore, it is recommended to measure bone conduction thresholds always with contralateral masking noise. For headphones, interaural attenuation amounts to about 40 dB. Masking should be used if the difference in air conduction in one ear and bone conduction in the other ear is 40 dB or greater. For insert earphones, interaural attenuation is in the range of 55 dB so that the use of insert earphones reduces the need for masking. A masking dilemma occurs when masking from the non-test ear crosses over to the test ear and affects threshold testing for the test ear. In this case, a reliable masked threshold cannot be obtained. This phenomenon generally occurs in the presence of a substantial conductive hearing loss component.

Automatic PTA assesses hearing thresholds following Békésy or Hughson-Westlake procedures. In contrast to the method of adjustment, for which the subject can vary the stimulus level until it is just audible, the subject can only control the direction in which the stimulus level varies. The subject increases and decreases the stimulus level around hearing threshold depending on its audibility (method of tracking). The available measuring procedures correspond to the algorithms suggested by ISO 8253-1.

Screening PTA examines whether the hearing threshold levels at different frequencies are better, equal to, or worse than the specified screening level. The stimulus is delivered to the ear at specific frequencies and levels that have been chosen for screening. The subject passes the screening test according to predefined criteria, e.g., a *pass* occurs if all stimuli are heard in each ear. Screening at certain selected levels and frequencies is typically faster but provides less information than measuring hearing thresholds using common PTA procedures.

Different PTA modules are available:

- **PTA 4** for assessing AC pure-tone thresholds at frequencies from 125 Hz to 6 kHz with stimulus levels from -10 to 70 dB HL in compliance with IEC 60645-1 class 4.
- **PTA 4 Advanced** for assessing AC pure-tone thresholds at frequencies from 125 Hz to 8 kHz with stimulus levels from -10 up to 110 dB HL in compliance with IEC 60645-1 class 4 with extended frequency and level range.
- **PTA 3** for assessing AC and BC pure-tone thresholds in compliance with IEC 60645-1 class 3. Frequencies and levels are the same as for PTA 4 Advanced. Due to the limited capacity of the electro-mechanical bone conductor, BC stimulus levels are lower. The module also provides a contralateral masking option, presentation of different stimulus types (pure-tone, pulsed pure-tone, warble tone), and automatic threshold detection procedures (Békésy, Hughson-Westlake).
- **PTA-HF** for assessing pure-tone thresholds at high frequencies from 9 to 16 kHz. PTA-HF is only available with specific headphones, e.g. Sennheiser HDA 300.

PRACTICAL USE

Select *Audio* from the module selection screen. Select the preset that you would like to use. If necessary, change the parameters (e.g. audiometry protocol). If *Expert* mode is selected, the user interface appears as shown in *Figure 1*. Press the *settings* button ③ (PTA3 only) in order to select the transducer and its placement (AC, BC mastoid, BC forehead), select the stimulus type (pure-tone, pulsed pure-tone, warble tone), switch to uncomfortable loudness (UCL) mode, activate automatic testing, or activate the bilateral audiogram view (plot data from both ears in one audiogram graph). Please note that you can also switch transducers by pressing the *AC* button ③ which toggles from AC to BC (mastoid) to BC-F (forehead) if an AC and a BC transducer are connected.

For some subjects (e.g. children) a pulsed pure-tone or a warble tone may be more interesting than a continuous pure-tone and hence may be preferred. For the Békésy procedure, pulsed stimuli are used. For subjects with a tonal tinnitus, a pulsed pure-tone or warble tone may distinguish the stimulus from the tinnitus more efficiently. Enable masking mode if required by pressing the *masking* button (12) (see *METHODOLOGY*). Masking mode can be left by pressing button (14).

Before the transducer is fitted and the test is started, the subject must be instructed about the task. The subject shall indicate if a tone is heard or not, e.g. by raising a hand or pressing the patient response button if a tone is heard or by putting down the hand or releasing the patient response button if *no* tone is heard. For pulsed tones, please make sure that the subject understands to respond to the pulsed tone once (i.e. indicate *heard* as long as the pulsed tone is audible) and not for each single stimulus pulse. The test is *not* intended to be self-controlled by the subject. Manual tests must be fully controlled by the examiner. Automatic tests do not need level/frequency control by the examiner, but should be supervised by qualified personnel nevertheless.

Select the test ear by pressing the ear button (1). The currently active ear is shown. If loudspeakers are connected or binaural measurement is enabled on the device, also binaural stimulation can be selected. Perform either a manual test or an automatic test at the selected frequencies. For further information about a typical workflow of pure-tone audiometry see e.g. ISO 8253-1 and ASHA (2005). During manual tests, stimulus frequency (8), level and/or masking level (7) (in masking mode masking level buttons are highlighted with a gray area) can be selected. The selected levels can be read from the crosshairs in the audiogram (red: right ear, blue: left ear, black: binaural, green: masking) and from the values between the selector elements (in Figure 1: stimulus level: 40 dB HL, making level: -10 dB HL). For PTA-HF, the audiogram shifts to higher frequencies when moving the crosshairs beyond 8 kHz. In masking mode, either the stimulus level or the masking level or both levels locked together can be controlled. The locked mode can be activated by pressing the lock button (5). The stimulus level buttons can be configured in the device Audio Preferences (up arrow: decrease level, down arrow: increase level or vice versa). The stimulus can be presented as long as the *loudspeaker* button 9 is pressed. If the stimulus is played, the orange status indicator 1 is on. If the masking noise is played, the green status indicator (2) is on. If the patient response button is pressed, the large green status indicator (3) is on. If a threshold is determined, set the symbol at the position of the red or blue crosshairs (6) by pressing the audiogram (5). The symbol can be toggled from *heard* to *not heard* to deleting the symbol. The test can be finished with the *stop* button ①.




Figure 1: Audio Expert user interface (left: standard mode; right: masking mode)

Hearing thresholds are plotted across test frequencies in an audiogram form. Right ear symbols are plotted in red, left ear symbols in blue, and binaural symbols in black. The used symbols refer to ISO 8253-1, Table 1. BC symbols can be plotted as proposed by ISO 8253-1 or mirrored (usual representation e.g. in Germany). The BC symbol and no response symbol representation can be configured in *Audio Preferences*.

The audiogram is divided into three colored zones: green (normal hearing \leq 20 dB HL), orange (mild hearing loss >20 to 40 dB HL), and gray (>40 dB HL: moderate to profound hearing loss).

For the interpretation of the result, the following hints may be considered (see Figure 2):

- Conductive hearing loss: Audiogram shows normal bone conduction thresholds, air conduction thresholds are poorer than normal by at least 10 dB. This air-bone gap occurs due to the damping of the air conduction stimuli on the way through the outer and middle ear. A conductive hearing loss typically occurs due to outer or middle ear problems including abnormalities of the tympanic membrane, occlusion of the auditory canal (e.g. by cerumen), middle ear infection or fluid (e.g. otitis media), perforation of the tympanic membrane, ossicular defects.
- Sensorineural hearing loss: Audiogram shows air and bone conduction thresholds that are both within a range of about 10 dB. Thresholds are higher than 20 dB HL. A sensorineural hearing loss typically occurs due to lesions of sensory cells, auditory nerve, or central auditory pathways (presbyacusis, noise-induced hearing loss, Menière's disease, vestibular schwannoma).
- **Mixed hearing loss**: Mixture of conductive and sensorineural hearing loss.




Figure 2: Audio result examples (left: conductive hearing loss; middle: sensorineural hearing loss; right: mixed hearing loss with conductive and sensorineural components)

Screen, Diagnostic and Auto Audio protocols provide adapted user interfaces as shown in Figure 3. The user interfaces of these protocols are in many parts similar to the Expert user interface. In the following the differences are explained. The Screen user interface provides a result table ① instead of an audiogram graph. The table contains either the tested level (heard) or NR = no response (not heard). The currently active frequency/level combination is highlighted in the result table. The Auto Audio user interface provides a play button ③ instead of the loudspeaker button, which allows starting the automatic threshold detection procedure. The patient's response, which is read from the patient response button, influences the direction of the stimulus level. If a threshold is detected by the automatic algorithm, it is automatically set in the audiogram graph. In all protocols, there is no settings and no masking button. Instead, in the footer you have direct access for changing the stimulus type (pure-tone, pulsed pure-tone, warble tone) by pressing the stimulus type button ②.



Figure 3: Audio user interfaces for different protocols (left: Screen; middle: Diagnostic; right: Auto Audio)

LITERATURE

- American Speech-Language-Hearing Association (2005): Guidelines for manual pure-tone threshold audiometry, online: http://www.asha.org/policy/GL2005-00014.htm
- IEC 60645-1: Electroacoustics Audiometric equipment Part 1: Equipment for pure-tone audiometry
- ISO 8253-1: Acoustics Audiometric test methods Part 1: Pure-tone air and bone conduction audiometry

3.1.2 Multiple-Choice Auditory Graphical Interactive Check (MAGIC)

METHODOLOGY

Play-audiometry is a well-established method in pediatric audiology. However, the test procedure, in which the child must be continuously observant, is exhausting, and the task (usually placing a peg in a pegboard) varies only slightly. In contrast, image-based and self-paced test procedures can considerably enhance the child's attentiveness. MAGIC is an image-based, self-controlled test which provides a playful alternative to conventional pure-tone audiometry. The test is especially designed for pre-school and school children but can be used for other cooperative subjects as well. Please note that for subjects who are not expected to be able to perform the test or are not able to react to the sound in a proper way, a DPOAE or ASSR Threshold measurement is advised.

Two MAGIC workflows are available:

- **MAGIC Screen** allows quick hearing screening at user-selectable frequencies and levels. If multiple levels are selected, the algorithm starts at the lowest level and proceeds to the next higher level each time the tone is not heard.
- **MAGIC Audio** allows hearing threshold determination at user-selectable frequencies with an automatically controlled bracketing algorithm. Mute tones are included for detecting incorrect inputs of the subject.

MAGIC can be performed at standard audiometric frequencies from 250 Hz to 8 kHz. Different animals represent different frequencies (cow = 250 Hz, bear = 500 Hz, elephant = 1 kHz, cat = 2 kHz, sheep = 3 kHz, mouse = 4 kHz, bird = 5 or 6 kHz, dolphin = 8 kHz – see *Figure 4*).




Figure 4: MAGIC animals (from left to right: cow, bear, elephant, cat, sheep, mouse, bird, and dolphin)

Each animal is presented in three variants: *neutral* to start the sound presentation, *healthy* and *sick* for indicating the two conditions *heard* and *not heard*, respectively (see *Figure 5*). The self-controlled workflow needs extensive explanation to the subject. An introductory story is recommended for explaining the task and for motivating the subject to be tested. The subject shall take the role of a doctor whose task is to diagnose healthy and sick animals based on whether or not they make a sound, respectively. This can also help to ease the pressure on the subject, because it is not him/her who is being diagnosed, but he/she will diagnose the animals' ability to produce sounds.

Via touch-screen the subject controls the sound presentation and the response if the sound was heard (healthy animal) or not (sick animal). The examiner explains the test procedure during the instruction phase, supervises the measurement, and assists the subject if indicated. The examiner can also intervene in the measurement procedure at any time or can do a retest at any frequency if necessary (MAGIC Audio).

For general information about conventional pure-tone audiometry please refer to section 3.1.1: Pure-Tone Audiometry (PTA).

PRACTICAL USE

Select *MAGIC* from the module selection screen. Select the preset that you would like to use. If necessary, change the parameters (e.g. test mode, frequency, screening levels (MAGIC Screen) / test level range (MAGIC Audio), stimulus type) and preset name as required. For some subjects (e.g. children) a warble tone may be more interesting than a pure-tone and hence may be preferred. Also, for subjects with a tonal tinnitus a warble-tone may distinguish the stimulus from the tinnitus more efficiently.

Make sure that a valid transducer is connected (e.g. headphone, insert earphones, bone conductor) and select the test ear (*Right+Left*: sequential measurement of right and left ear for each frequency). The measurement can be started with the introduction phase (if enabled in the settings) or the test phase.

It is recommended to start with the instruction phase. For instruction tell the subject that his/her task is to be a doctor now. Some animals will produce a sound because they are healthy and happy (shown cheering with hands up) and some animals will not produce a sound because they are sick (shown wearing a scarf). The subject in the role as a doctor shall find out which animals are healthy and which are sick. To find out whether the animal produces a tone or not the neutral animal shall be pressed (see *Figure 5* - ①). The response shall be given by pressing the respective healthy ③ or sick ② animal. A response is only accepted *after* the belly of the neutral animal has been pressed, i.e. after tone presentation. The instruction phase starts with a stimulus at a user-definable initial level (MAGIC Audio) or at a fixed level determined by the device (MAGIC screen). The second stimulus is presented with a level 10 dB lower (if first response was heard). Finally, a mute tone is presented in order to accustom the subject to the situation that also no sound may occur.




Figure 5: MAGIC user interface (left: tone presentation and response; middle/right: animal selection)

A typical instruction may follow the workflow presented in *Table 2*. A short instruction sheet is also available for download on the PATH MEDICAL website's *Learning* section.

Instruction of examiner	Task for examiner
Here you see an animal. Press the belly of the animal and you can hear if the animal makes a sound or not. The first animal makes a sound.	Show the subject how to press the belly of the neutral animal.
Press the animal's belly. Do you hear the sound?	Ask the subject to press the belly several times.

Instruction of examiner	Task for examiner
Each time you press the animal's belly it makes a sound.	The subject should notify that the sound is audible.
	If not a higher stimulus level needs to be selected (pressing <i>not heard</i> : stimulus level increases by 20 dB automatically). Repeat the instruction from the first step.
	Show the subject how to press the healthy animal corresponding to <i>heard</i> .
The next animal also makes a sound, but it is a quieter sound. Press the belly and listen.	The subject should hear this sound as well. Make sure that the subject presses the healthy animal corresponding to <i>heard</i> .
The next animal does not make a sound. You won't hear anything. Press the belly and listen.	The subject should notify that no sound is audible. If not, repeat the instruction. Show the subject how to press the sick animal corresponding to <i>not heard</i> .

Table 2: MAGIC instruction

After the instruction phase the test phase starts. Depending on age and physical/mental abilities, the test phase may be conducted completely self-controlled (i.e. subject enters response) or with assistance of the examiner. Supervision by a qualified examiner is recommended at all times.

The progress of the test procedure may be visualized by a shelf, from which the current test animal is selected (option for MAGIC Audio). During the test, please check if the subject conducts the test appropriately. If there is any hint that the subject does not understand the workflow (e.g. pressing healthy animal after mute tone presentation or pressing the screen randomly) you may restart the test or repeat the instruction phase in-between (option in settings ⑤) (MAGIC Audio). Also a previous patient response can be undone when pressing the *undo* button ④ in the footer (press the power on/off switch to make the footer visible).

In the header, information about the current test status including stimulus levels and subject's response is provided (see *Table 3*). For example, the header in *Figure 5* shows that the test is in the instruction phase with a tone presented at 50 dB HL when pressing the animal's belly. The previous stimulus, which was presented at 60 dB HL, was heard by the subject. If you suspect that the subject might take advantage of the header information, turn off the current level information in the module settings (also possible during the test for MAGIC Audio via the *settings* button (5)).

Stimulus	Test phase	Current level	MAGIC mode:	Previous level	Response
S: sine F: warble (FM tone)	I: Instruction T: Test	Value [dB HL], mute	Audio, Screen	Value [dB HL], mute	heard not heard

Table 3: MAGIC header information

At the beginning of the test phase and after finishing a test run at a frequency, the subject is allowed to select a new animal (see *Figure 5 - right*). The order in which the different animals are selected is up to the operator (i.e. tested subject or examiner).

After the test is finished, for MAGIC Screen the result is shown as a table (see *Figure 6*) with heard ④, not heard ⑥, and skipped/not tested ⑤ symbols. For MAGIC Audio, the result is shown as an audiogram. For MAGIC Audio it is possible to retest specific frequencies. A retest may be indicated if an incorrect response of the subject is shown. An incorrect response is represented by a "?" in the audiogram below the affected frequency. The number of incorrect responses is shown below the "?" ⑦. Another indication for retesting a frequency may be an outlier in the audiogram. In the MAGIC Audio result screen you can see test details by pressing in the audiogram the respective frequency. The responses at the selected frequency are shown (in this case example: 25 dB HL: heard, 5 dB HL: not heard, 15 dB HL: heard, mute: not heard, 10 dB HL: heard; test run finished and threshold determined to be 10 dB HL). The gray line marks the minimum level configured in the module settings.

Please note that during this test the footer is hidden. You may activate the footer by pressing the power on/off switch briefly.




Figure 6: MAGIC result (left: MAGIC Screen; middle: MAGIC Audio; right: MAGIC Audio detail)

STUDY RESULTS

The feasibility and reliability of MAGIC was investigated in children in a multi-center study conducted at the Department for Communication Disorders, ENT department of the University of Mainz Medical School's hospital (Germany), at the ENT clinic of Klinikum rechts der Isar, Technische Universität München (Germany), and at Cyprus Audiology Centre Nicosia (Cyprus) (Schirkonyer et al., 2010b, 2011). MAGIC tests were conducted in 108 children aged between 3;6 and 11;11 years at frequencies 0.5, 1, 2, 4 kHz (n = 82) and additionally at 0.25, 3, 6, 8 kHz (n = 26). Additionally, hearing threshold estimation by means of DPOAE I/O functions were performed in 36 children from the collective. DPOAE I/O functions were recorded at frequencies f_2 = 1, 1.5, 2, 3, 4 and 6 kHz with primary tone levels L₂ in the range from 10 to 65 dB SPL. Thresholds were estimated by linear regression analysis. Primary tone levels L_1 and L_2 were set according to the scissor paradigm. The frequency ratio f_2/f_1 was 1.2. Both MAGIC and DPOAE I/O functions were measured with Sentiero. For comparison, playaudiometry pure-tone thresholds were determined at the corresponding frequencies. There was a highly significant (p<0.001) correlation between image-based (MAGIC) and play-audiometry thresholds. The correlation coefficient (n = 1247) was r = 0.73. The histogram of the difference between image-based and play-audiometry thresholds showed normal distribution. The mean difference amounted to -1.5 ± 9.6 dB. Test time per frequency (one ear) was on average about 30 s, ranging from 14 to 91 s. Test time decreased with increasing age. Also, MAGIC pure-tone thresholds and estimated DPOAE thresholds were highly significant (p<0.005) with a correlation coefficient of r = 0.60.

LITERATURE

- Bohnert A, Schirkonyer V, Müller S, Brantzen P, Keilmann A, Janssen T (2010): Hearing screening in children using DPOAE thresholds and interactive, self-paced behavioural pure-tone thresholds (MAGIC), NHS annual meeting.
- Heller O, Opp-Enzinger AHA (2007): Self-paced in-situ audiometry, international patent WO2007009287 A2
- Janssen T, Bohnert A, Heller O, Schirkonyer V, Oswald H, Zoth P, Keilmann A (2009): Hearing diagnostics in children using DPOAE and pure-tone thresholds, EFAS annual meeting.
- Niedermeyer HP, Schirkonyer V, Lochmahr S, Janssen T (2010): Bildgestützte Tonschwellenaudiometrie und DPOAE-Schwellen in der Pädaudiologie (German: Image-based pure-tone audiometry and DPOAE thresholds in paedaudiology), DGPP annual meeting.
- Schirkonyer V, Bohnert A, Heller O, Keilmann A, Janssen T (2010a): Bildgestützte Tonschwellenaudiometrie und DPOAE-Schwellenschätzung bei Kindern (German: Image-based puretone audiometry and estimated DPOAE thresholds for children), DGA annual meeting.
- Schirkonyer V, Bohnert A, Niedermeyer HP, Thodi C, Keilmann A, Janssen T (2010b): Validierung eines Handgerätes zur bildgestützten Tonschwellenaudiometrie und DPOAE-Schwellenbestimmung bei Kindern (German: Validation of a handheld device for image-based pure-tone audiometry and DPOAE thresholds for children), DGPP annual meeting.
- Schirkonyer V, Bohnert A, Thodi C, Niedermeyer HP, Keilmann A, Janssen T (2011): Image-based and self-controlled test procedure for assessing pure-tone thresholds in children, Journal of Hearing Science 1(1).
- Thodi C, Zoth P (2010): MAGIC interactive, self-paced pure-tone audiometry for children, NHS annual meeting.
- Zoth P, Lodwig A, Schirkonyer V, Haszprunar B, Janssen T (2009): Hearing diagnostics by means of extrapolated DPOAE I/O-functions and pure-tone thresholds using a unique hand-held device, IERASG annual meeting.

3.1.3 Overview: Speech Intelligibility Tests

Analytic speech tests assess the phonetic level of speech perception. They reflect the auditory processing of individual speech sounds. These tests aim at evaluating the perception of consonants and vowels in simple words or word-like contexts. Consonants are e.g. presented within a nonsense vowel-consonant-vowel (VCV) word (e.g. SUN: aga, afa, asa). Vowels are presented in a CVC format, e.g. by using voiced plosive consonants, with the first and second consonant fixed (e.g. bid, bed, bad in an English test).

Word tests determine speech recognition ability. These tests aim at assessing the maximum achievable recognition score at an optimum speech level or test the loss in recognition score at a given speech level. In this respect they determine the effective hearing impairment for realistic speech stimuli that represent commonly used speech as e.g. monosyllabic words (e.g. bed, cow, row), spondees/two-syllable words (e.g. pancake, playground), or two-digit numbers (e.g. 21, 45, 97). Individual speech recognition values, i.e. the percentage of words which are correctly heard and repeated by the subject for each presented speech level, are displayed in the format of a speech audiogram. A specific value, e.g. hearing loss at 50 % intelligibility may be provided as speech recognition threshold. Initial speech levels may be derived from the individual hearing threshold at a specific frequency (e.g. 0.5 kHz, see e.g. PTA, MAGIC) plus a fixed offset. The average across the puretone thresholds 0.5, 1, and 2 kHz (Fletcher index) is expected to be close to the speech recognition threshold (for spondees). If the average pure-tone threshold is significantly better than the speech recognition threshold (for spondees), the possibility of a central disorder should be considered.

Speech tests may be conducted with headphones, insert earphones, and loudspeakers. Loudspeakers are especially used when testing subjects with hearing-aids or other amplifying devices which prevent the use of headphones or insert earphones.

In general, speech tests may be used to determine speech intelligibility in quiet or in noise (i.e. with ipsilateral noise) or for checking the fitting performance of hearing aids, i.e. by testing the subject with and without wearing the hearing aid. Compared to pure-tone audiometry, speech audiometry examines auditory processing on a higher intellectual level: not only hearing but also understanding.

Especially for speech tests with realistic words, word lists must be appropriate for the target subject and hence must consider the subject's age and language abilities (i.e. the expected knowledge of words). Please also consider that speech is evolving and that some previously commonly used words may get out of use (maybe dependent on age). Moreover, test results depend on the speaker (e.g. male, female voice) and the pronunciation and accentuation of words. Word tests should in general be provided in the subject's native language. If the above recommendations are not followed, the outcome may not reflect the actual hearing status. Please note that if the examiner evaluates the correctness of the response, the examiner's hearing and language abilities also influence the outcome of the test.

Senti and Sentiero devices provide various speech tests such as SUN and MATCH, and the speech test platforms UST, Speech CD Player, and Live Speech which are explained in the following sections.

3.1.4 Speech Understanding in Noise (SUN)

METHODOLOGY

SUN is a quick and self-controlled speech in noise test especially designed for screening adults and older adults for hearing disability using a set of intervocalic consonants (VCVs). Understanding speech in noise is the most common listening difficulty experienced by adults (Kramer et al., 1998). However, the test can be used for other cooperative subjects as well. The test is designed in a way that it is self-convincing, fast, and of low cognitive load.

The test aims at logatom recognition in noise, specifically at the identification of consonants and hence at the early assessment of a high frequency hearing loss. The test is fully automated. Three alternative VCVs (e.g. aFa, aGa, aSa) are displayed on the touch screen. VCVs are presented via headphones, insert earphones or loudspeakers at different signal-to-noise ratios. The subject has to select the correct item from the screen. Recordings are available with speakers of different languages (e.g. Italian, German, English²) and with different character representations (e.g. Latin, Greek, Farsi, Hindi, and Cyrillic).

Two SUN workflows are available:

- SUN Predefined allows for conducting the test with pre-determined signal-to-noise ratio (SNR) groups, i.e. the test starts at a certain SNR whereupon during the test the SNR is reduced after a fixed number of words. The sequence of words is fixed and optimized with regard to the intelligibility of the VCVs. The speech level remains constant during the entire test. The test determines a score dependent on the number of correct answers (Paglialonga et al., 2011a,b).
- SUN Adaptive allows for conducting the test with an adaptive SNR (similar to Kaernbach, 1991). The speech level is decreased after a correct answer and increased after an incorrect answer. The noise level remains constant during the entire test. An additional question mark button is present which may be pressed if the subject did not understand the logatom. The test provides an SNR threshold as result.

PRACTICAL USE

Select SUN from the module selection screen. If more than one speech test is licensed, SUN can be found in the Speech section. If necessary, change the parameters (e.g. test mode, test level, masking noise type, language, character set) as required. Make sure to select a language and character set that is familiar to the tested subject.

Before starting the test, instruct the subject about the task. On the display, three different VCV options are visualized. However, a voice pronounces only one of them. The subject shall listen to a speech sample. After the speech sample is played, the subject is asked to press the button, which corresponds to the word that has been understood. If the word has not been understood the subject shall guess and press any button (SUN Predefined) or the subject may press the question mark button (SUN Adaptive) (see Figure 7). Make sure that the subject understands the task.

² Please ask PATH MEDICAL for further information on available languages.

In order to accustom the subject to the *words* you may start with a training phase (*Training* mode needs to be enabled in settings). In the training phase, some of the logatoms from the test are played without noise. Proceed to the test phase if the subject is familiar with the task.

Before starting the training or test phase, make sure that a valid transducer is connected (headphone, insert earphones or loudspeakers) and select the test ear. The training or the test starts. Please note that during the *SUN Predefined* test, the initial responses are not counted for the final result.

Dependent on age and physical/mental abilities, the test may be conducted completely self-controlled (i.e. subject enters response) after instruction or with assistance of the examiner. Supervision by a qualified examiner is recommended at all times.



Figure 7: User interfaces for different SUN test modes (left: SUN Predefined; right: SUN Adaptive)

After the test is finished, for SUN Predefined the result is shown as a score with a traffic light status (see *Figure 8 - left*) or for SUN Adaptive as an SNR threshold together with a score and a traffic light status (see *Figure 8 - middle/right*).



Figure 8: SUN result (left: SUN Predefined; middle: SUN Adaptive with traffic light view; right: SUN Adaptive with SNR time course view)

The time course of the SNR is additionally shown in a graph when pressing on the result screen. The three-stage status light refers to the following definitions: green for hearing ability within normal

range, yellow for hearing ability which might be below normal range, red for hearing ability well below normal range. The limits between hearing range groups depend on the selected language.

Please note that during this test the footer is hidden. You may activate the footer by pressing the power on/off switch briefly.

STUDY RESULTS

A study on the effectiveness of SUN (predefined mode) for adult hearing screening was performed at the Institute of Biomedical Engineering (IsIB), Milan, Italy (Paglialonga et al., 2011a). The study was conducted on 1273 adolescents and adults (13 to 89 years) with varying degrees of audiometric thresholds including SUN (predefined mode) in Italian language and pure-tone audiometry at 1, 2, and 4 kHz as reference. Tests were performed both in an environment with low and high ambient noise. After the test, all subjects were asked to fill a questionnaire for evaluating the difficulty of the task, the test duration and the overall rating. The main outcome of the study was that SUN is suited for adult hearing screening because of the following detailed findings: The overall SUN result was in line with pure-tone audiometry with a good correlation between the three SUN categories (traffic light status) and three specifically defined PTA classes. The test performance was similar across all subjects independent of age. The test performance was also not influenced by ambient noise levels up to 65 dB(A) so that the test may also be performed in a non-clinical setting, where ambient noise is not typically controlled (e.g. hearing aid providers). The test time was very short with two minutes on average for both ears. Even older adults typically managed to conduct the test within one minute per ear. This is deemed as an important factor for a hearing screening test because inattentiveness and fatigue are likely to increase with increasing test time. The cognitive load was low and the acceptance of the test was very high. About 85 % of the subjects considered the test to be easy or slightly difficult. 95 % of the subjects judged the test duration to be short or fair. More than 90 % of the subjects rated SUN pleasant or neutral.

A further study (Paglialonga *et al.*, 2013) extended the above data by testing SUN (predefined mode) in more than 6000 subjects. SUN was found to be very reliable with a good correlation to pure-tone audiometry and self-reported hearing handicap. The test's sensitivity and specificity to identify disabling hearing impairment were 84 % and 75 %, respectively. Test time was on average lower than 1 minute per ear.

LITERATURE

- Grandori F, Tognola G, Paglialonga A (2010): The SUN test Speech Understanding in Noise: A new test for screening hearing ability in adults, EUHA annual meeting.
- Kaernbach C (1991): Simple adaptive testing with the weighted up-down method, Perception & Psychophysics 49, p. 227-229.
- Kramer SE, Kapteyn TS, Festen JM (1998): The self-reported handicapping effect of hearing disabilities, International Journal of Audiology 37(5), p. 302-312.
- Paglialonga A, Tognola G, Grandori F (2011a): SUN-test (Speech Understanding in Noise): a method for hearing disability screening, Audiology Research 1:e13.
- Paglialonga A, Tognola G, Grandori F (2011b): Pilot initiatives of adult hearing screening in Italy, Audiology Research 1:e17.

- Paglialonga A, Tognola G, Grandori F, Chisolm TH, Deggouj NAD, Fenech AJ, Janssen T, Li X, Walger M (2012): The SUN test (Speech Understanding in Noise) in different languages: test development and fine-tuning, AHS annual meeting.
- Paglialonga A, Grandori F, Tognola G (2013): Using the speech understanding in noise (SUN) test for adult hearing screening, American Journal of Audiology 22(1), p. 171-174.

3.1.5 Mainzer Audiometric Test for Children (MATCH)

METHODOLOGY

MATCH is an interactive, image-based speech-intelligibility test, which can be performed in quiet and in noise. The test is based on the Mainzer Kindersprachtest I (German speech test for children), which uses one- and two syllable words presented on image cards. Disadvantages of the Mainzer Kindersprachtest I were that the number of words was too low, words and images were not up to date, and words did not match the typical vocabulary of a two-year-old child.

MATCH consists of 26 words corresponding to the typical vocabulary of a two-year-old child (see Figure 9). Although the test is especially designed for children, it can be used for other cooperative subjects as well. The words are selected from recent studies on the vocabulary of young children (Suchodoletz and Sachse, 2008). The phoneme distribution of the German test shows a good correspondence with conversational German (Schiel, 2010; Schirkonyer et al., 2014). The test items are figuratively presentable and not gender-specific. Simple pictures appropriate for children represent the words (e.g. duck, milk, car, etc.) so that the subject does not need reading skills. The subject just needs to point at the suspected correct image. The test is designed as closed-set test (i.e. limited number of alternatives available for selection) because an open-set test is expected to be too difficult for smaller children. Four pictures are presented simultaneously on the touch-screen together. One of them belongs to the spoken word. If the word has not been understood, the subject does not need to guess but may press the "?" button. The sound presentation is started by pressing on a "magic box". After the subject has given the response the test continues until all items have been presented. The test is completely randomized, i.e. the sequence of words and the alternatives are randomized so that the test can be repeated with the same subject. It is possible to conduct the test with a subset of the complete word list. For example, if the child does not know some words, these words can be removed from the test.



Figure 9: MATCH images (from left to right: 1) monkey, apple, eye, car, bear, ball, tree, 2) bed, boat, butter, eggs, duck, cucumber, hair, 3) bunny, trousers, light, milk, mouth, nose, ear, 4) grandpa, horse, door, clock, water)

MATCH is also available with additional ipsilateral noise with a fixed noise level of 65 dB. This mode assesses speech intelligibility in noise following the same procedure as without noise. It was shown that the threshold of speech intelligibility in noise is more sensitive for assessing a hearing loss (Leensen *et al.*, 2011). The test with ipsilateral noise is suitable for older children and adults.

Two MATCH workflows are available:

- **MATCH Fixed** allows for conducting the test with a fixed speech level. The test determines a score dependent on the number of correct answers.
- **MATCH Adaptive** allows for conducting the test with an adaptive speech level presentation. The speech level is decreased after a correct answer and increased after an incorrect answer. The test determines a 71.4 % speech recognition threshold in quiet (without ipsilateral noise) or in noise (with ipsilateral noise).

Normative data for children have been collected and presented by Bohnert *et al.* (2013), Zoth *et al.* (2013) and Schirkonyer *et al.* (2014). A more detailed explanation is available in the MATCH Quick Guide (100800-24).

PRACTICAL USE

Select *MATCH* from the module selection screen. If more than one speech test is licensed, *MATCH* can be found in the *Speech* section. If necessary, change the parameters (e.g. test mode, test level (MATCH Fixed) / test level range (MATCH Adaptive), masking noise, language, items) as required. Make sure to select a language that is familiar to the tested subject³. Select or deselect masking noise if you prefer to test speech intelligibility in noise or in quiet, respectively.

Make sure that a valid transducer is connected (headphone, insert earphones, or loudspeakers) and start the test by selecting the test ear. When starting the test, instruct the subject about the task. Tell the child that there is a magic box on the screen of the device and that the task is to find out what is inside the box. A typical instruction may follow the workflow presented in *Table 4*. The instruction is also explained in the online help of the device.

Instruction of examiner	Task for examiner	
The magic box may contain different things.	Show the subject how to press the box with the green question marks (see <i>Figure 10</i> - ①).	
Press the magic box and listen carefully what is inside the box this time.	Ask the subject to press the magic box. Afterwards, the response alternatives appear (see Figure 10 - ③).	
Did you hear what is inside the box? If you do not hear what is inside the box press the green question mark at the bottom of the screen. If you hear it, show me what is inside the box.	The subject shall notify if the word has been understood or not. Show the subject the green question mark button and show how to press the correct response image.	

Table 4: MATCH instruction

-

³ Please ask PATH MEDICAL for further information on available languages.

Please note that the response input is only possible after the word has been played. During the word presentation the *question mark* button ② is grayed out and no response is possible.



Figure 10: MATCH user interfaces (left: elicit stimulus; right: response alternatives)

In the header, information about the current test status including stimulus levels and subject's response is provided (see *Table 5*). For example, the header in *Figure 10* (right image) shows that the word is presented at 40 dB HL. The correct answer is at the bottom right image (door). One item has already been tested. The previous word has been correctly selected by the subject.

Speech level	Correct word	Number of items	Response
Value [dB HL] 4	Image position index (5): 0 2 1 3	Number of tested items 6	correct item selected incorrect item selected

Table 5: MATCH header information

If you suspect that the subject might take advantage of the header information, turn off the header information in the settings before starting the test.

Dependent on age and physical/mental abilities, the test may be conducted completely self-controlled (i.e. subject elicits stimulus and enters response) after instruction or with assistance of the examiner. Supervision by a qualified examiner is recommended at all times.

During the test, there are some specific features available when connecting the patient response button. The patient response button shall in this case be operated by the examiner. A speech sample may be repeated by the examiner by pressing the button briefly and a previous response can be undone (e.g. if the subjects would like to correct the response after response input) when pressing the button for about 2 s.

After the test is finished, for MATCH Fixed the result is shown as a score (see *Figure 11*) and for MATCH Adaptive the result is shown as a speech audiogram (estimated discrimination function) with a speech recognition threshold (see *Figure 12*). Additionally, the time course of speech level and reaction time as well as the test protocol are shown when pressing on the result screen. The test protocol shows a list of all presented words, the presented alternatives, and the subject's response.

The time course graph can help to evaluate the reliability of the test, whereas the test protocol may provide information on problems e.g. with discrimination of specific phonemes.



Figure 11: MATCH Fixed result (left: score in quiet; right: score in noise)



Figure 12: MATCH Adaptive in quiet result (left: speech audiogram and discrimination function; middle: time course of speech level (blue symbols) and reaction time (green line); right: test protocol with word selections)



Figure 13: MATCH Adaptive in noise result (left: speech audiogram and discrimination function; middle: time course of speech level (blue symbols) and reaction time (green line); right: test protocol with word selections)

STUDY RESULTS

A study (Bohnert et al., 2013; Schirkonyer, 2013; Schirkonyer et al., 2014) for evaluating the test performance of MATCH with children was conducted at the Department for Communication Disorders, ENT-department of the University of Mainz Medical School's hospital (Germany). In the study, pure-tone audiometry and MATCH was performed with Senti Desktop in 111 children (213 ears) from 2;6 to 6;9 years. Thereof, 157 ears exhibited normal hearing, whereas 56 ears were hearing-impaired. Some children exhibited a unilateral hearing impairment. Normal hearing was defined in the study as screening pass at 30 dB HL at 0.5, 1, 2, 4, and 6 kHz. If a fail occurred, hearing thresholds were determined by extrapolated DPOAE I/O-functions. All tests were performed in a quiet room with ambient noise < 50 dB(A) or in a booth. Data from all subjects were split up into three age groups (AG1: 2;6-4;3 years, AG2: 4;3-5;6 years, AG3: 5;5-6;9 years). Speech recognition thresholds at 50 % were 29.4 \pm 1.2 (AG1), 27.8 \pm 0.9 (AG2) and 25.4 \pm 0.8 dB HL (AG3), whereas speech recognition thresholds at 71.4 % were 38.2 ± 1.0 (AG1), 34.9 ± 0.7 (AG2) and 31.1 ± 1.0 dB HL (AG3), i.e. the speech recognition thresholds slightly decreased with increasing age. The median of the test duration amounted to 4.0 (AG1), 3.4 (AG2), and 3.1 minutes per ear (AG3), i.e. test duration decreased with increasing age. Test-retest stability (n = 79) was good with a correlation of r = 0.89. When comparing average pure-tone thresholds at 0.5, 1, and 2 kHz to MATCH 71.4 % speech recognition thresholds (n = 55) there was a good correlation of r = 0.81.

LITERATURE

- Bohnert A, Schirkonyer V, Harmuth C, Leonhard K, Janssen T, Keilmann A (2013): Was ist neu am Neuen Mainzer Kindersprachtest? (in German: What is new with the Mainzer speech test for children), ADANO annual meeting.
- Leensen MC, de Laat JA, Snik FM, Dreschler WA (2011): Speech-in-noise screening tests by internet, part 2: improving test sensitivity for noise-induced hearing loss, International Journal of Audiology 50(11), p. 835-848.
- Schiel F (2010): BAStat: New statistical resources at the Bavarian Archive for Speech Signals, In: Proc. of LREC, paper 277.
- Schirkonyer V, Bohnert A, Keilmann A, Janssen T (2014): MATCH ein neuer Kindersprachtest (in German: MATCH a new speech test for children), DGA annual meeting.
- Suchodoletz W, Sachse S (2008): Handbuch: Sprachbeurteilung durch Eltern Kurztest für die U7 (SBE-2-KT) (German: Manual: language assessment by parents quick test for the U7 (SBE-2-KT)), online: http://www.kjp.med.uni-muenchen.de/download/SBE-2-KT-Handbuch.pdf.
- Zoth P, Bohnert A, Schirkonyer V, Keilmann A, Janssen T (2013): A new interactive & self-paced speech test for children, EFAS annual meeting.

3.1.6 Universal Speech Test (UST)

METHODOLOGY

The Universal Speech Test (UST) provides a platform for speech tests which follow a typical workflow of presenting words (organized in word lists) at different speech levels with the aim of determining a 50 % speech recognition threshold and optionally a maximum discrimination level, which amounts at the best to 100 %. The UST offers three test phases: (1) voice hearing phase (optional) for determining the level at which a word or in general human voice can be heard but not understood, (2) word understanding phase (optional) for determining a level at which the subject begins to understand a word and (3) the test phase which allows for determining speech recognition scores at different speech levels. The optional pre-test phases can be used for determining a starting level for the test phase. Alternatively, audiometric threshold data (e.g. average hearing threshold across 0.5, 1 and 2 kHz) may be used. UST is available with optional ipsilateral noise for assessing speech intelligibility in noise. Contralateral noise is available for masking purposes, e.g. for patients with unilateral hearing loss or distinct difference in hearing threshold between left and right ear. The clinical application and target patient group depend on the speech material used. Typical speech material includes single syllable words, two syllable words/spondees, and two-digit numbers⁴.

As an example the Freiburger speech test (Halbrock, 1953), a well-established German speech test for adults, is explained in the following. Please note that other speech tests may use different word types and may aim at different target speech recognition values. The Freiburger speech test uses monosyllabic nouns that represent commonly used speech (e.g. Ring, Spott, Farm, Hang, ...) and twodigit numbers with mainly four-syllables (e.g. 98, 22, 54, ...). Word lists are based on the most commonly used words at the time at they were recorded. Lists were prepared in such a way that the words exhibit a similar phonetic distribution. The test comprises 20 word lists with 20 monosyllabic nouns each and 10 word lists with 10 two-digit numbers each (see DIN 45621-1). The Freiburger speech test can be split up in two separate subtests. Both are typically conducted without ipsilateral noise in order to determine speech understanding in quiet. The monosyllabic word test is intended for determining individual speech recognition values, i.e. the percentage of words which are correctly heard and repeated by the test person for multiple speech levels (65, 80, 95, 110 dB SPL), which are displayed in the format of a speech audiogram. The intended use of the two-digit number test is to determine a hearing loss for numbers at 50 % intelligibility. Numbers are presented at two speech levels. The lower one corresponds to the hearing loss of the subject at 4 kHz determined by puretone audiometry plus 20 dB. The higher one has to be set to a level that is 5 dB higher than the initial level. The initial speech level corresponds to the individual hearing loss at 500 Hz from the pure-tone audiogram plus 20 dB.

For normal-hearing subjects, monosyllabic nouns usually require 10 to 20 dB higher levels to achieve the same speech recognition score as for two-digit numbers (Brinkmann and Richter, 1997). Hearing-impaired subjects may not reach a 100 % speech recognition score even at high speech levels. Normative data for monosyllabic words and two-digit numbers is presented in DIN 45626-1.

PRACTICAL USE

Select the UST subtest (e.g. *Freiburger*) from the module selection screen. If more than one speech test is licensed, the respective UST subtest can be found in the *Speech* section. If necessary, change

⁴ Please contact PATH MEDICAL for an overview of available speech tests under the UST framework.

the parameters (e.g. ipsilateral noise or contralateral masking noise, noise level, voice hearing / word understanding phase) as required. Make sure to select a UST subtest with language and word lists that are matching the tested subject's language abilities.

Make sure that a valid transducer is connected (headphone, insert earphones, or loudspeakers) and start the test by selecting the test ear. Before the test, the subject must be instructed about the task. During the voice hearing phase the subject shall respond if a voice is heard or not. During the word understanding phase the subject shall respond if a word has been understood or not. The determined word understanding level is marked in the speech audiogram as lowest speech level with a score of zero. During the test phase, the subject shall repeat the word. The examiner needs to evaluate if the word was understood correctly or not. Please note that the test must be conducted by a qualified and normal-hearing examiner. The examiner needs to elicit the stimulus and enter the response. The test is *not* intended to be self-controlled by the subject.




Figure 14: UST user interface (top left: word list selection; top right: main user interface during voice hearing phase; bottom left: level selection during test phase; bottom middle: main user interface during test phase; bottom right: main user interface with visible footer)

For each pre-test phase and during the test for each test run (i.e. for each tested speech level), a word list needs to be selected (see *Figure 14* - ①). Already selected word lists are grayed out but can be selected again. Make sure that previously presented words are not already known to the subject. Therefore, during the test phase, select for each test run (i.e. for each tested speech level) a new word list. After selecting the word list, adjust the speech level to your needs by pressing the *up/down*

buttons ②. In most tests, typically multiple levels are selected for determining speech recognition scores from very low (around 0 %) to very high (around 100 %).

A word is played by pressing the large *word* button ④. If you would like to select another word from the list, you can browse through the word list via the triangular *up/down* buttons above/below the list. The smaller word buttons show the previous ③ and next ⑤ word in the list. After the word has been played and the subject has answered, enter the patient response (correctly repeated: ⑦, incorrectly repeated: ⑧). Already tested words are colored green (correctly repeated) or red (incorrectly repeated). The test can be stopped by pressing the *stop* button ⑥. Any pre-test phase or test run at the current speech level can be cancelled by pressing the *forward* button ⑨. The current result screen can be displayed by pressing the *score* button ⑩. During the test masking can be turned on or off by pressing the *masking* button ⑪ in the footer.

After the test is finished, the result (see *Figure 15*) is shown as a graphical and tabular speech audiogram together with a 50 % speech recognition threshold. Toggle the views by pressing the result screen. There are specific result adaptations for specific tests (e.g. *Freiburger*: see *Figure 16*).




Figure 15: Typical UST result (left: speech audiogram graph; right: speech audiogram table)




Figure 16: Specific Freiburger result (left: speech audiogram graph; right: speech audiogram table)

Please note that during this test the footer is hidden. You may activate the footer by pressing the power on/off switch briefly.

There are two adaptations of the basic UST procedure as described above for performing specific speech tests.

The first adaptation offers a user interface which allows presenting words from one or multiple word lists at different speech levels (see *Figure 17*) for a simplified determination of a speech threshold. This workflow is available e.g. for the *Spondees Children* or *Spondees Adult* speech tests. After selecting the word list ①, the user changes the speech level until a speech threshold is obtained. The determined speech threshold ③ is the level which is set by the examiner when finishing the test with the *stop* button ②.



Figure 17: UST user interface adaptation for simplified determination of speech thresholds (left: list selection; middle: test with level selection for determination of the speech threshold; right: result view)

The second adaptation is similar to the standard UST procedure but provides level selection within the main test screen (and not between word list selection and test screen as for the standard UST procedure) so that the level can be changed while going through a single word list (see *Figure 18*). This workflow is available e.g. for the *NU-6*, *CID W-22*, *PBK-50*, *Maryland CNC*, and *Mots français CAD* speech tests. The test result is initially presented as a table but can be switched to a graph by pressing the result screen.




Figure 18: UST user interface adaptation for determination of speech thresholds with level control during testing a single list (left: list selection; middle: test with level selection and scoring; right: initial result view as table)

LITERATURE

- American Speech-Language-Hearing Association (1988): Determining threshold level for speech, online: http://www.asha.org/policy/GL2005-00014.htm
- Brinkmann K, Richter U (1997): Ensuring reliability and comparability of speech audiometry in Germany, In: Speech Audiometry (editor: Martin M), Whurr Publishers Ltd., p. 106-130.
- Causey GD, Hood LJ, Hermanson CL, Bowling LS (1984): The Maryland CNC Test: Normative Studies, Audiology 23, p. 552-568
- DIN 45621-1: Sprache für Gehörprüfung Teil 1: Ein- und mehrsilbige Wörter (in German: Word lists for recognition tests part 1: monosyllabic and polysyllabic words).
- DIN 45626-1: Tonträger mit Sprache für Gehörprüfung Teil 1: Tonträger mit Wörtern nach DIN 45621-1 (in German: Sound carrier with speech for recognition tests part 1: sound carrier with word lists in accordance with DIN 45621-1 (recording 1969)).
- Hahlbrock KH (1953): Über Sprachaudiometrie und neue Wörterteste (in German: On speech audiometry and new word-based tests). European Archives of Oto-Rhino-Laryngology 162, p. 394-431.
- Heckendorf AL, Wiley TL, Wilson RH (1997): Performance Norms for the VA Compact Disc Versions of CID W-22 (Hirsh) and PB-50 (Rush Hughes) Word Lists, J Am Acad Audiol 8, p. 163-172
- Hillock-Dunn A (2015): Pediatric Speech Recognition Measures What's Now and What's Next!, http://www.audiologyonline.com
- Picard M (1997): Speech Audiometry in French speaking Quebec, Revue d'Orthophonie et d'audiologie 21(4), p. 301-311
- Schoepflin JR (2012): Back to Basics: Speech Audiometry, http://www.audiologyonline.com
- Stoppenbach DT, Craig JM, Wiley TL, Wilson RH (1999): Word Recognition Performance for Northwestern University Auditory Test No. 6 Word Lists in Quiet and in Competing Message, Journal of the American Academy of Audiology 10, p. 429-435
- Thomson S (2012): Comparison of Word Familiarity: Conversational Words vs. NU-6 list Words, http://www.audiologyonline.com

3.1.7 Speech CD Player

METHODOLOGY

The Speech CD Player provides an alternative platform for speech tests which follow a typical workflow of presenting words organized in word lists that are grouped in a CD track at one or multiple speech levels. The test can be performed in quiet or in noise. Contralateral masking noise can be applied optionally. The Speech CD Player provides the ability to select a CD track and play the words from this track with start, stop, and pause functions. Scoring can be performed for each presented word. The clinical application and target patient group depend on the speech material used⁵.

PRACTICAL USE

Select the Speech CD Player subtest from the module selection screen. If more than one speech test is licensed, the respective subtest can be found in the *Speech* section. Make sure to select a subtest with language and word lists that are matching the tested subject's language abilities. In the following the workflow is explained as an example for the Dutch NVA speech test (see *Figure 19*).




Figure 19: NVA user interface (top left: list selection; top middle: track selection; top right: stimulus output configuration; bottom left to right: test user interface at different times within a track)

-

⁵ Please contact PATH MEDICAL for an overview of available speech tests under the Speech CD Player framework.

Make sure that a valid transducer is connected (headphone, insert earphones, or loudspeakers). Before starting the test, the subject must be instructed about the task. During the test, the subject shall repeat the word. The examiner needs to evaluate if the word has been understood correctly or not. Please note that the test must be conducted by a qualified and normal-hearing examiner. The examiner needs to control the stimulus playback if required and enter the response. The test is *not* intended to be self-controlled by the subject.

Initially, the appropriate speech material needs to be selected by the examiner. Each button ① selects a group of tracks from the CD (e.g. with words for different age groups). Afterwards, the respective track (2) from the track group can be selected. In order to set up the stimulus output configuration, press the settings (3) button. In the settings menu, the user can configure the output for the two stimulus channels. For each channel, the speech level can be set by pressing the up/down buttons 4. If you wish to change both levels simultaneously you can lock the speech levels of channel 1 and 2 by pressing the lock button (5). For each channel, the stimulus type (6) and the ear side (7) can be configured. After the track selection, the test starts automatically by playing the words from the chosen track as the CD is running. During the test, the currently played word (8) and the elapsed time are displayed. As soon as activated, at any time the scoring buttons (10) may be pressed (for NVA to score the number of correct phonemes). The presentation can be interrupted by pressing the pause button (2). If you wish to repeat the current word or proceed to the next word press the back 11 or forward 13 buttons. By pressing the back button twice you can rewind to the previous word. The score (14) is displayed according to the user input. At the end of a track, a single data point for the speech audiogram is stored. When leaving the test the result screen is shown as a graphical and tabular speech audiogram together with a 50 % speech recognition threshold and a maximum score (see Figure 20). Toggle between the views by pressing the graph / table.




Figure 20: NVA result (left: speech audiogram graph; right: speech audiogram table)

3.1.8 Live Speech

METHODOLOGY

The Live Speech module provides a basic platform for speech tests with live speech presented by the examiner via microphone. The output level can be balanced by means of a volume unit (VU) meter. The user is responsible for the correct selection, pronunciation and balancing of the words as well as for the result evaluation. In general, there are known deficits of this method related to reproducibility across presentations (e.g. due to varying pitch, voice quality, speed of production of words by the examiner). The Live Speech module offers three modes: Speech Detection Threshold (SDT), Speech Recognition Threshold (SRT), and Word Recognition (WR).

The speech recognition threshold or speech reception threshold is the minimum speech level at which the patient recognizes 50 % of the speech material. Recognition means that the subject can reproduce the presented speech material correctly (either by repeating the presented speech material or by selecting the test item from a closed set of choices). The speech detection threshold or speech awareness threshold is the minimum speech level at which the patient is able to notice the presence of speech 50 % of the time. The patient does not need to identify the presented material as speech, but must indicate awareness of the presence of sound. Word recognition is similar to speech recognition threshold, but allows setting a speech threshold at a recognition rate differing from 50 %.

The clinical application of Live Speech is the determination of speech recognition or detection thresholds in quiet in cooperative subjects for audiological diagnostics. The clinical application and intended patient target group depend on the applied speech material.

PRACTICAL USE

Select the *Live* speech module from the module selection screen. If more than one speech test is licensed, the respective speech module can be found in the *Speech* section.

Make sure that a valid transducer is connected (headphone or insert earphones). Before starting the test, the subject must be instructed about the task. During the test, the subject shall repeat the presented test item (SRT, WR) or indicate that a sound has been detected (SDT). The examiner needs to evaluate if the presented test item has been understood correctly (SRT, WR). Please note that the test must be conducted by a qualified and normal-hearing examiner. The examiner needs to produce live speech and enter the response. The test is *not* intended to be self-controlled by the subject.

The *Live* speech module provides several options for configuring the test (see *Figure 21*). First, please select the test mode (SDT: Speech Detection Threshold, SRT: Speech Recognition Threshold, WR: Word Recognition) by pressing the *test mode* button ② until the favored test mode appears. The test ear (right, left, both) can be selected by pressing the respective *ear* checkboxes ①.

Before starting the live speech presentation, the microphone input level may need to be adjusted: make sure that the *input/output* button (§) shows *MIC*. If necessary, turn on the microphone by pressing the *microphone on/off* button (7) and speak into the microphone with a speech intensity that you will use for the live speech presentation. Adjust the microphone input with the *level* buttons (4) so that the VU meter indicator (3) settles close to 0 and does not continuously reach the red area. Press the *level up* button to increase and the *level down* button to decrease the sensitivity of the microphone. You may skip the microphone level adjustment procedure if your speech intensity is

appropriate for the given microphone level setting. If the microphone level is correct you can proceed with setting the speech level: press the *input/output* button (8) to show the headphone icon and use the *level* buttons (4) to change the speech level as needed.



Figure 21: Live Speech user interface (left: SDT user interface with microphone turned off; middle: SDT user interface with microphone turned on and microphone level adjustment; right: WR user interface after accepting the determined threshold)

During the test, make sure that when you speak into the microphone the VU meter indicator mainly remains in the blue area close to zero. If the indicator continuously reaches the red area during the presentation of live speech, please adjust the microphone level as described above. During SDT or SRT mode the speech level needs to be adjusted by the examiner until the respective threshold is determined. In WR mode, correct and incorrect answers can be additionally counted by pressing the OK (9) or not OK button (11), respectively. The score and percentage are displayed (10). Please note that these values are not remembered when changing the speech level. If the speech threshold has been determined, press the accept threshold button (5). The currently selected speech level and for WR also the score are taken over and displayed in the test mode button (2). Press the stop button (6) to finish the test.

Results of all three modes are saved. The result screen shows a speech audiogram including the word recognition score and the determined speech detection and recognition threshold (see *Figure 22*).




Figure 22: Live Speech result (with word recognition score in speech audiogram, speech detection threshold and speech recognition threshold)

3.1.9 Munich Auditory Screening Test for Processing Disorders (MAUS)

METHODOLOGY

MAUS is a German screening test for identifying subjects with central auditory processing disorder (APD). The test was especially designed and validated for children from 6 to 11 years but may be used for other cooperative subjects as well. The test is supposed to be conducted after exclusion or treatment (e.g. hearing aid, surgery) of a peripheral hearing disorder. In case of a positive result regarding APD, further diagnostics may be advised for the differentiation to other dysfunctions, such as cognitive dysfunctions, attention deficit disorders, general dysfunctions of short-term memory. MAUS is only available in German language (recordings: © Westra Elektroakustik GmbH 2003/2004). MAUS includes three subtests: sequence of syllables, words in noise, and phoneme differentiation and identification. These subtests allow for maximum sensitivity for detecting APD (see section STUDY RESULTS). The sequence of syllables subtest includes 18 items with six three-syllable-words (e.g. muwage = mu-wa-ge), six four-syllable-words and six five-syllable-words. All items are meaningless words. The occurrence of vowels and consonants noticeable for subjects with ADP is balanced. The word in noise test includes 12 items per ear with speech-simulating ipsilateral noise. In contrast to the other subtests, for which items are presented binaurally, all word in noise items are presented monaurally. The phoneme differentiation and identification test includes 13 consonantvocal and 10 consonant-consonant-vocal pairs of meaningless syllables. The entire speech material is presented at a fixed level of 65 dB HL. The test duration is about 15 minutes. For further details please refer to the section STUDY RESULTS and to the MAUS manual (Nickisch et al., 2004).

PRACTICAL USE

Select *MAUS* from the module selection screen. If more than one speech test is licensed, *MAUS* can be found in the *Speech* section. Make sure that a valid transducer is connected (e.g. headphone, insert earphones). The test starts with the sequence of syllables subtest and proceeds to the subtests words in noise and phoneme differentiation and identification (see *Figure 23*). For each subtest, instruct the subject about the tasks as described in *Table 6*. The instructions are also displayed on the device screen before the subtest starts.

Instruction of examiner	Task for examiner				
Subtest 1: Sequence of syllables					
You will hear magic words. Listen carefully and repeat what you heard.					
Subtest 2: Words in noise					
Now you will hear words while there is noise. Listen carefully and repeat what you heard.					
Subtest 3: Phoneme differentiation and identification					
Now I will say two meaningless words. Listen carefully and tell me if these words are exactly the same or not! Both cases may occur. After that, please repeat the words in the correct order. Let's give it a try. Mi-Mi, these words are exactly the same. Which words did I say? [if correct:] Great.	Wait for subject's response and check if the subject understands the task for similar syllables.				

[if wrong:] Listen carefully again, I will tell you the words once more [repeat instruction]	
Now listen again: Bo-Sa, these words sound different. Can you hear that? Which words did I say? Try it yourself now: Li-Pa. Is that similar or different?	Wait for subject's response and check if the subject understands the task for different syllables.
[if correct:] Great.	
[if wrong:] Listen carefully again [repeat instruction]	

Table 6: MAUS instruction

Please note that for the first subtest (sequence of syllables) the subject has two tries, i.e. if the sequence of syllables was not repeated correctly the first time, the item is automatically replayed and the subject may repeat the word a second time. For the other subtests, all test items are played only once. For further qualitative analysis, it is recommended to write down incorrect responses on a separate protocol sheet. The test is *not* intended to be self-controlled by the subject.




Figure 23: MAUS user interface for the different subtests (left: sequence of syllables; middle: words in noise; right: phoneme discrimination and identification)

The current test item is displayed in a box on the screen ①. The box containing the test item is highlighted green as long as the test item is presented. The response status can be entered by selecting the OK (correctly repeated) ② or the Error (incorrectly repeated) ③ radio button. For subtest phoneme discrimination and identification additionally the correctness of the discrimination must be entered (Yes: correctly discriminated ⑦; No: incorrectly discriminated ⑧). The initial test item for each subtest is played automatically. After the response status is entered, the next test item can be started by pressing the Play ⑤ button. You can cancel the current subtest and proceed with the next subtest by pressing the Play ⑤. The test can be stopped by pressing the Play button ⑥. The test can be stopped by pressing the Play button ④.

The result shows a list of scores with color code (see *Figure 24*). The scores are calculated automatically from the input of the examiner. The maximum scores are for sequence of syllables 36 points, for 'words in noise' 12 points per ear, and for 'phoneme discrimination and identification' 22 points each. The score is mapped to an age-dependent T-value range which is represented by colored areas (<30: extremely substandard, 30-40: substandard, 41-60: normal, 61-70: surpassing, >70

extremely surpassing). Please note that even if all responses were correct, for some age groups, the score is not mapped to the best T-value range. If all results are normal or better, the probability for an APD is very low. If at least one result is substandard or worse, an APD may possibly be present. Please note that the test is a screening test, which may deliver false positive results. So, if MAUS delivers a conspicuous result, further diagnostics (e.g. intelligence diagnostics, neuro-pedaudiologic diagnostics, objective audiometric diagnostics as e.g. ABR) is advised.



Figure 24: MAUS result

STUDY RESULTS

In a pre-study by Nickisch and Oberle (2002), 79 children from 6 to 12 years (average: 8,8 years) with normal peripheral hearing (pure-tone thresholds ≤ 20 dB HL) and clinical indication for APD including two or more APD-positive results with the following psycho-acoustic tests: Mottier test from Zurich reading test, Heidelberg sound differentiation test (HLAD) test with sound differentiation, analysis and identification, speech audiometry in noise, number sequences, dichotic listening test, directional hearing, time-compressed speech, binaural fusion, Patsy test, gap detection, and joining phonemes. The incidence rate of results atypical for normal-hearing subjects was analyzed. A combination of tests was investigated with the aim of getting maximum sensitivity regarding APD detection. A combination of three tests was found to deliver maximum sensitivity: speech audiogram in noise, Mottier test, and phoneme differentiation from HLAD. Normative data (Nickisch *et al.*, 2004b) was derived from 356 school children (189 male, 167 female) with an average age of 8.76 years and with normal peripheral hearing.

A clinical study (Nickisch *et al.*, 2004a) with 52 children (36 with APD, 16 without APD) from 6 to 12 years (average age of 9.5 years) with normal peripheral hearing was performed. All children exhibited normal peripheral hearing. APD was diagnosed if at least two of the other tests mentioned above exhibited an APD-positive result. MAUS was conducted for comparison. The correlation between the MAUS subtests and the other test methods was examined and amounted to r = 0.78 (MAUS: sequence of syllables \leftrightarrow Mottier test), r = 0.51 (MAUS: words in noise \leftrightarrow speech audiometry in noise), r = 0.76 (MAUS: phoneme differentiation \leftrightarrow HLAD: sound identification). All correlations were highly significant. The agreement rate of APD screening results was compared between the MAUS subtests and the other test methods and amounted to 71 % (MAUS: sequence of syllables \leftrightarrow Mottier test), 75 % (MAUS: words in noise \leftrightarrow speech audiometry in noise), 73 % (MAUS: phoneme differentiation \leftrightarrow HLAD: sound differentiation), and 87 % (MAUS: phoneme identification \leftrightarrow HLAD: sound

identification). The sensitivity of MAUS amounted to 97 %. In an extended study (Heuckmann *et al.*, 2006) the sensitivity was validated with 132 school children (80 with APD, 52 without APD) from 6 to 12 years (average age: 9,5 years) and amounted to 96 %.

Retest reliability of MAUS was investigated in a study (Heuckmann *et al.*, 2005) with 62 school children from 8 to 10 years. Therefore, MAUS was conducted twice with a time lag of two to six weeks from the first to the second test. The overall retest reliability for MAUS was r = 0.75 with the subtests ranging from r = 0.64 (sequence of syllables) to r = 0.75 (phoneme differentiation).

LITERATURE

- Heuckmann C, Massinger C, Burger T, Nickisch A (2005): Münchner Auditiver Screeningtest für Verarbeitungs-und Wahrnehmungsstörungen (MAUS) - Retest-Reliabilität und Sensivität (in German: Munich Auditory Screening test (MAUS) for Perception Disorders - reliability and sensitivity), DGPP annual meeting.
- Heuckmann C, Massinger C, Burger T, Kunze S, Nickisch A (2006): Münchner Auditiver Screeningtest für Verarbeitungs- und Wahrnehmungsstörungen (MAUS) – Untersuchungen zur Kriteriumsvalidität (in German: Munich Auditory Screening test (MAUS) for Perception Disorders – Analysis of vriteria validity), DGPP annual meeting.
- Nickisch A, Oberle D (2002): Analyse von Testprofilen bei auditiven Verarbeitungs- und Wahrnehmungsstörungen (in German: Analysis of test profiles for auditory processing disorders), In: Aktuelle phoniatrisch-pädaudiologische Aspekte (in German: Current phoniatric-paedaudiological aspects) (editors: Kruse E, Gross M), Median-Verlag, p. 327-331.
- Nickisch A, Heuckmann C, Burger T (2004): MAUS manual, Westra Elektroakustik.
- Nickisch A, Heuckmann C, Burger T, Massinger C (2006): Münchner Auditiver Screeningtest für Verarbeitungs- und Wahrnehmungsstörungen (MAUS) (in German: Münchner screening of auditory perception disorders (MAUS)), Laryngorhinootologie, 85, pp.253–259.

3.1.10 Bochum Auditory Speech Discrimination Test (BASD)

METHODOLOGY

The Bochum Auditory Speech Discrimination Test (BASD) provides an audiometric test battery for detecting speech development deficits and related auditory processing disorders (APD), i.e. functional deficits that are not caused by peripheral hearing disorders or attention deficit disorders. BASD examines auditory processing of speech and non-speech stimuli by measuring just noticeable differences (JND) in frequency, level, duration, and amplitude modulation of a presented sinusoidal signal. The non-speech-based part of BASD is based on the Leipzig inventory for assessing centralauditory processing disorders during child development (LIPP: see Ludwig, 2008). Additional test items are speech stimuli for testing the discrimination between different consonants ba/ga (contrast in voicing) and ga/ka (contrast in place of articulation). The test is designed as a playful task, i.e. a memory game, and is especially intended for children from 4 years on. However, the test may also be used for other cooperative subjects. The test is conducted as a three alternatives forced choice procedure. The task is to find two similar items (i.e. with same consonant, frequency, level, duration) from three alternatives. The two similar items are randomly spread across the three memory cards. When pressing the cards images appear as in a memory game. If the two selected cards show the same image the response was correct. This provides an immediate feedback to the tested subject and hence increases motivation.

Start parameters are set as follows: stimulus level differences from 5 to 20 dB (5 dB steps), differences in duration from 100 to 300 ms (50 ms steps), and differences in frequency are set to 0.05, 0.1, 0.2, 0.4, 0.8, and 1.6 octaves. Basic frequencies are 500 Hz und 1 kHz. For the amplitude modulation subtest the reference modulation frequency is 20 Hz. The deviating tone has a modulation frequency of 20 Hz + x Hz with x varying from 0 to 80. Stimulation modes are: monaural (left or right ear), dichotic (signal ipsilateral, broadband noise contralateral), binaural (same signal in both ears), as well as inter-aural (same stimulus on contralateral ear, different stimulus on test ear).

For determining the JND an adaptive algorithm is performed by changing the step size dependent on the duration of the measuring time, i.e. in the beginning steps are larger compared to those at the end of the procedure. For level, frequency and duration tests, in the beginning, the target property (i.e. level difference, frequency difference, duration difference) is increased/decreased after each trial, later it is increased/decreased after two trials (decreased if both trials were correct, increased if at least one trial was incorrect). JNDs are calculated on the basis of the reversal points.

PRACTICAL USE

Select *BASD* from the module selection screen. If more than one speech test is licensed, the test can be found in the *Speech* section. If necessary, change the parameters (e.g. presentation level, presentation mode, frequency for tonal subtests, initial difference for level, frequency, and duration, number of trials and reversals) as required. Select the subtest.

Make sure that a valid transducer is connected (e.g. headphone, insert earphones) and start the test by selecting the test ear. The available test ear configurations depend on the selection of the presentation mode (monaural, binaural, dichotic, interaural) in the module settings. If *monaural*, dichotic or interaural is selected you may test either the right or left ear. If dichotic is selected, all test stimuli are presented together with contralateral noise. If interaural is selected, all test stimuli are

presented with three identical contralateral stimuli (e.g. right ear selected for subtest ba/ga - right ear: ga - ba - ga, left ear: ga - ga - ga). If binaural is selected, both ears are tested together.

Before starting the test, the subject must be instructed about the task. The task is to listen to three stimuli represented by three memory cards. Two of these three stimuli are identical, i.e. they feature the same consonant, level, frequency, duration or amplitude modulation. The two identical stimuli shall be indicated by the subject by uncovering the corresponding two cards.

Depending on age and physical/mental abilities, the test may be conducted completely self-controlled (i.e. subject elicits stimulus and enters response) after instruction or with assistance of the examiner. Supervision by a qualified examiner is recommended at all times.

When starting the test, the main user interface is shown (see *Figure 25*). The stimulus can be elicited by pressing the *play* button ②. The stimulus presentation can be repeated by pressing the *play* button again. The subject responds by pressing two out of three memory cards ④ corresponding to the stimuli that are considered similar regarding the target property (consonant, level, frequency, duration, amplitude modulation). If the selected cards are correct, the test proceeds to the next trial. If the selected cards are not correct, the third card is automatically uncovered in order to provide feedback regarding which combination would have been correct. The blue bars on the left ① and right ③ of the *play* button represent card decks showing the test progress, i.e. the left bar decreases with each trial and the right bar increases with each correct response.




Figure 25: BASD user interface (left: ba/ga differentiation; right: frequency differentiation)

The header provides various information as explained in *Table 7*.

Frequency	Ear	Current differ. / Current position of different item	Test type	Previous differ. / Previous not- pressed item	
500 Hz, 1 kHz for tonal tests	R (right), L (left), B (binaural)	Difference value for tonal tests, otherwise position of different item	ba/ga or ga/ka not available for tonal tests	Difference value for tonal tests, otherwise position of not pressed item	≠ correct (difference could be detected) = incorrect

Table 7: BASD header information

For example, the left screenshot from *Figure 25* provides the following information: test performed on right ear, current item that is different to the other two items is the third item (i.e. item to the right of the screen), the performed subtest is ba/ga, the previously not pressed button is the button on the left of the screen (i.e. the middle and right button were pressed), the response was correct. The right screenshot from *Figure 25* provides the following information: the test is performed with a 1 kHz tone at the right ear, the current frequency difference is 16.6 Hz, the previous frequency difference was 34 Hz, the response was correct.

BASD results are shown as a score for consonant tests and as a discrimination threshold for tonal tests (see *Figure 26*). For tonal tests, the time course of the responses can be displayed by pressing the result screen. Preliminary normative data is available (see *STUDY RESULTS*).



Figure 26: BASD result (left: ba/ga differentiation; middle, right: frequency differentiation result data and time course)

STUDY RESULTS

In a study performed at the University of Leipzig (Freigang *et al.*, 2011), the LIPP test (similar to non-speech-based part of BASD) was performed in 59 adults from 65 to 89 years. Tests were performed at a stimulus level of 35 dB SL with monaural, dichotic, and interaural presentation. The measurements yielded significant age-dependent deteriorations in the ability to discriminate acoustic properties. Most prominently, interaural frequency and duration discrimination at low test frequencies was elevated which could be explained by a deterioration of time- and phase-dependent processing at brain stem and cortical levels.

In a project supported by the German Federal Ministry for Education and Research (BMBF) a speech screening program for pre-school children was developed (Neumann, 2012). The screening program should allow to differentiate between normal speech development (group A), conspicuous speech development with need for educational speech training (group B), and speech development disorders with need for therapy (group C). Moreover, cases of APD within each group should be detected. A previously developed speech screening test for children called KiSS (Euler *et al.*, 2010) was found to exhibit good sensitivity and specificity for identifying conspicuous speech development, but insufficient discrimination between groups B and C. The task of the project was to develop an enhanced psycho-acoustical test battery for examining auditory processing (with differentiation between groups A to C) in children from 4 to 5 years. Moreover, the electrophysiological correlates

of processing the same stimuli with measuring event-related potentials, i.e. mismatch negativity (MMN), was investigated.

In 189 pre-school children from 4 to 4;5 years BASD was tested along several reference tests. Peripheral hearing was tested by means of pure-tone audiometry, tympanometry, stapedius reflex threshold, and Göttinger speech intelligibility test (in quiet) for children. Speech development was tested by means of AWST-R vocabulary test for children, SET-K 3-5 speech development test for children (with subtests morphological rule construction, speech understanding, phonological commemoration of non-words, sentence commemoration), PLAKSS articulation test, and KiSS.2 speech test for children. Auditory processing was tested by means of Göttinger speech intelligibility test (in quiet/in noise) for children, Uttenweiler dichotic test, Mottier test, parental questionnaire. Cognitive abilities were tested by means of Colored Progressive Matrices (CPM), language-free intelligence test, and a phoniatric/audiological examination. Based on these tests, experts classified general speech development and suspicion of APD for each child. From the 189 children, 120 (63.5 %) were classified to group A, 37 (19.6 %) to group B, and 32 (16.9 %) to group C. From 188 children, 13 (6.9 %) were suspected of APD. From the APD-positive children, 2 (1.7 %) were from group A, 4 (10.8 %) from group B, and 7 (22.6 %) from group C. The sensitivity of KiSS.2 for detecting a conspicuous speech development (including group B and C) amounted to 89.7 %, the specificity to 89. 1%. In 71 of the 189 children (with similar shares of each group A to C) MMN was applied. There were no significant differences between the test groups.

Normative BASD data were established in a study with 293 children from 4;0 to 5;11 years. An age-dependent development effect could be observed. Medians of discrimination thresholds are presented in *Table 8* for the different age groups. For frequency differentiation age-dependent dropouts occurred. When considering a frequency discrimination exceeding 500 Hz as invalid, 35 % of all data from age group 4;0 to 4;5 years had to be removed. This share decreased with increasing age to 9 % for age group 5;6 to 5;11 years. This means that the frequency discrimination subtest may not be suited for very young children.

Age (years)	4;0-4;5	5	4;6-4;1	11	5;0-5;5	5	5;6-5;1	11
Level [dB]	8.7	[n=143]	8.6	[77]	7.6	[50]	6.0	[23]
Frequency [Hz]	67.8	[103]	69.0	[59]	62.3	[41]	40.0	[21]
Duration [ms]	146.0	[143]	146.7	[77]	113.4	[50]	65.3	[23]
Ga-Ka [score/12]	6	[122]	5	[73]	6	[50]	10	[23]
Ba-Ga [score/12]	6	[125]	8	[74]	9	[50]	11	[23]

Table 8: BASD discrimination threshold medians for different age groups

Normative data for the 90 % percentile is shown for different age groups in *Table 9*. The correlation between BASD subtests and expert classification was significant for the subtests frequency (if invalid results mainly from young children were removed), place of articulation, and voicing. Also, the differentiation of APD was investigated. Significant correlations could be found between voicing and expert decision as well as between duration and suspicion of parents or of child care workers.

Age (years)	4;0-4;5	4;6-4;11	5;0-5;5	5;6-5;11
Level [dB]	13.4	13.4	12.7	12.7
Frequency [Hz]	296	344	223	70
Duration [ms]	249	254	158	157
Ga-Ka [score/12]	2	2	2	4
Ba-Ga [score/12]	2	3	5	6

Table 9: BASD 90 % percentile for different age groups

LITERATURE

- Biedermann F (2004): Zentrale Verarbeitung akustischer Informationen nach unilateralen Läsionen im auditorischen Di- und Telenzephalon (in German: Central processing of acoustic information after unilateral lesions in the auditory di- and telencephalon), dissertation, University of Leipzig.
- Biedermann F, Bungert P, Dorrscheidt GJ, von Cramon DY, Rübsamen R (2008): Central auditory impairment in unilateral diencephalic and telencephalic lesions, Audiology and Neurootology 13, p. 123-144.
- Bungert-Kahl P, Biedermann F, Dörrscheidt GJ, von Cramon DY, Rübsamen R (2004): Psychoacoustic test tools for the detection of deficits in central auditory processing: Normative data, Audiological Acoustics 43, p. 48-71.
- Euler HA, Holler-Zittlau I, van Minnen S, Sick U, Dux W, Zaretsky Y, Neumann K (2010).
 Psychometrische Gütekriterien eines Kurztests zur Erfassung des Sprachstandes vierjähriger Kinder, (in German: Validity criteria of a short test to assess speech and language competence in 4-year-olds),
 HNO 58, p. 1116–1123.
- Freigang C, Schmidt L, Wagner J, Eckardt R, Steinhagen-Thiessen E, Ernst A, Rübsamen R (2011): Evaluation of central auditory discrimination abilities in older adults, Frontiers in Aging Neuroscience 3(6)
- Ludwig A (2008): Psychoakustische Untersuchungen zu zentral-auditiven Verarbeitungsstörungen während der Kindesentwicklung (in German: Psycho-acoustical examinations regarding central-auditory processing disorders during child development), dissertation, University of Leipzig.
- Neumann K (2012): Ein Screening-Verfahren zur flächendeckenden Erfassung des Sprachstandes vierbis viereinhalbjähriger Kinder: Optimierung, Validierung, Erweiterung, elektrophysiologische Fundierung (in German: A screening procedure for a comprehensive acquisition of the speech status in 4- to 4,5-year-old children: optimization, validation, extension, electro-physiological foundation), final report of BMBF project DLR 01GJ0982.
- Neumann K, Moliadze V, Bieck S, Oswald J, Darquea M, Euler H, Zeratsky Y, Dazert S (2012a): Auditive Verarbeitung und Wahrnehmung vierjähriger Kinder: ein neues Testverfahren mit sprachfreien und sprachhaltigen Stimuli (in German: Auditory processing and perception of 4-year-old children), DGA annual meeting.
- Neumann K, Zaretsky Y, Euler HA, Kunz L, Looschen A, Beste C, Ocklenburg S, Lücke T, Schirkonyer V,
 Oswald J (2012b): A new test to assess central auditory processing of 4-year-old children: A psychoacoustic and MMN study, NHS annual meeting.

3.1.11 Hearing Handicap Inventory for the Elderly (HHIE-S)

METHODOLOGY

In screening questionnaires, subjects are asked to rate their hearing status by answering one or more specific questions. The HHIE questionnaire was developed and standardized by Ventry and Weinstein (1982) and included 25 questions. The abbreviated screening version HHIE-S (Weinstein, 1986) includes 10 questions which concern activity limitations, participation restrictions and associated feelings and emotions of daily life. For example, one of the questions of the emotional category is: Does a hearing problem cause you to feel embarrassed when meeting new people? For example, one of the questions of the participation restriction category is: Do you have difficulty hearing when someone speaks in a whisper? The answer can be Yes, Sometimes, or No. The answers are matched to a multiple-point response scale: Yes = 4 points, Sometimes = 2 points, and No = 0 points. The total score is calculated by adding the points and is then subsequently compared to the defined cut-off to decide whether there is a pass or a fail. Normally, the questionnaire is performed in a face-to-face administration format. The interviewer subsequently reads each question and possible response categories out loud, awaits the subject's response and fills out the corresponding answer on the form.

The questionnaire is especially designed for hearing screening in the elderly, however, may be used for other cooperative subjects as well. Please note that the test is a self-assessment test and therefore depends on the subject's cooperation and ability to responds honestly and realistically and may hence exhibit a result that is deviating from the real hearing status. Further audiologic diagnostics is advised, especially if the HHIE-S screening result suggests hearing problems.

PRACTICAL USE

Select *HHIE-S* from the module selection screen. Select the language in which you would like to perform the questionnaire. Inform the subject about the task: a number of questions from daily life situations are posed (see *Figure 27* left). The subject shall respond with *Yes, Sometimes*, or *No* dependent on which answer seems most appropriate. Dependent on age and physical/mental abilities, the test can be performed completely self-controlled after instruction or with assistance of the examiner (i.e. examiner reads questions and/or enters responses). The result is shown as score with a traffic light status (see *Figure 27* right).




Figure 27: HHIE-S user interface (left) and result (right)

Please note that during this test the footer is hidden. You may activate the footer by pressing the power on/off switch briefly.

STUDY RESULTS

Jupiter and DiStasio (1998) investigated the performance of HHIE-S in 50 elderly subjects from 65 to 85 years with no history of otologic disease, high-level noise exposure or ototoxicity. In the study, the HHIE-S questionnaire was performed together with determination of a three-frequency (0.5, 1, 2 kHz) pure-tone threshold average (Fletcher index) in the better ear. The average HHIE-S score amounted to 9.8 whereas the average pure- tone threshold amounted to 36.8 dB HL. The correlation between the HHIE-S score and the pure-tone threshold average amounted to r = 0.67 for the situational questions, to r = 0.52 for the emotional questions and to r = 0.63 for all questions.

LITERATURE

- Jupiter T, DiStasio D (1998): An evaluation of the HHIE-S as a Screening Tool for the elderly homebound population, Journal of the Academy of Rehabilitative Audiology 31, p. 11-21
- Ventry IM, Weinstein BE (1982): The hearing handicap inventory for the elderly: a new tool, Ear and Hearing 3(3), p. 128-134.
- Weinstein BE (1986): Validity of a screening protocol for identifying elderly people with hearing problems, American Speech-Language-Hearing Association 28, p. 41-45.

3.2 Physiological Test Procedures

Physiological test procedures include otoacoustic emissions (OAE), auditory evoked potentials (AEP), acoustic immittance tests (tympanometry), and other tests that rely on the measurement of physiological properties.

Physiological test procedures do not require active participation of the tested subject during the test. In contrast to psycho-acoustical tests, the ability and willingness of the subject to cooperate only play a subordinate role. Nevertheless, some cooperation is necessary. This may include (depending on the test procedure) that the subject keeps calm, does not swallow, does not move, or does not remove the transducer or any other relevant test equipment from its required position during the test. In contrast to psycho-acoustical test procedures, physiological test procedures may be conducted in subjects that are not able or not willing to perform a psycho-acoustical task, i.e. subjects who are too young or otherwise incapable of responding behaviorally. This may be neonates and small children, mentally or physically handicapped subjects, or people simulating a hearing loss. The examiner must ensure that the subject is capable of fulfilling the respective prerequisites (e.g. being quiet during the test). Also, older subjects should be informed about the test procedure, so that they know what to expect during the test. This may positively contribute to patient relaxation.

Dependent on the test procedure, the test environment must be appropriate. This may include low ambient noise conditions (OAE) and low electromagnetic radiation (AEP). Non-compliance may deteriorate the result.

Physiological test procedures rely on specific physiological processes that are active at different stages of the auditory pathway and hence reflect the mechanical or the neurological functioning of the auditory system. Therefore, physiological tests (single test or combination of tests) can be used for topological hearing diagnostics, i.e. for detecting the site of impairment along the auditory pathway. Tympanometry, OAE and AEP are reportedly able to differentiate between middle-ear, cochlear, and neural disorders. This is important since an adequate therapy for a hearing disorder can only be developed if it is known which stage of the auditory pathway is impaired. With psychoacoustical tests, a differentiation between sound-conductive and sensorineural hearing loss is possible by evaluating the difference between air and bone conduction pure-tone thresholds. However, the discrimination of a sensorineural hearing loss, i.e. the differentiation between a sensory (cochlear) and a neural disorder with subjective testing may be not reliable, because the validity of relevant tests (e.g. Short Increment Sensitivity Index (SISI), Fowler, Carhart) is limited. Moreover, in infants and other non-cooperative patients psycho-acoustical tests cannot reliably be performed.

As for psycho-acoustical tests, proper transducer placement and calibration is important for a reliable outcome. Incorrect placement of a transducer or usage of an incorrectly calibrated transducer may adulterate the result. A suitable transducer should be connected and properly fitted.

If not specified otherwise all physiological tests can be performed in subjects of all ages as long as appropriate participation and cooperation can be expected.

3.2.1 Overview: Otoacoustic Emissions

The discovery of otoacoustic emissions (OAE) (Kemp, 1978) has produced a fast, powerful, and versatile tool for diagnosing cochlear integrity. OAE measurements are today a standard part of the audiometric test battery. OAEs are measured by means of a highly sensitive low-noise microphone within an ear probe that is placed in the outer ear canal. There are spontaneous and evoked OAEs. OAEs are the by-product of the non-linear sound amplification process in the cochlea (Davis, 1983; Dallos, 1992). OAEs are low-level sound emissions generated by the outer hair cells (OHCs) of the inner ear. OAE levels depend on the number of functioning outer hair cells given a normal middle-ear function. Also, OAE levels depend on the ear canal volume. Because of the smaller ear canal volume, OAE amplitude in newborns is usually higher compared to that in adults (Norton, 1992; Lasky, 1998a,b; Abdala, 2000). As a consequence, OAEs in newborns are easier to measure.

Spontaneous otoacoustic emissions (SOAE) are sinusoidal signals and appear without any sound stimulation in a healthy cochlea and seem to be a direct consequence of the cellular force generation of OHCs (Zwicker and Schloth, 1984; Burns *et al.*, 1998; Jülicher *et al.*, 2003). SOAEs do not appear in all normal hearing subjects and are present in about a half of the normal hearing population, with a distinctly higher prevalence in women than in men (Bilger *et al.*, 1990; Penner *et al.*, 1993; Penner and Zhang, 1997). SOAEs are therefore not suited for audiological diagnostics.

Evoked OAEs are triggered by external sounds, either by transient (clicks or tone bursts) or stationary stimuli (sine tones) delivered via loudspeakers within the ear probe. Commonly used types of evoked OAEs are:

Stimulus frequency otoacoustic emissions (SFOAE) are elicited by one continuous, low-level sinusoidal signal. Recording of SFOAEs is difficult because stimulus and response superimpose. SFOAEs are typically derived as the difference in sound pressure in the ear canal with and without a suppressor tone added to the probe tone (Kalluri and Shera, 2001; Neely *et al.*, 2005).

Transient evoked otoacoustic emissions (TEOAE) are elicited by clicks or tone bursts (see *Figure 28*). TEOAEs represent the sum of the pulse responses of OHCs along the cochlea. They already disappear at mild hearing losses and are therefore commonly used in hearing screening programs. For more information on TEOAEs, please refer to section *3.2.2: Transient Evoked Otoacoustic Emissions (TEOAEs)*.

Distortion product otoacoustic emissions (DPOAE) are elicited by two sine tones with a specific frequency and level ratio (see *Figure 28*). DPOAEs represent cubic distortions of OHCs and arise directly from the frequency-selective compressive non-linearity of OHCs. DPOAEs can be applied as a test for a frequency-specific assessment of cochlear dysfunction. DPOAEs are reported to be measurable at a cochlear hearing loss of up to 40 to 50 dB HL, corresponding to the range of amplification of OHCs. For more information on DPOAEs, please refer to section *3.2.3: Distortion Product Otoacoustic Emissions (DPOAEs)*.

Sentiero devices provide both TEOAE and DPOAE test procedures. TEOAEs and DPOAEs are present in essentially every normal-hearing subject. TEOAEs give a rapid overview of cochlear function, whereas DPOAEs provide more quantitative information about sound processing at distinct cochlear sites.

OHCs are reported to be impaired by sound overexposure, ototoxic drugs (e.g. therapeutic antibiotics), infections (e.g. meningitis, mumps, materno-fetal infection), and anoxia (e.g. birth

trauma), or to be partly missing in genetic hearing loss. OHC impairment results in a loss of sensitivity and frequency selectivity of the hearing organ (Liberman and Dodds, 1984). OAEs, as a by-product of cochlear non-linear sound amplification, then appear with reduced amplitude or disappear (Mills and Rubel, 1994).

Since OAEs are a by-product of the non-linear sound amplification process of OHCs in the cochlea they can only serve as a measure for evaluating OHC integrity, i.e. lesions of inner hair cells or retro-cochlear defects (e.g. neural defects, auditory processing disorders) are not detectable by means of OAEs. In sound-conductive hearing loss, both the stimulus and the response amplitude are reduced, so that OAEs are not present, even with a mild sound-conductive hearing loss (Margolis, 2002), as e.g. due to Eustachian tube dysfunction or amniotic fluid in the tympanic cavity. So, if there is no OAE response detectable the following pathologies are possible: sound conductive hearing loss, cochlear hearing loss exceeding 20-30 dB HL (TEOAE) or 40-50 dB HL (DPOAE), or retro-cochlear hearing loss. In these cases, tympanometry, auditory brainstem responses (ABR), and auditory steady-state responses (ASSR) should be performed to determine type and degree of the hearing loss.

As a general rule, if there is a suspicion of a hearing disorder, OAEs should be used first. The test method is fast and helps to confirm normal middle-ear and cochlear function. This is the case if OAEs are present over a wide frequency range. If OAEs are absent, the presence of a middle-ear or cochlear (OHC) pathology is likely. OAEs then should be followed by tympanometry. If the tympanogram is normal and OAEs are absent, then a cochlear disorder is likely. If the tympanogram is abnormal, a sound-conductive hearing loss is likely. If there is an indication for a hearing disorder and both the tympanogram and OAEs are normal, ABR or ASSR may reveal if there is a cochlear (inner hair cell) or neural pathology. For example, in auditory neuropathy, where synchronization of neural activity is malfunctioning (either due to inner hair cell synaptic or neural dysfunction), normal OAEs and abnormal ABRs occur (Doyle *et al.*, 1998; Starr *et al.*, 1996).

Typical clinical applications of OAEs are: hearing screening, follow-up diagnostics after newborn hearing screening, confirmation of cochlear hearing loss (together with tympanometry and ABR), quantitative evaluation of hearing loss and recruitment for providing hearing aid fitting parameters, early detection and monitoring of OHC impairment after noise over-exposure or ototoxic drug administration, topological diagnostics, as well as identifying subjects simulating a hearing loss.

Hearing screening is a selection procedure used to decide whether further diagnostics is advised or not. Consequently, a screening decision is binary, i.e. *pass* (negative finding, no diagnostics necessary) or *refer* (positive finding, follow-up diagnostics advised). The requirements for screening are different from those for diagnostics. Because screening is performed in large populations, the typically used devices provide automatic evaluation rather than rely on an expert's judgement. A screening test should be performed as quickly as possible. Therefore, the respective methods must avoid long preparation times and the test should stop automatically if the desired quality of the result is achieved. OAEs (especially TEOAEs) are widely regarded as being suitable for screening in newborns and infants, as they are quickly measured and not present in the case of OHC dysfunction (e.g. Kemp and Ryan, 1991; Norton *et al.*, 2000a,b). The premise for this approach is that cochlear hearing loss always includes OHC damage or malfunction. It should be emphasized that a sound conductive loss due to Eustachian tube dysfunction and/or amniotic fluid in the tympanic cavity also cause *refer* results under screening conditions mainly due to the attenuation of an existing OAE signal.

A disadvantage of using OAEs in screening protocols is a lower validity as compared to ABR methods (Barker *et al.*, 2000; Norton *et al.*, 2000a). This is especially true in populations with a high prevalence of minor threshold elevation due to a temporary conductive hearing loss, as it is found in term neonates in the first 36 hours of life due to Eustachian tube dysfunction or amniotic fluid in the tympanic cavity or due to a persisting sensory hearing loss in premature and neonatal intensive care unit infants. In order to maintain a high sensitivity, the specificity may be reduced dramatically, making a screening procedure inefficient. To avoid high referral rates, OAE referrals should be followed up with an ABR screening before proceeding with further diagnostics, i.e. two-stage screening (Rhodes *et al.*, 1999; Norton *et al.*, 2000b). DPOAE audiograms may provide an alternative means for revealing a temporary hearing loss. When applying DPOAE audiograms before ABR screening, time and costs can be saved in those babies where DPOAEs are measurable and thus no additional ABR is needed. DPOAE audiograms have an advantage over TEOAEs or click-evoked ABRs because they can quantitatively assess cochlear hearing thresholds at distinct test frequencies.

Especially for hearing-aid adjustment in infants, a quantitative evaluation of the hearing loss and recruitment is necessary. When elicited at high stimulus levels (which is common in clinical practice), TEOAEs are absent at a cochlear hearing loss exceeding 20-30 dB HL, whereas DPOAEs are absent at a cochlear hearing loss exceeding 40-50 dB HL. Thus, when applying both TEOAEs and DPOAEs, a rough estimate of the hearing loss is possible.

OAE measures are stable through time and hence are capable of monitoring recovery from OHC impairment. Therapeutic drugs such as antibiotics (e.g. aminoglycosides) and anti-tumor chemotherapeutic agents (e.g. cisplatin) are reported to induce an irreversible hearing loss, that typically affects the highest frequencies first, with hearing loss systematically progressing to the lower frequencies (e.g. Kopelman et al., 1988; Fausti et al., 1994; Berg et al., 1999; Stavroulaki et al., 2001). Early detection of ototoxicity is important for providing effective management options such as substitution of medications, change of dosage, and mode of administration (Lonsbury-Martin and Martin, 2001). Because TEOAEs are less effective above 4 kHz, DPOAEs are the test of first choice for detecting and monitoring OHC dysfunction due to ototoxic drugs. Moreover, DPOAEs have an additional advantage over TEOAEs, in that they can give information about compression of the OHC amplifiers. If OHC function is disturbed during the toxic process then not only DPOAE level, but also DPOAE growth should be altered. Like antibiotic and chemotherapeutic drugs, salicylate is also known to affect hearing sensitivity and to induce tinnitus (Myers and Bernstein, 1965; McFadden and Plattsmier, 1984; Wier et al., 1988; Long and Tubis, 1988; Boettcher and Salvi, 1991; Brown et al., 1993; McFadden and Pasanen, 1994). However, most importantly, impairment due to salicylate toxicity is reversible. Assuming that a loss of OHC stiffness is responsible for distortions within cochlear micromechanics, the corresponding change in inner hair cell activity may be one potential correlate of tinnitus (Janssen et al., 2000).

Microphone noise, physiological noise (breathing, blood flow), and ambient noise do not allow proper OAE measurements at very low stimulus levels. At low frequencies, OAE measurements are not reliable even at high stimulus levels mainly due to technical noise. Because of the limited frequency range of the ear probe's electro-acoustic transducers, OAE measurements are also limited towards higher frequencies with increasing variability of OAE levels (Dreisbach *et al.*, 2006). To achieve low noise floor levels, OAE measurements have to be done in a sound-attenuating booth or any other quiet environment especially when close-to-threshold stimulus levels are used. Alternatively, a sound insulation headphone, that is covering the ear probe, may be used. Moreover,

the subject needs to be calm, i.e. not heavily moving, breathing, or swallowing. The ear probe cable needs to hang loosely, so that it does not touch any material as e.g. clothes. Touching or rubbing the cable may add noise.

Automated measuring and evaluation procedures guarantee test consistency and simplify the interpretation of OAE recordings. Fourier transformation computations from the time domain signal allow automatic evaluation of OAE signals. To minimize the influence of unwanted external signals, algorithms for noise reduction and artefact rejection are applied. The noise floor level is usually higher at low frequencies due to microphone properties and low-frequency body sounds such as breathing. Artefact rejection can be performed by elimination of high noise level buffers or by weighting each buffer dependent on its noise content. In addition, the noise floor level is reduced by time domain averaging of the recorded signal. The idea of averaging is that the signal is constant and hence the same in each buffer whereas the noise is random and hence changes in each buffer. Adding all the buffers increases the signal and reduces the noise. Theoretically, the improvement in SNR is proportional to the square root of the number of samples that are averaged, limited by technical properties (e.g. quantization noise).

OAEs can only be tested with a special ear probe that typically contains one (for TEOAE) or two (for TEOAE or DPOAE) loudspeakers and a microphone (see *Figure 28*). For DPOAE recording, separate loudspeakers are commonly used for each primary tone in order to exclude technically generated distortion components.




Figure 28: Schematic overview of TEOAE and DPOAE measurement (A/D = Analog-Digital-Converter, D/A = Digital-Analog-Converter, PC = Personal Computer, L = sound pressure level, p = sound pressure, f = frequency, t = time)

Before each OAE measurement, a calibration of the ear probe is automatically performed in order to adapt the stimulus output to the ear canal volume. However, due to standing wave effects, the

estimated ear canal volume may not exactly reflect the real ear canal volume in all cases (Siegel, 1994; Whitehead *et al.*, 1995). The ear probe needs to be inserted properly with a tight fit and without any leakage between ear probe and ear canal. If there is leakage, low frequency sound components cannot be delivered to the cochlea properly and hence no apical and medial OHCs will contribute to the OAE response. Furthermore, proper occlusion of the ear canal by the ear tip diminishes the influence of external sounds. Ear probe channels must not be blocked by e.g. cerumen.

For a review of methodology, technical and clinical aspects of OAEs see Janssen and Müller (2008).

PRACTICAL USE

In general, it is recommended to check the functionality of the ear probe regularly (e.g. if the ear probe is frequently used: once a day before starting to test patients or if any failure of the ear probe is suspected). The ear probe can be tested by placing the ear probe together with an adequate probe tip in the correct test cavity (small probe tip PT-S fits into blue test cavity, large probe tip PT-A fits into red test cavity). Start the probe test (see *Figure 29*) and check if the probe test passes. If not check the ear probe cable, the probe tip placement on the ear probe, and the ear probe and probe tip channels e.g. regarding contamination). If applicable, clean the probe tip. Do not use sharp items with the ear probe or probe test does not pass or if you suspect any dysfunction, please retry with another ear probe or probe tip and/or contact your distributor.




Figure 29: Probe test

LITERATURE

- Abdala C (2000): Distortion product otoacoustic emission (2f1-f2) amplitude growth in human adults and neonates, Journal of the Acoustical Society of America 107, p. 446-456.
- Barker SE, Lesperance MM, Kileny PR (2000): Outcome of newborn hearing screening by ABR compared with four different DPOAE pass criteria, American Journal of Audiology 9, p. 142-148.
- Berg AL, Spitzer JB, Garvin JH Jr. (1999): Ototoxic impact of cisplatin in pediatric oncology patients, Laryngoscope 109, p. 1806-1814.
- Bilger RC, Matthies ML, Hammel DR, Demorest ME (1990): Genetic implications of gender differences in the prevalence of spontaneous otoacoustic emissions, Journal of Speech and Hearing Research 33, p. 418-432.
- Boettcher FA, Salvi RJ (1991): Salicylate ototoxicity: review and synthesis, American Journal of Otolaryngology 12, p. 33-47.

- Brown AM, Williams DM, Gaskill SA (1993): The effect of aspirin on cochlear mechanical tuning, Journal of the Acoustical Society of America 93, p. 3298-3307.
- Burns RM, Keefe DH, Ling R (1998): Energy reflectance in the ear canal can exceed unity near spontaneous otoacoustic emission frequencies. Journal of the Acoustical Society of America 103, p. 462-474.
- Dallos P (1992): The active cochlea, Journal of Neuroscience 12, p. 4575-4585.
- Davis H (1983): An active process in cochlear mechanics. Hearing Research 9, p. 79-90.
- Doyle KJ, Fujikawa S, Rogers P, Newman E (1998): Comparison of newborn hearing screening by transient otoacoustic emissions and auditory brainstem response using ALGO-2, International Journal of Pediatric Otorhinolaryngology 43, p. 207-211.
- Dreisbach LE, Long KM, Lees SE (2006): Repeatability of high-frequency distortion-product otoacoustic emissions in normal-hearing adults. Ear and Hearing 27(5), p. 466-479
- Fausti SA, Larson VD, Noffsinger D, Wilson RH, Phillips DS, Fowler CG (1994): High-frequency audiometric monitoring strategies for early detection of ototoxicity, Ear and Hearing 15, p. 232-239.
- Janssen T, Müller J (2008): Active processes and otoacoustic emissions, In: Springer handbook of auditory research (editors: Manley AG, Fay RR, Popper AN), Springer, p. 421-460.
- Janssen T, Boege P, Oestreicher E, Arnold W (2000): Tinnitus and 2f1-f2 distortion product otoacoustic emissions following salicylate overdose, Journal of the Acoustical Society of America 107, p. 1790-1792.
- Jülicher F, Camalet S, Prost J, Duke TAJ (2003): Active amplification by critical oscillations, In: Biophysics of the cochlea: from molecule to models (editor: Gummer A), World Scientific.
- Kalluri R, Shera CA (2001): Distortion-product source unmixing: A test of the two-mechanism model for DPOAE generation, Journal of the Acoustical Society of America 109, p. 622-637.
- Kemp DT (1978): Stimulated acoustic emissions from within the human auditory system. Journal of the Acoustical Society of America 64, p. 1386-1391.
- Kemp DT, Ryan S (1991): Otoacoustic emission tests in neonatal screening programmes. Acta Otolaryngologica Supplementum 482, p. 73-84.
- Kopelman J, Budnick AS, Sessions RB, Kramer MB, Wong GY (1988): Ototoxicity of high-dose cisplatin by bolus administration in patients with advanced cancers and normal hearing, Laryngoscope 98, p. 858-864.
- Lasky RE (1998a) Distortion product otoacoustic emissions in human newborns and adults. I. Frequency effects, Journal of the Acoustical Society of America 103, p. 981-991.
- Lasky RE (1998b) Distortion product otoacoustic emissions in human newborns and adults. II. Level effects, Journal of the Acoustical Society of America 103, 992-1000.
- Liberman MC, Dodds LW (1984): Single-neuron labeling and chronic cochlear pathology. III. Stereocilia damage and alterations of threshold tuning curves. Hearing Research 16, p. 55-74.
- Long GR, Tubis A (1988): Modification of spontaneous and evoked otoacoustic emissions and associated psychoacoustic microstructure by aspirin consumption, Journal of the Acoustical Society of America 84, p. 1343-1353.
- Lonsbury-Martin BL, Martin GK (2001): Evoked otoacoustic emissions as objective screeners for ototoxicity, Seminars in Hearing 22, p. 377-392.
- Margolis RH (2002): Influence of middle-ear disease on otoacoustic emissions, In: Otoacoustic emissions: clinical applications (editors: Robinette MS, Glattke TJ), Thieme.
- McFadden D, Pasanen EG (1994): Otoacoustic emissions and quinine sulphate, Journal of the Acoustical Society of America 95, p. 3460-3474.

- McFadden D, Plattsmier HS (1984): Aspirin abolishes spontaneous oto-acoustic emissions, Journal of the Acoustical Society of America 76, p. 443-448.
- Mills DM, Rubel EW (1994): Variation of distortion product otoacoustic emissions with furosemide injection, Hearing Research 77, p. 183-199.
- Myers EN, Bernstein JM (1965): Salicylate ototoxicity; a clinical and experimental study, Archives of Otolaryngology 82, p. 483-493.
- Norton SJ (1992): The effects of being a newborn on otoacoustic emissions, Journal of the Acoustical Society of America 91, p. 2409-2409.
- Norton SJ, Gorga MP, Widen JE, Folsom RC, Sininger Y, Cone-Wesson B, Vohr BR, Mascher K, Fletcher K (2000a): Identification of neonatal hearing impairment: evaluation of transient evoked otoacoustic emission, distortion product otoacoustic emission, and auditory brain stem response test performance, Ear and Hearing 21, p. 508-528.
- Norton SJ, Gorga MP, Widen JE, Folsom RC, Sininger Y, Cone-Wesson B, Vohr BR, Fletcher KA (2000b): Identification of neonatal hearing impairment: summary and recommendations, Ear and Hearing 21, p. 529-535.
- Penner MJ, Glotzbach L, Hunag T (1993): Spontaneous otoacoustic emissions: measurement and data, Hearing Research 68, p. 229-237.
- Penner MJ, Zhang T (1997): Prevalence of spontaneous otoacoustic emissions in adults revisited, Hearing Research 103, p. 28-34.
- Rhodes MC, Margolis RH, Hirsch JE, Napp AP (1999): Hearing screening in the newborn intensive care nursery: comparison of methods, Otolaryngology Head Neck Surgery 120, p. 799-808.
- Siegel JH (1994): Ear-canal standing waves and high-frequency sound calibration using otoacoustic emission probes, Journal of the Acoustical Society of America 95, p. 2589-2597.
- Starr A, Picton TW, Sininger Y, Hood LJ, Berlin CI (1996): Auditory neuropathy, Brain 119(3), p. 741-753.
- Stavroulaki P, Apostolopoulos N, Segas J, Tsakanikos M, Adamopoulos G (2001): Evoked otoacoustic emissions an approach for monitoring cisplatin induced ototoxicity in children, International Journal of Pediatric Otorhinolaryngology 59, p. 47-57.
- Whitehead ML, Stagner BB, Lonsbury-Martin BL, Martin GK (1995): Effects of ear-canal standing waves on measurements of distortion product otoacoustic emissions, Journal of the Acoustical Society of America 98, p. 3200-3214.
- Wier CC, Pasnanen EG, McFadden D (1988): Partial dissociation of spontaneous otoacoustic emissions and distortion products during aspirin use in humans, Journal of the Acoustical Society of America 84, p. 230-237.
- Zwicker E, Schloth E (1984): Interrelation of different oto-acoustic emissions, Journal of the Acoustical Society of America 75, p. 1148-1154.

3.2.2 Transient Evoked Otoacoustic Emissions (TEOAEs)

METHODOLOGY

Transient otoacoustic emissions (TEOAEs) elicited by clicks represent the sum of the pulse responses of OHCs along the cochlea (see *Figure 28*). Almost all OHCs along the cochlear partition or a part of them (the site in the cochlea depending on the carrier frequency of the tone burst) are excited. Due to frequency dispersion in the cochlea, a specific component of the TEOAE response can be directly traced to a specific frequency component of the transient signal. As the basilar membrane at basal sites moves faster than at more apical sites, high-frequency TEOAE components stem from basal cochlear sites, whereas low-frequency TEOAE components come from more apical ones. Also, basal responses appear at the beginning and apical responses at the end of the TEOAE time function.

TEOAE responses are typically evoked by one of the following stimulus trains: (i) four clicks of equal magnitude (linear protocol), (ii) three clicks of positive polarity followed by a fourth click of inverse polarity with a relative magnitude of 9.5 dB higher than the corresponding positive clicks (non-linear protocol) (Kemp *et al.*, 1986; Bray, 1989). Under the hypothesis that the TEOAE recordings originate from saturated cochlear generators, it is assumed that the non-linear protocol removes stimulus artefacts of a linear nature (i.e. the stimulus itself), because stimulus signals increase linearly with the stimulus level, while response signals (i.e. emission emerging from the non-linear operation of OHCs) increase non-linearly with the stimulus level. It is generally accepted that the non-linear protocol is a practical compromise to maximize the reliability of a TEOAE recording (Kemp *et al.*, 1990a,b; Grandori and Ravazzani, 1993; von Specht *et al.*, 2001; Hatzopoulos *et al.*, 2003).

TEOAE signals are recorded between the short pauses between the stimuli within the stimulus trains. There are several objective methods for separating the TEOAE signal from the background noise and for automatically evaluating the validity of a recorded emission. The first method is based on the calculation of the buffer correlation of the time domain averaged signals between two separate signal buffers (Kemp *et al.*, 1990a). If the two buffers are completely identical, the correlation coefficient is 1 and thus the reproducibility 100 %. A signal is commonly accepted as valid for a reproducibility exceeding a minimum of 60 %. The second method relies on the computation of the spectral power ratio of the sum and the difference of the two signal buffers, denoted as the signal-to-noise ratio (SNR). The pass criterion for a valid signal is typically set to an SNR of 6 dB. The third signal validation procedure is based on a binomial statistical test, which determines the statistical probability that an emission has been recorded. Binomial statistics reduces the recorded signal to binary events, and uses knowledge on the expected distribution of these events (binomial distribution) (Giebel, 2001).

The inter-individual variance of the TEOAE level is high with a standard deviation exceeding 10 dB (Kemp *et al.*, 1986; Probst *et al.*, 1987; Bonfils and Uziel, 1989; Smurzynski and Kim, 1992). However, the intra-individual variance is quite low with a standard deviation around 1 dB (Harris *et al.*, 1991).

When recording TEOAEs, stimulus artefacts may generate signals being in phase in two averaging buffers resulting in a pseudo response of high reproducibility. By means of windowing functions, the stimulus artefact can be excluded so that the reproducibility of the overall signal is restricted to the signal section of interest (Kemp *et al.*, 1990a, Kemp *et al.*, 1990b). It should be emphasized that since the stimulus artefact always appears in the early recording period, the high frequency TEOAE components from the basal site of the cochlea get lost as a result of the windowing procedure. Due

to the fact that the stimulus and the high-frequency TEOAE components superimpose and therefore have to be cancelled during TEOAE recording, TEOAEs fail to measure OHC functionality in the basal region of the cochlea, i.e. above 4 kHz. Therefore, a high-frequency hearing loss cannot be detected by means of TEOAEs. Also, valid TEOAE responses can be present in ears with a mid-frequency hearing loss (e.g. congenital hearing loss). The problem with TEOAEs here is that infants with congenital hearing loss may be overlooked and the hearing loss is detected later in life.

The wide-band TEOAE stimulus calibration is not influenced as much as DPOAEs by standing wave problems. Moreover, TEOAEs are usually stimulated with a relatively high click level, i.e. where cochlear compression already saturates basilar membrane displacement. As opposed to DPOAEs, no level ratio between primary tones needs to be fulfilled. This results in TEOAEs being less susceptible to stimulus calibration errors.

TEOAEs already disappear at low conductive hearing losses (> 10-15 dB HL) or mild cochlear hearing losses (> 20-30 dB HL) and are therefore a suitable tool for (newborn) hearing screening (Robinette and Glattke, 2002; Janssen, 2009; Janssen and Müller (2009)).

In order to get more frequency-specific and quantitative information on hearing loss, DPOAE or ASSR threshold measurements are advised.

Two TEOAE workflows are available:

- TEOAE Quick allows for simple TEOAE testing with an automated statistical algorithm for response detection with fixed detection criteria. TEOAE Quick assesses cochlear outer hair cell function more qualitatively. The test can be performed binaurally if two ear probes are connected. In doing this, test time is reduced by a factor of two.
- **TEOAE Diagnostic** allows for a more profound testing for assessing responses in five half-octave frequency bands around 1, 1.5, 2, 3, and 4 kHz and at different stimulus levels. User-definable criteria can be set, i.e. SNR criteria (6 or 9 dB) for each frequency band and an overall pass criterion (number of passed frequency bands for overall pass: 3/5, 4/5, 5/5).

PRACTICAL USE

Select the TEOAE test that you would like to perform from the module selection screen, i.e. *TEOAE Quick* or *TEOAE Diagnostic*. If more than one OAE test is licensed, *TEOAE Quick* and *TEOAE Diagnostic* can be found in the *OAE* section. For TEOAE Diagnostic, if necessary, change the parameters (e.g. single band SNR criterion, overall pass criterion) as required. Please note that the selection of SNR and pass criteria influence the overall test result, i.e. different settings may result in different overall results. The stricter the selected criteria, the more prominent the TEOAE response must be yielding a valid result. The criteria are saved together with the result data and cannot be changed afterwards.

The subject should be calm and sitting comfortably in a chair or lying on a bed. For babies, try to test the subject during sleep. Make sure that a valid ear probe (e.g. EP-TE or EP-DP) is connected and that the cable of the ear probe does not rub against any material as e.g. clothes. Select an ear tip with appropriate size matching the probe tip size and the subject's ear canal size. Make sure that the ear probe is inserted without any leakage between ear probe and ear canal. In cases where the ear probe is not inserted properly or ear probe channels are blocked a warning will appear.

Select the test ear (Right+Left: simultaneous measurement of right and left ear – only available when two ear probes are connected). Before the measurement starts, the calibration screen is shown (see Figure~30). In expert calibration view mode (configuration in TEOAE~Preferences), during calibration the parameters stimulus level, symmetry between both channels (not for one channel ear probe EPTE), and leakage are checked (green check mark: ok, gray \emptyset : not ok) ①. If the stimulus and symmetry criteria pass but the leak check does not pass, after a while the calibration can be skipped with the forward button ②. If all three criteria are met, the calibration graph turns green. In simple calibration view mode (configuration in TEOAE~Preferences), during calibration a traffic light status is shown. As long as not all calibration criteria are met, the traffic light remains red or yellow. If all criteria are met, the traffic light turns green. After successful ear probe calibration, the measurement starts automatically.



Figure 30: TEOAE ear probe calibration (top left: expert mode monaural; top middle: expert mode monaural with skip option; top right: expert mode binaural; bottom left/middle: simple mode monaural; bottom right: simple mode binaural)

A non-linear broadband click stimulus sequence is presented for TEOAE Quick at a fixed level (85 dB peSPL) and for TEOAE Diagnostic at a user-defined level. The response is detected via the ear probe microphone.

During the measurement (see *Figure 31*) the overall measurement progress (5) and the noise floor (4) are displayed. The noise floor is an indicator for the adequacy of measurement conditions. With good measurement conditions, the noise floor bar should remain in the lower quarter. If the noise floor bar is higher, ambient noise levels may be too high, the patient may be not calm enough (e.g.

swallowing, moving), or the ear probe cable may rub e.g. on clothes. If the leak check is enabled in *TEOAE Preferences* and the criterion fails during the test (e.g. when the probe slips out of the ear canal) at the place of the green noise bar the warning information "unstable" is shown. The overall measurement time depends on test conditions, i.e. the higher the noise floor the lower the progress bar speed.

The time-domain response (averaged into two buffers, time window from 5 to 13 ms) ② and the frequency spectrum (signal: red area (right ear), blue area (left ear), noise: gray area) are shown during the measurement. For *TEOAE Quick* a continuous spectrum ① and a TEOAE validity bar ③ are shown. For *TEOAE Diagnostic*, the spectrum is shown separated for the different frequency bands together with the current SNR values ⑧. A cartoon mode (plane or other object flying through a comic style landscape ⑨) can be activated by pressing the *plane* button ⑦. The cartoon mode is especially meant for focusing a child's interest and hence improving noise conditions during a test. The test can be aborted by the examiner with the *stop* button ⑥.




Figure 31: TEOAE measurement (top: TEOAE Quick; bottom: TEOAE Diagnostic; left: monaural; right: binaural; bottom far right: cartoon mode)

After the test is finished, the result screen (see *Figure 32*) shows the overall result status ① (green symbol if the test criterion is reached or red symbol if not; termination of the measurement by the user is indicated by a yellow question mark symbol). For further analysis, the spectrum of the response signal and noise (TEOAE Quick: continuous spectrum ②; TEOAE Diagnostic: frequency bands ⑦) and the time signal ③ are shown together with quality data, i.e. noise level ⑤, artefact rate, and stimulus stability ⑥. For TEOAE Quick the statistical validity can be read from the TEOAE

validity bar ④. For TEOAE Diagnostic additional result views can be displayed by pressing the result screen, i.e. SNR and overall pass criteria ⑧, continuous spectrum ②, and test parameters (stimulus type, level, and averaged) ⑨.

In case no valid result is detected and the artefact rate exceeds about 20 % or the stimulus stability is below about 80 %, please try to eliminate possible causes (e.g. ambient noise, inappropriate ear probe placement) and restart the measurement.




Figure 32: TEOAE result (top left: TEOAE Quick; top right: TEOAE Diagnostic frequency band view; bottom left: SNR view, bottom middle: spectrum view, bottom right: parameter view)

LITERATURE

- Bonfils P, Uziel A (1989): Clinical applications of evoked acoustic emissions: results in normally hearing and hearing-impaired subjects, Annals of Otology, Rhinology, and Laryngology 98, p. 326-31.
- Bray P (1989): Click evoked otoacoustic emissions and the development of a clinical otoacoustic hearing test instrument, dissertation, University College and Middlesex School of Medicine, London.
- Giebel A (2001): Applying signal statistical analysis to TEOAE measurements, Scandinavian Audiology Supplementum 52, p. 130-132.
- Grandori F, Ravazzani P (1993): Non-linearities of click-evoked otoacoustic emissions and the derived non-linear technique, British Journal of Audiology 27, p. 97-102.
- Harris FP, Probst R, Wenger R (1991): Repeatability of transiently evoked otoacoustic emissions in normally hearing humans, Audiology 30, p. 135-141.

- Hatzopoulos S, Petrucelli J, Morlet T, Martini A (2003): TEOAE recording protocols revised: data from adult subjects, International Journal of Audiology 42, p. 339-347.
- Janssen T (2009): Otoakustische Emissionen (in German: Otoacousic emissions), In: Praxis der Audiometrie (editors: Lehnhardt E, Laszig R), Thieme, p. 113-136.
- Janssen T, Müller J (2008): Otoacoustic emissions as a diagnostic tool in a clinical context. In: Manley GA, Fay RR, Popper AN. Springer Handbook of Auditory Research (Vol 30). Heidelberg, Germany
- Kemp DT, Bray P, Alexander L, Brown AM (1986): Acoustic emission cochleography practical aspects, Scandinavian Audiology Supplementum 25, p. 71-95.
- Kemp DT, Ryan S, Bray P (1990a): A guide to the effective use of otoacoustic emissions, Ear and Hearing 11, p. 93-105.
- Kemp DT, Ryan S, Bray P (1990b): Otoacoustic emission analysis and interpretation for clinical purposes, In: Cochlear mechanisms and otoacoustic emissions (editors: Grandori F, Ciafrone G, Kemp D), Karger.
- Probst R, Lonsbury-Martin BL, Martin GK, Coats AC (1987): Otoacoustic emissions in ears with hearing loss, American Journal of Otolaryngology 8, p. 73-81.
- Robinette MS, Glattke TJ (2002): Transient evoked otoacoustic emissions, In: Otoacoustic emissions: clinical applications (editors: Robinette MS, Glattke TJ), Thieme.
- Smurzynski J, Kim DO (1992): Distortion-product and click-evoked otoacoustic emissions of normally-hearing adults, Hearing Research 58, p. 227-240.
- von Specht H, Ganz M, Pethe J, Leuschner S, Pytel J (2001): Linear versus non-linear recordings of transiently-evoked otoacoustic emissions – methodological considerations, Scandinavian Audiology 30, p. 116-118.

3.2.3 Distortion Product Otoacoustic Emissions (DPOAEs)

METHODOLOGY

Distortion product otoacoustic emissions (DPOAEs) represent cubic distortions of outer hair cells (OHC) when stimulated simultaneously by two tones f_1 (lower frequency) and f_2 (higher frequency) (see Figure 28). DPOAEs arise directly from the frequency-selective compressive nonlinearity of OHCs (Brownell et al., 1985; Kemp et al., 1986). The two primary tones interact in the cochlea within the region of overlap of the traveling waves of the two primary tones close to the characteristic place of f_2 . Thus, DPOAEs can be applied for frequency-specific assessment of cochlear dysfunction at the f_2 place. In humans, both quadratic (e.g. f_2 - f_1) and cubic distortion products (e.g. $2f_1$ - f_2) can be detected. The cubic distortion component $2f_1$ - f_2 yields the highest amplitude and is therefore primarily used in audiological diagnostics (Gorga et al., 2000). DPOAE amplitudes typically range from about 20 dB SPL down to the limiting noise floor level, i.e. about -20 dB SPL. DPOAEs provide quantitative and frequency-specific information about the range and operational characteristics of the cochlear amplifier, i.e. sensitivity, compression, and frequency selectivity of the hearing organ. Especially, extrapolated DPOAE I/O functions (see Figure 33) allow for assessing loss of cochlear sensitivity and compression (Janssen and Müller, 2008). The number of OHCs contributing to DPOAE generation depends on the size of the overlapping region, which is determined by the primary tone levels L_1 and L_2 , and the frequency ratio f_2/f_1 . A frequency ratio of about 1.2 was found to be optimal. Different primary tone level ratios have been suggested. A primary tone level setting, which accounts for the different compression of the primary tone traveling waves at the f_2 place, is the scissor paradigm (Whitehead et al., 1995; Kummer et al., 2000; Boege and Janssen, 2002). Due to the steep slope of the traveling wave towards the cochlear apex, the maximum interaction site is close to the f_2 place in the cochlea. To preserve optimum overlap of the primary tone traveling waves at a constant frequency ratio, the primary tone level difference has to be increased with decreasing stimulus level. This results in a decrease in L_1 being lower than the decrease in L_2 (see Figure 33).




Figure 33: Schematic overview of DPOAE generation with different primary tone level settings: scissor paradigm (solid lines), $L_1 = L_2$ (dotted lines)

Despite the fact that DPOAEs are supposed to primarily reflect OHC activity at the f_2 place, there is evidence that DPOAEs are generated by two distinct cochlear sources (Whitehead *et al.*, 1992; Brown *et al.*, 1996; Shera and Guinan, 1999). As already mentioned, the first source, the effect of which is actually intended to be measured, is located at the region of overlap of the traveling waves of the two primary-tones near the f_2 place and is due to intermodulation distortion. The second source, which is unintentionally adding constructively or destructively to the first source emission, is located at the characteristic frequency place of the emission at $2f_1$ - f_2 and is due to reflection of energy that has travelled apically from the overlap region near f_2 . Thus, energy from both interacting sources yields the composite DPOAE signal which is actually recorded in the outer ear canal. The influence of the second DPOAE source may be observed when monitoring the DPOAE level across frequency with narrow frequency spacing of f_2 (DPOAE fine-structure).

There are several objective methods for separating the DPOAE signal from the background noise and for automatically evaluating the validity of a recorded emission. The first method is based on the calculation of the noise floor level by averaging the levels of several adjoining frequency components around the DPOAE frequency component, with the SNR being indicated by the difference between the emission level and the noise floor level. The SNR criterion is usually set to 6 dB. The second data validation procedure is based on phase statistics, which checks the coupling of the DPOAE component phase to the phase of the primary tones. The phase statistics average normalized phase vectors of the signal received at the known DPOAE frequency. Like the binomial statistics, the vector sum can be scaled in probability terms, providing defined and very high sensitivities. A typical level of significance exceeds 99 % per single frequency test.

The intra-individual variance of the DPOAE level is quite low with a standard deviation below 2 dB (Johnsen and Elberling, 1982a,b). Repetitive DPOAE measurements with unchanged sound probe position exhibited an exponentially increasing standard deviation of DPOAE level with increasing SNR (Janssen *et al.*, 2005a). For example, at an SNR of 10 dB, the standard deviation amounts to 1.8 dB, at an SNR of 20 dB to 0.7 dB, and at an SNR of 40 dB to 0.1 dB. This means that the higher the SNR, the higher the reliability of the DPOAE measurement. This finding is important with respect to the evaluation of small DPOAE changes. For clinical practice, repetitive DPOAE measurements with changed sound probe position are relevant. The standard deviation of DPOAE level amounted to about 1.6 dB (Müller *et al.*, 2005) when changing sound probe position between measurements.

DPOAE grams plot the DPOAE level $L_{\rm dp}$ as a function of f_2 (the main DPOAE generation site) for a selected combination of primary-tone levels L_1 and L_2 . It should be emphasized that DPOAE grams reflect the sensitivity of the cochlear amplifier (CA) best when recorded at close-to-threshold stimulus levels (Janssen *et al.*, 1998; Kummer *et al.*, 1998; Dorn *et al.*, 2001). In normal hearing (normal CA), DPOAE grams are close to each other at high and more separated at low stimulus levels reflecting cochlear non-linear sound processing. In cochlear hearing loss ears (impaired CA), DPOAE grams are more separated even at high stimulus levels, revealing loss of CA compression (Janssen *et al.*, 1998, Kummer *et al.*, 1998, Neely *et al.*, 2003).

DPOAE level I/O-functions plot the DPOAE level $L_{\rm dp}$ as a function of primary-tone level L_2 for a selected f_2 and thus reflect CA dynamics at the f_2 place in the cochlea (Dorn *et al.*, 2001). In normal hearing, in response to low-level stimuli, DPOAE level I/O-functions exhibit steep slopes, while at high stimulus levels slopes decrease, thus mirroring the strong amplification at low and decreasing amplification (saturation) at moderate sound levels. However, this is only true when a specific

stimulus level setting is used which accounts for the different compression of the primary-tones at the f_2 place (scissor paradigm: Kummer *et al.*, 2000).

DPOAE pressure I/O-functions plot the DPOAE pressure p_{dp} (instead of the DPOAE level L_{dp}) as a function of the primary-tone level L_2 . Due to the logarithmic dependency of the DPOAE level on the primary tone level there is a linear dependency between DPOAE pressure p_{dp} and primary tone level L_2 (Boege and Janssen, 2002). Thus, DPOAE data can easily be fitted by linear regression analysis.

The intersection point of the linear regression line with the L_2 -axis at $p_{dp} = 0$ Pa can then serve as an estimate of the stimulus level at the DPOAE threshold (see Figure 34; Boege and Janssen, 2002; Gorga et al., 2003). The estimated threshold level L_{dpth} when plotted across frequency f_2 provides a measure for estimating CA threshold at the f_2 place. Due to standing waves in the outer ear canal, stimulus and response cannot always be reliably determined and thus DPOAE thresholds do not always match pure-tone thresholds. This is true especially for large ear canal volumes. If this is the case a discrepancy between DPOAE thresholds and PTA thresholds mainly occurs for adults (dependent on ear canal size) in the mid-frequency region around 3 kHz and at the higher test frequencies (> 6 kHz). Standing wave effects are less important in newborns and infants due to the shorter ear canal length (Keefe et al., 1993). The relation between OAE level and auditory threshold – or rather the lack of it - is strongly debated. Formerly, it was common to define confidence limits to determine the degree of certainty with which any measured response could be assigned to either normal or impaired hearing (Gorga et al., 1996; Gorga et al., 2000), or to define a DPOAE detection threshold as the stimulus level at which the response equaled the noise present in the instrument (Dorn et al., 2001). However, since the noise is of technical origin (e.g. microphone noise) the threshold evaluated in this way does not match the behavioral threshold. A more relevant measure is the intersection point between the extrapolated DPOAE I/O-function and the primary tone level axis at which the response's sound pressure is zero and hence at which OHCs are inactive. A linear dependency between DPOAE sound pressure and primary tone sound pressure level is present when using the scissor paradigm for eliciting DPOAEs. The estimated DPOAE threshold $L_{\rm doth}$ is independent of noise and seems to be more closely related to behavioral threshold than the DPOAE detection threshold (Boege and Janssen, 2002; Gorga et al., 2003; Janssen et al., 2006).




Figure 34: Schematic overview of a DPOAE audiogram derived from DPOAE threshold estimation

When plotting the DPOAE threshold in dB hearing level (HL), the estimated DPOAE thresholds can be plotted in an audiogram form (**DPOAE audiogram**) (see *Figure 34*). DPOAE audiograms can be applied in babies with a refer result in newborn hearing screening to reveal a transitory conductive hearing loss due to Eustachian tube dysfunction and/or amniotic fluid in the tympanic cavity or to confirm a

persisting cochlear hearing loss in follow-up diagnostics. In case of mild and moderate hearing loss DPOAE audiograms are an alternative method to behavioral audiometry or frequency-specific evoked response audiometry (auditory brainstem responses (ABR) with narrowband stimuli, auditory steady state responses (ASSR)). Especially in infants where the conditioned free-field audiogram does not reflect the real threshold, DPOAE audiograms may assess cochlear hearing loss more precisely than behavioral tests. Moreover, unilateral hearing loss can be detected. DPOAE audiograms are able to quantitatively assess the hearing loss at distinct frequencies in a couple of minutes. Predicting hearing loss at five frequencies by tone burst ABR or ASSR may take half an hour and more. This is an essential advantage of DPOAE over tone burst ABR or ASSR. Thus, DPOAE audiograms can serve as an advanced tool for bridging the gap between screening and audiological testing in pediatric audiology.

DPOAE slope, calculated from DPOAE level I/O-functions (e.g., between stimulus levels L_2 of 40 and 60 dB SPL), indicates CA compression. When plotted across frequency, a slope profile can be established. In ears with cochlear hearing loss, the slope s of the DPOAE level I/O-function increases with increasing hearing loss indicating loss of CA compression (Janssen $et\ al.$, 1998; Kummer $et\ al.$, 1998; Müller and Janssen, 2004; Neely $et\ al.$, 2003). DPOAE slope differs significantly between hearing loss classes (Janssen $et\ al.$, 2005b). DPOAE slope is reported to be related to the slope of the loudness functions (Neely $et\ al.$, 2003, Müller and Janssen, 2004). Thus, the slope of DPOAE I/O-functions is suggested to allow a quantitative assessment of CA compression and hence provide an objective recruitment test. Especially for hearing aid adjustment in children, a quantitative evaluation of hearing loss and recruitment is necessary. With the help of DPOAE audiograms and DPOAE growth characteristics of the cochlear impaired ear, additional parameters for a non-cooperative hearing aid adjustment can be provided (Müller and Janssen, 2004).

The influence of the second DPOAE source may be observed when plotting the DPOAE level across frequency with narrow frequency spacing, i.e. ≤ 100 Hz (DPOAE fine-structure). Due to either destructive or constructive superposition of the second source across frequency, a pattern of dips and peaks in the DPOAE fine-structure can be observed in subjects with normal or near-normal hearing (He and Schmiedt, 1993, 1996, 1997; Talmadge et al., 1999; Mauermann et al., 1999a,b). DPOAE fine-structure may give information about the fine-structure of behavioral pure-tone thresholds. However, due to the superposition of the second source, the correlation between the two measures is not clear. Moreover, DPOAE fine-structure is reportedly able to reveal OHC impairment in the very early stage, i.e. early stage hearing loss e.g. due to noise over-exposure or ototoxic drug administration. Peaks and dips of the DPOAE fine-structure and with that apparently the impact of the second source are reported to disappear in patients with increasing hearing loss (Mauermann et al., 1999b). In ears with early stage hearing loss the second source is active at least at some cochlear places. For clinical evaluation of DPOAE I/O functions, the interference of the second DPOAE source deteriorates the interpretability and the accuracy of deduced measures such as DPOAE threshold and compression, so that an elimination of the second source is expected to improve their reliability. This can be done to a certain extent by suppressing the second DPOAE source by stimulating the ear simultaneously with a supra-threshold level tone with a frequency close to $2f_1-f_2$ (Heitmann et al., 1998), by applying a windowing technique (Mauermann and Kollmeier, 2004), by an onset decomposition technique (Dalhoff et al., 2013), or as proposed by PATH MEDICAL by frequency-modulated primary tones (Lodwig, 2012, 2013 a,b). When applying a DPOAE suppression technique high suppressor tone levels have to be used to reduce the impact of the second source. As a consequence, OHCs in the neighborhood of the $2f_1-f_2$ place are affected that lead to a less frequency-specific evaluation of the hearing loss. The disadvantage of the windowing technique is that this method is very time consuming. The onset decomposition technique seems to be able to predict cochlear hearing threshold with high accuracy in a short time. However, this technique has only been applied in normal-hearing subjects and in a limited frequency range (1.4 to 2.6 kHz) so far.

The **FMDPOAE**TM measurement technique was developed by PATH MEDICAL (Lodwig, 2012, 2013 a,b, 2014 a,b) as a means of suppressing the impact of the second DPOAE source. Primary tone frequencies are varied over time, following $f_1(t) = f_{1\text{norm}} + d_1(t)$ and $f_2(t) = f_{2\text{norm}} + d_2(t)$ resulting in $f_{dp}(t) = 2*f_1(t) + f_2(t)$. Primary tone frequencies are shifted between ±50 Hz at 1 kHz and ±100Hz at 4 kHz with a modulation rate of 1.4 to 1.6 Hz. Because of the associated phase shift the impact of the second source, which may deteriorate DPOAE detection and with that also DPOAE threshold estimation especially at fine structure minima, is reduced. Moreover, due to frequency modulation the number of stimulated OHCs is increased resulting in a higher DPOAE level. Thus, FMDPOAE does not need any additional stimuli and does not extend test time. FMDPOAE seems to have the potential to improve both screening and diagnostic DPOAE testing performance significantly.

DPOAEs are reported to be measurable at a conductive hearing loss up to 20 to 25 dB HL and a cochlear hearing loss of up to 40 to 50 dB HL, representing approximately the range of the cochlear amplifier (Davis, 1983; Ruggero *et al.*, 1997).

In general, the following are the main clinical applications with DPOAE: follow-up diagnostics after (newborn) hearing screening, assessment of loss of sensitivity (cochlear threshold) and compression (recruitment) of OHC amplifiers, and detection of early stage cochlear impairment during noise over-exposure or ototoxic drug administration, and hearing aid fitting especially in young children.

DPOAE newborn hearing screening is usually performed in the mid-frequency region (e.g. between 1.5 and 4 kHz). High-frequency testing is less reliable due to the standing-wave problem and limitations of the electro-acoustic transducer. Below 1 kHz, the SNR is insufficient. The objective of hearing screening in childhood is to identify hearing impairments which are not apparent and may cause significant handicap for the child concerned. Preschool hearing screening tests should provide more frequency-specific and quantitative information on the hearing loss as available in a DPOAE audiogram.

Two DPOAE workflows are available:

- DPOAE Quick/Diagnostic allows for a frequency-specific measurement of DPOAEs at one or multiple stimulus levels. A response detection status (valid/invalid response) is delivered for each selected frequency/stimulus level combination. DPOAE Quick refers to a measurement at one screening level, whereas DPOAE Diagnostic refers to a measurement at multiple levels. User-definable criteria can be set, i.e. SNR criteria (6, 9 or 12 dB) and for DPOAE Quick an overall pass criterion (number of valid responses: x out y).
- **DPOAE Threshold** allows for a frequency-specific determination of cochlear hearing thresholds by measuring DPOAEs with adaptive level setting. Estimating cochlear hearing loss by means of extrapolated DPOAE I/O-functions provides a quantitative assessment of outer hair cell function/dysfunction. An automatic frequency and level optimization procedure is conducted before each measurement at a specific frequency.

DPOAE tests can be performed binaurally if two ear probes are connected. In doing this, test time is reduced by a factor of two. With the above standard modules, DPOAEs can be measured at frequencies f_2 from 1 to 8 kHz (for DPOAE Diagnostic with *DPHIRES* the frequency range is extended).

For Sentiero Desktop devices and *DPOAE Diagnostic* pressurized DPOAEs are available. This feature allows measuring DPOAEs with static pressure offset in order to compensate for a shift in maximum middle ear compliance. The use of pressurized DPOAEs may improve DPOAE detectability in subjects with abnormal tympanogram (see e.g. Zebian *et al.*, 2013; Beck *et al.*, 2016).

Several module extensions are available, which provide specific additional features:

- **Multichannel DPOAE** allows for a simultaneous measurement of DPOAEs at multiple primary tone pairs (Zurek and Rabinowitz, 1993; Lodwig, 2013b). In order to avoid overlap of the traveling waves of the tone pairs (f_{11} : f_{21} , f_{12} : f_{22} , f_{13} : f_{23} , ...) on the basilar membrane, tone pair frequencies f_{1i} : f_{2i} should have a distinct distance of at least one octave. Frequency distance of the primary tones is controlled automatically.
- **FMDPOAE**TM allows measurement of DPOAEs with frequency-modulated primary tones, which is meant to reduce the influence of the second DPOAE source and therefore is expected to improve reliability of DPOAE detection and hence hearing threshold estimates.

Extension to *DPOAE Diagnostic* only:

- **DPOAE High Resolution** allows DPOAE measurement at user-definable start and stop frequencies from 0.8 to 10 kHz and user-definable linear or logarithmic step size. This allows on the one hand to measure DPOAE grams in a user-definable frequency range and on the other hand to measure DPOAE fine structure with narrow frequency spacing.

PRACTICAL USE

Select the DPOAE test that you would like to perform from the module selection screen, i.e. DPOAE Quick, DPOAE Diagnostic or DPOAE Threshold. If more than one OAE test is licensed, DPOAE Quick, DPOAE Diagnostic, and DPOAE Threshold can be found in the OAE section. Select the preset that you would like to perform. If necessary, change the parameters (e.g. f_2 frequency; for DPOAE Quick/Diagnostic: L2 level, SNR criterion; for DPOAE Quick: overall pass criterion; for DPOAE Diagnostic: L_2/L_1 setup, minimum L_{dp} , timeout) and the preset name as required. For L_2/L_1 level setup, PATH MEDICAL recommends to use the Auto setting which provides optimized primary tone levels similar to the scissor paradigm by Kummer et al. (2000). You can either select adaptive timeout (timing is controlled by the device) or manual timeout with a minimum and maximum measurement time. Please consider that a doubling of measurement time corresponds to an increase in SNR of up to 3 dB. If licensed, FMDPOAE and Multichannel DPOAE mode may be activated in order to improve DPOAE reliability and measurement time, respectively. For DPOAE Diagnostic, the DPOAE High Resolution option can be used with linear or logarithmic step size for measuring DPOAE grams in a user-definable frequency range or for measuring DPOAE fine structure. If fine structure properties are meant to be investigated, FMDPOAE should be deactivated. In all other cases, it is recommended to activate FMDPOAE.

The subject should be calm and sitting comfortably in a chair or lying on a bed. For babies, try to test the subject during sleep. Make sure that a valid ear probe (e.g. EP-DP, EP-VIP) is connected and that the cable of the ear probe does not rub against any material as e.g. clothes. Select an ear tip with appropriate size matching the probe tip size and the subject's ear canal size. Make sure that the ear

probe is inserted without any leakage between ear probe and ear canal. In cases where the ear probe is not inserted properly or ear probe channels are blocked a warning will appear.

Select the test ear (*Right+Left*: simultaneous measurement of right and left ear). After successful ear probe calibration (see *Figure 30*), the measurement starts. For further information on ear probe calibration please refer to section 3.2.2: Transient Evoked Otoacoustic Emissions (TEOAEs). The ear probe calibration view mode can be configured in DPOAE Preferences. The DPOAE stimuli are presented according to the parameter settings and the response is detected via the ear probe microphone. During the DPOAE Quick/Diagnostic measurement, the following screen items are displayed:

DPOAE Quick (see Figure 35, top):

- DPOAE validity bar ① for all selected f_2 at the currently tested L_2 (the bar color corresponds to the ear: red: right ear, blue: left ear)

DPOAE Diagnostic (see Figure 35, bottom):

- DPOAE validity 1 and timeout bar 4 for currently tested f_2/L_2 combinations
- Response status matrix 5 (displayed for right or left ear when pressing the screen on the left or right screen side, respectively): indicates for each f_2/L_2 if a DPOAE is valid (green check mark), invalid (red \emptyset) or skipped (gray \emptyset).
- DPOAE gram (10) (DPOAE high resolution only see *Figure 36*) (displayed for right or left ear when pressing the screen on the left or right screen side, respectively)

DPOAE Quick/Diagnostic

- Overall progress bar ②
- Noise bar (3)

DPOAE Threshold (see Figure 37):

- DPOAE validity (1) and timeout bar (2) for currently tested f_2/L_2 combinations
- Noise bar (3)
- Test status level-frequency-matrix 1: filled vertical bars indicate at which minimum level a valid DPOAE was detected, open box symbols indicate at which f_2/L_2 a DPOAE measurement is currently running
- DPOAE audiogram (12) (displayed for right or left ear when pressing the screen on the lower or upper screen side, respectively)

Please note that the *FMDPOAE* option does not influence the measurement screen. For binaural and multichannel measurements, the user interface is adapted so that all data recorded simultaneously can be seen on a single screen.

If the DPOAE validity bar achieves full amplitude, a valid DPOAE is detected. If the progress bar achieves full amplitude, the measurement timeout is reached. The noise floor is an indicator for the adequacy of test conditions. With good test conditions, the noise floor bar should remain in the lower quarter. If the noise floor bar is higher (the noise bar gradually turns red), ambient noise levels are too high, the patient may be not calm enough (swallowing, moving), or the ear probe cable may rub e.g. against clothes. During the measurement the stimulus stability is observed by means of a low frequency probe tone. If the stimulus stability deteriorates, e.g. due to movement of the ear probe, the ear probe calibration procedure restarts. After successful re-calibration, the measurement continues.

Tests can be paused and continued after a pause ⑦ or stopped ⑥. Also, current f_2/L_2 measurements can be skipped ⑧. For *DPOAE Quick*, all currently ongoing f_2/L_2 measurements are skipped. If a binaural measurement is performed, the user can decide at which ear to skip the currently ongoing f_2/L_2 measurements. For *DPOAE Diagnostic*, a single running f_2/L_2 measurement can be selected. A cartoon mode (see Figure 31) is available by pressing the *plane* button ⑨. The cartoon mode is especially meant for focusing a child's interest and hence improving noise conditions during a test.




Figure 35: DPOAE Quick/Diagnostic measurement (top: DPOAE Quick, middle row: DPOAE Diagnostic; left: monaural test, middle: binaural test, right: binaural, skip current f_2/L_2 measurement; bottom left: DPOAE Diagnostic with current response status matrix; bottom right: DPOAE Quick retest views)

After the *DPOAE Quick* test is finished, DPOAE measurements at specified frequencies can be repeated by pressing the *retest* button 3 in the bar graph view. By default the frequencies with invalid responses are marked for retest but any frequency can be selected. The retest is started by pressing the *play* button 4. If the *retest failed frequency* button 5 is pressed in the overall result view, DPOAEs are immediately retested at the frequencies with invalid responses. *DPOAE Diagnostic* can be configured to automatically repeat f_2/L_2 combinations with invalid responses after the *DPOAE Diagnostic* test is finished.




Figure 36: DPOAE High Resolution measurement (left: monaural; middle: binaural; right: DPOAE gram)




Figure 37: DPOAE Threshold measurement (left: monaural; middle: binaural; right: audiogram)

The following screen items are displayed on the result screen for the respective modules:

DPOAE Quick (see Figure 38):

- Overall result (dependent on selected x out of y overall result criterion)
- DPOAE and noise level bars for each tested f_2 for the given L_2
- DPOAE gram: plots DPOAE and noise levels for all f_2 at the given L_2
- Result data table (including L_2 , f_2 , DPOAE level L_{dp} , noise floor level L_{nf} , and SNR)

Please note that the different result views can be switched by sliding a finger on the screen horizontally. Relevant test parameters (SNR criterion, FMDPOAE, multichannel) are shown in gray at the bottom of the result screen together with the currently selected result view index (filled circle in

row of circles: ••••• = first page of four pages). The preferred initial result screen can be selected in the *DPOAE Preferences* section in the device settings.




Figure 38: DPOAE Quick result (top left: overall result; top right: DPOAE and noise level bars for each tested f_2 ; bottom left: DPOAE gram; bottom right: data table)

DPOAE Diagnostic (see Figure 39):

- DPOAE validity matrix: indicates for each f_2/L_2 if a DPOAE is valid (green check mark), invalid (red \emptyset) or skipped (gray \emptyset)
- DPOAE gram bar graph: plots DPOAE and noise levels as vertical bars and shows their numeric values together with the resulting SNR for each f_2 at the selected L_2 . The DPOAE gram bar graph can be run through all L_2 by pressing the result screen.
- DPOAE gram: plots DPOAE levels for all f_2 at all L_2 (different color for each L_2). The DPOAE gram can be run through all L_2 by pressing the result screen. For a single L_2 , the DPOAE and noise levels for all f_2 at the selected L_2 are plotted.
- Result data table (including L_2 , f_2 , DPOAE level L_{dp} , noise floor level L_{nf} , and SNR)



Figure 39: DPOAE Diagnostic result (top left: DPOAE validity matrix; top middle: DPOAE and noise level bars for each tested f_2 ; top right, bottom left, bottom middle: DPOAE gram; bottom right: data table)

DPOAE High Resolution (see Figure 40):

- DPOAE gram: plots DPOAE levels for all f_2 at all L_2 (different color for each L_2). The DPOAE gram can be run through all L_2 by pressing the result screen. For a single L_2 , the DPOAE and noise levels for all f_2 at the selected L_2 are plotted.
- Result data table (including L_2 , f_2 , DPOAE level L_{dp} , noise floor level L_{nf} , and SNR)



Figure 40: DPOAE High Resolution result (left, middle: DPOAE gram; right: data table)

DPOAE Threshold (see Figure 41):

- DPOAE audiogram, i.e. estimated DPOAE thresholds $L_{\rm dpth}$ plotted across f_2 . The type of DPOAE threshold estimation is indicated by means of different symbols. Estimated thresholds are limited to 50 dB HL. If one or two valid points are available, thresholds are estimated by subtracting an experimental offset. If more than two valid points are available, thresholds are estimated by means of linear regression analysis (see *Figure 34*). When pressing the result screen, also the valid DPOAE f_2/L_2 combinations are shown in the audiogram graph as red (right ear) or blue (left ear) dots.
- DPOAE gram: plots DPOAE levels for all f_2 at all L_2 (different color for each L_2). The DPOAE gram can be run through all L_2 by pressing the result screen. For a single L_2 , the DPOAE and noise levels for all f_2 at the selected L_2 are plotted.
- DPOAE I/O function: plots DPOAE levels for all L_2 at a selected f_2 . The DPOAE I/O function can be run through all f_2 by pressing the result screen.
- DPOAE threshold table (including DPOAE threshold L_{dpth} and minimum L_2)
- Result data table (including L_2 , f_2 , DPOAE level L_{dp} , noise floor level L_{nf} , and SNR)




Figure 41: DPOAE Threshold result (top left/middle: DPOAE thresholds; top right: DPOAE gram; bottom left: DPOAE I/O function; bottom middle: threshold table; bottom right: data table)

Please note that the *FMDPOAE* and *Multichannel DPOAE* options do not influence the layout of the result screen.

LITERATURE

- Beck DL, Speidel S, Arrue Ramos J, Schmuck C (2016): Otoacoustic Emissions and Pressurized OAEs, Hearing Review 23(7), p. 30-33
- Boege P, Janssen T (2002): Pure-tone threshold estimation from extrapolated distortion product otoacoustic emission I/O-functions in normal and cochlear hearing loss ears, Journal of the Acoustical Society of America 111, p. 1810-1818.
- Brown AM, Harris FP, Beveridge HA (1996): Two sources of acoustic distortion products from the human cochlea, Journal of the Acoustical Society of America 100, p. 3260-3267.
- Brownell WE, Bader CR, Bertrand D, de Ribaupierre Y (1985): Evoked mechanical responses of isolated cochlear outer hair cells, Science 227, p. 194-196.
- Dalhoff E, Turcanu D, Vetesnik A, Gummer AW (2013): Two-source interference as the major reason for auditory-threshold estimation error based on DPOAE input-output functions in normal-hearing subjects, Hearing Research 296, p. 67-82.
- Davis H (1983): An active process in cochlear mechanics, Hearing Research 9, p. 79-90.
- Dorn PA, Konrad-Martin D, Neely ST, Keefe DH, Cyr E, Gorga MP (2001): Distortion product otoacoustic emission input/output functions in normal-hearing and hearing-impaired human ears, Journal of the Acoustical Society of America 110, p. 3119-3131.
- Gorga MP, Stover L, Neely ST, Montoya D (1996): The use of cumulative distributions to determine critical values and levels of confidence for clinical distortion product otoacoustic emission measurements, Journal of the Acoustical Society of America 100, p. 968-977.
- Gorga MP, Nelson K, Davis T, Dorn PA, Neely ST (2000): Distortion product otoacoustic emission test performance when both 2f1-f2 and 2f2-f1 are used to predict auditory status, Journal of the Acoustical Soceity of America 107, p. 2128-2135.
- Gorga MP, Neely ST, Dorn PA, Hoover BM (2003): Further efforts to predict pure-tone thresholds from distortion product otoacoustic emission input/output functions, Journal of the Acoustical Society of America 113, p. 3275-3284.
- He NJ, Schmiedt RA (1993): Fine structure of the 2f1-f2 acoustic distortion product: changes with primary level. Journal of the Acoustical Society of America 94, p. 2659-2669.
- He NJ, Schmiedt RA (1996): Effects of aging on the fine structure of the 2f1-f2 acoustic distortion product. Journal of the Acoustical Society of America 99, p. 1002-1015.
- He N, Schmiedt RA (1997): Fine structure of the 2 f1-f2 acoustic distortion products: effects of primary level and frequency ratios, Journal of the Acoustical Society of America 101, p. 3554-3565.
- Heitmann J, Waldmann B, Schnitzler HU, Plinkert PK, Zenner HP (1998): Suppression of distortion product otoacoustic emissions (DPOAE) near $2f_1$ - f_2 removes DP-gram fine structure Evidence for a secondary generator, Journal of the Acoustical Society of America 103, p. 1527-1531.
- Janssen T, Müller J (2008): Otoacoustic emissions as a diagnostic tool in a clinical context. In: Manley GA, Fay RR, Popper AN. Springer Handbook of Auditory Research (Vol 30). Heidelberg, Germany
- Janssen T, Kummer P, Arnold W (1998): Growth behavior of the 2 f1-f2 distortion product otoacoustic emission in tinnitus, Journal of the Acoustical Society of America 103, p. 3418-3430.
- Janssen T, Boege P, Mikusch-Buchberg J, Raczek J (2005a): Investigation of potential effects of cellular phones on human auditory function by means of distortion product otoacoustic emissions, Journal of the Acoustical Society of America 117, p. 1241-1247.
- Janssen T, Gehr DD, Klein A, Müller J (2005b): Distortion product otoacoustic emissions for hearing threshold estimation and differentiation between middle-ear and cochlear disorders in neonates, Journal of the Acoustical Society of America 117, p. 2969-2979.

- Janssen T, Niedermeyer HP, Arnold W (2006) Diagnostics of the cochlear amplifier by means of distortion product otoacoustic emissions. ORL: Journal for Otorhinolaryngology and its Related Specialties 68, p. 334-339.
- Johnsen NJ, Elberling C (1982a): Evoked acoustic emissions from the human ear. I. Equipment and response parameters, Scandinavian Audiology 11, p. 3-12.
- Johnsen NJ, Elberling C (1982b): Evoked acoustic emissions from the human ear. II. Normative data in young adults and influence of posture 11, p. 69-77.
- Keefe DH, Bulen JC, Arehart, KH, Burns EM (1993): Ear-canal impedance and reflection coefficient in human infants and adults, Journal of the Acoustical Society of America 94, p. 2617-2638.
- Kemp DT, Bray P, Alexander L, Brown AM (1986): Acoustic emission cochleography practical aspects, Scandinavian Audiology Supplementum 25, p. 71-95.
- Kummer P, Janssen T, Arnold W (1998): The level and growth behavior of the 2 f1-f2 distortion product otoacoustic emission and its relationship to auditory sensitivity in normal hearing and cochlear hearing loss, Journal of the Acoustical Society of America 103, p. 3431-3444.
- Kummer P, Janssen T, Hulin P, Arnold W (2000): Optimal L(1)-L(2) primary tone level separation remains independent of test frequency in humans, Hearing Research 146, p. 47-56.
- Lodwig A (2012): Frequenzmodulierte DPOAE (in German: Frequency modulated DPOAE), AGERA annual meeting.
- Lodwig A (2013a): Frequency modulated DPOAE, IERASG annual meeting.
- Lodwig A (2013b): Mehrkanal-FMDPOAE (in German: Multichannel FMDPOAE), AGERA annual meeting.
- Lodwig A (2014a): Finding hearing threshold with Multichannel FMDPOAE, AudiologyNow annual meeting
- Lodwig A (2014b): Method and Apparatus for measuring Distortion Product Otoacoustic Emissions (DPOAE) by means of frequency modulated stimuli, US patent 13/657,158
- Lodwig A (2015): Multi-Stimulus suppression in DPOAE vs. ASSR, IERASG annual meeting
- Lonsbury-Martin BL, Martin GK (1990): The Clinical Utility of Distortion-Product Otoacoustic Emissions, Ear and Hearing 11(2), p. 144-154
- Mauermann M, Kollmeier B (2004): Distortion product otoacoustic emission (DPOAE) input/output functions and the influence of the second DPOAE source, Journal of the Acoustical Society of America 116, p. 2199-2212.
- Mauermann M, Uppenkamp S, van Hengel PWJ, Kollmeier B (1999a): Evidence for the distortion product frequency place as a source of DPOAE fine structure in humans I. Fine structure and higher order DPOAE as a function of the frequency ratio f2/f1, Journal of the Acoustical Society of America 106, p. 3473-3483.
- Mauermann M, Uppenkamp S, van Hengel PWJ, Kollmeier B (1999b): Evidence for the distortion product frequency place as a source of distortion product otoacoustic emissions (DPOAE= dine structure in humans. II. Fine structure for different shapes of cochlear hearing loss, Journal of the Acoustical Society of America 106, p. 3484-3491.
- Müller J, Janssen T (2004): Similarity in loudness and distortion product otoacoustic emission input/output functions: implications for an objective hearing aid adjustment, Journal of the Acoustical Society of America 115, p. 3081-3091.
- Müller J, Janssen T, Heppelmann G, Wagner W (2005): Evidence for a bipolar change in distortion product otoacoustic emissions during contralateral acoustic stimulation in humans, Journal of the Acoustical Society of America 118, p. 3747-3756.

- Neely ST, Gorga MP, Dorn PA (2003): Cochlear compression estimates from measurements of distortion-product otoacoustic emissions. Journal of the Acoustical Society of America 114, p. 1499-1507.
- Ruggero MA, Rich NC, Recio A, Narayan SS, Robles L (1997): Basilar-membrane responses to tones at the base of the chinchilla cochlea, Journal of the Acoustical Society of America 101, p. 2151-2163.
- Shera CA, Guinan JJ Jr. (1999): Evoked otoacoustic emissions arise by two fundamentally different mechanisms: a taxonomy for mammalian OAE, Journal of the Acoustical Society of America 105, p. 782-798.
- Smurzynski J, Janssen T (2015): Protocols of DPOAE Measurements Aimed at Reducing Test Time, American Auditory Society Scientific and Technology Meeting 2015-05, Poster #113 DX05
- Smurzynski J (2016): Frequency modulated distortion-product otoacous-tic emission (FMDPOAE) tests aimed for improving diagnostic perfor-mance, World Congress of Audiology annual meeting
- Talmadge CL, Long GR, Tubis A, Dhar S (1999): Experimental confirmation of the two-source interference model for the fine structure of distortion product otoacoustic emissions, Journal of the Acoustical Society of America 105, p. 275-292.
- Whitehead ML, Lonsbury-Martin BL, Martin GK (1992): Evidence for two discrete sources of 2f1-f2 distortion-product otoacoustic emission in rabbit: I. Differential dependence on stimulus parameters, Journal of the Acoustical Society of America 91, p. 1587-1607.
- Whitehead ML, McCoy MJ, Lonsbury-Martin BL, Martin GK (1995a): Dependence of distortion-product otoacoustic emissions on primary levels in normal and impaired ears. I. Effects of decreasing L2 below L1. Journal of the Acoustical Society of America 97, p. 2346-2358.
- Whitehead ML, Stagner BB, McCoy MJ, Lonsbury-Martin BL, Martin GK (1995b): Dependence of distortion-product otoacoustic emissions on primary levels in normal and impaired ears. II. Asymmetry in L1,L2 space, Journal of the Acoustical Society of America 97, p. 2359-2377.
- Williams DM, Brown AM (1995) Contralateral and ipsilateral suppression of the 2f1-f2 distortion product in human subjects. Journal of the Acoustical Society of America 97, p. 1130-1140
- Zebian M, Schirkonyer V, Hensel J, Vollbort S, Fedtke T, Janssen T (2013): Distortion product otoacoustic emissions upon ear canal pressurization, Journal of the Acoustical Society of America 133(4), p. 331-337
- Zurek PM, Rabinowitz WM (1993): Method for testing adequacy of human hearing, U.S. patent #5,267,571.

3.2.4 Overview: Auditory Evoked Potentials

Auditory evoked potentials (AEP) are a tool for diagnosing auditory processing along the ascending pathway from the outer ear up to the auditory cortex including the diagnosis of neural disorders. AEPs are generated in the cochlea, pass on to the cochlear nerve, through the cochlear nucleus, superior olivary complex, lateral lemniscus, to the inferior colliculus in the midbrain, on to the auditory cortex. AEPs are very small event-related electrical voltage potentials that are triggered by an acoustic stimulus and that are recorded as far field potentials via electrodes from the scalp or as near field potentials e.g. with an electrode close to the eardrum. Far field potentials are recorded with some distance from their source, whereas near field potentials are recorded close to their source. A recording channel requires one electrode as ground and two electrodes to measure the desired potential.

AEP recording is a derivative of electroencephalography (EEG). EEG recordings include various wave patterns: *delta* (4 Hz, seen in babies and in adults during slow wave sleep), *theta* (4-7 Hz, seen in young children and in adults e.g. during relaxation and wake-up), *alpha* (7-14 Hz, seen during relaxation and while closing the eyes), *beta* (15-30 Hz, seen during active movements and concentration), and *gamma* (30-100 Hz, seen during movements and cognitive processes).

AEP recordings deviate from standard EEG recordings in the fact that a response signal is averaged time-locked to an acoustic stimulus. The acoustically triggered firing of neurons results in very low electrical potentials. The small amplitudes require signal amplification and cancelling of unwanted external signals (e.g. via differential amplification, i.e. common mode rejection (CMR); filtering). In order to further minimize the influence of unwanted external signals, algorithms for artefact rejection are applied. Artefact rejection can be performed by elimination of high noise level buffers or by weighting each buffer dependent on its noise content. In addition, with enough repetitions of an acoustic stimulus, signal averaging allows the AEP response to emerge from the superimposing background noise, i.e. spontaneous neural firing and other interferences as biological artefacts (e.g. due to eye blinks and movements, swallowing, and other cardiac, muscle, and brain activity; see also EEG wave patterns), and environmental artefacts (e.g. due to electromagnetic radiation of 50/60 Hz depending on the local power system's frequency). The idea of averaging is that the signal is constant and hence the same in each buffer whereas the noise is random and hence changes in each buffer. Adding all the buffers increases the signal and reduces the noise. Theoretically, the improvement in SNR is proportional to the square root of the number of samples that are averaged, limited by technical properties (e.g. quantization noise). Besides environmental conditions, the response detection is dependent on electrode placement, electrode type (surface electrode, needle electrode), electrode cable shielding, and technical aspects such as the amount of amplification and filter properties.

AEPs can be differentiated according to the response latency, i.e. the time gap between stimulus onset and response (see *Figure 42*). With increasing latency the neural generator moves to more central areas.

Early latency responses (= Brainstem Evoked Response Audiometry, BERA) include potentials that can be recorded in a time frame of about 0 to 10 ms after stimulus onset. The physiological source of these potentials is located between the cochlea, the auditory nerve, and various regions of the auditory brainstem (superior olivary complex, lateral lemniscus, and inferior colliculus). Typical clinical application areas are: (newborn) hearing screening, topological diagnostics (together with

tympanometry and otoacoustic emissions), neuro-monitoring, and determination of hearing thresholds especially in small children and other non-cooperative subjects. The following test methods are commonly used for clinical applications:

Auditory Brainstem Responses (ABR) are recorded with surface electrodes on the scalp. Clicks, chirps or tone bursts are typically used. The response latency depends on the stimulus type and level, technical conditions (e.g. transducer type), and the subject's age and hearing status. The recorded potentials contain five to seven waves roughly associated with specific sites along the auditory pathway. ABRs are not affected by sleep. For more information on ABR, please refer to section 3.2.5: Auditory Brainstem Reponses (ABRs).

Auditory Steady-State Responses (ASSR) are recorded with surface electrodes on the scalp. As stimulus an amplitude-modulated signal (tone burst, chirp) is typically used. The repetition rate of the stimulus (e.g. 40 Hz, 80 Hz) is selected in such a way that the transient responses overlap in the brainstem, hence delivering a steady-state response. A response is detected according to statistical properties so that no visual analysis is required by the examiner. Depending on the repetition rate ASSRs may be affected by sleep (40 Hz) or not (80 Hz). For more information on ASSR, please refer to section *3.2.6: Auditory Steady-State Responses (ASSR)*.

Electrocochleography (ECochG) is a procedure where neural activity is recorded by placing a near-field electrode in the ear canal close to the eardrum (extratympanic recording) or with a needle electrode that is perforating the eardrum to rest on the cochlear promontory (transtympanic recording). The response latency amounts to about 1 to 3 ms. For more information on ECochG (module in preparation), please refer to section 3.2.7: Electrocochleography (ECochG).

Middle latency responses (MLR) are potentials that can be recorded from the scalp in a time frame of about 10 to 50 ms after stimulus onset. The physiological source of these potentials is located at the thalamus and the primary auditory cortex. The signal is comprised of multiple maxima with positive (P) and negative (N) polarity called PO, Na, Pa, Nb, and P1. Late latency responses (= Cortical Evoked Response Audiometry, CERA) are potentials that can be recorded in a time frame exceeding 50 ms after stimulus onset. The physiological source of these potentials is the auditory cortex. The signal is typically much larger than for early and middle latency responses and comprises waves P1, N1, P2, and N2. Middle and late latency responses are highly susceptible to the subject's vigilance and attention.

Other AEP methods include mismatch negativity (MMN), which is a procedure detecting potentials elicited by regular acoustic changes (e.g. in level or frequency) evoked by a sequence of identical sounds that is occasionally interrupted by a differing sound. MMN procedures may be independent of cognitive processing or not, i.e. the subject needs to attend to the stimulus difference (e.g. P300). For more complex stimulus differences as e.g. semantic incongruities in speech samples, the subject's language processing skills are required (e.g. N400).




Figure 42: Overview of short, middle and late latency AEP waves

PRACTICAL USE

In general, it is recommended to check the functionality of the electrode cable regularly. The electrode cable can be tested by attaching all electrode clips at the electrode testing device or any available conducting metal bolt. When starting an AEP test (e.g. ABR or ASSR), the impedance should be $0 \text{ k}\Omega$ for both the red and white electrode.




Figure 43: AEP electrode cable check with electrode testing device (left: all clips connected; middle: red clip removed; right: black clip removed or red and white clip removed)

When detaching either the red or white electrode clip from the electrode testing device, the impedance of the red or white electrode, respectively, changes to *open*, i.e. there is no conducting connection between the red or white electrode, respectively, against the black electrode. When detaching the black electrode clip, both impedances are *open*. If there is any deviation from this behavior or if you suspect any dysfunction, please retry with another electrode cable and/or contact your distributor.

LITERATURE

- American Speech-Language-Hearing Association (1987): Short latency auditory evoked potentials, online: http://www.asha.org/policy/RP1987-00024/.
- Kraus N, Nicol T (2009): Auditory evoked potentials, In: Encyclopedia of Neuroscience (editors: Binder MD, Hirokawa N, Windhorst U), Springer, p. 214-218.

3.2.5 Auditory Brainstem Reponses (ABRs)

Auditory brainstem responses (ABRs) recorded from electrodes placed on the scalp represent far field potentials generated by the fiber tracts and nuclei of the ascending auditory pathway between the cochlea and the brainstem. The ABR latency epoch consists of up to seven prominent wave peaks which typically appear within the first 10 ms (Jewett and Williston, 1971). Commonly ABRs are elicited by transient stimuli (click, chirp, tone burst). ABRs represent the sum of synchronized neural activity of the ascending auditory pathway. The waves are commonly referred to as Jewett/wave I to VII. In clinical diagnostics mainly waves I, III, and V are analyzed. Wave I and II stem from the auditory nerve, wave III from the cochlear nucleus, wave IV and V from the lateral lemniscus and inferior colliculus (lower brainstem), and wave VI and VII from sub-cortical regions. Due to the fact that the velocity of the traveling wave running into the cochlea is highest in the basal region, click-evoked ABRs mirror mainly the activity of basal (high-frequency) cochlear regions. This is also true for low frequency tone bursts and chirps when high stimulus levels are applied. Typical properties are wave amplitude and latency. Wave amplitudes represent the amplitude difference between a positive peak and a following negative peak typically amounting to less than 2 µV (Picton et al., 1981). Absolute wave latencies represent the time from stimulus onset until the wave peak occurs. Interpeak latencies (IPL) describe the time between the wave peaks. Wave V has the highest, wave I the lowest amplitude. Amplitudes also decrease with decreasing stimulus level. Wave latencies are shorter at higher stimulus levels due to an increasing basalward spread of excitation (Folsom, 1984). ABRs are quite independent of a subject's vigilance (i.e. sleep, attentiveness) (Picton and Hillyard, 1974) and can therefore be performed during sedation or anaesthetization.




Figure 44: Overview of ABR waves

Commonly, ABRs are elicited by stimulating the hearing organ via air conduction (AC) using headphones, insert earphones or loudspeakers. However, ABRs can also be elicited by stimulating the cochlea directly by bypassing the middle ear via bone conduction (BC). In this case, a bone conductor is placed on the mastoid (or forehead) for stimulating the sensory cells in the cochlea. A

hearing loss caused by middle-ear pathologies (middle ear effusion, Eustachian tube dysfunction, otosclerosis) can be detected, especially in children with doubtful behavioral pure-tone audiograms (Mauldin and Jerger, 1979). Conductive and cochlear hearing loss exhibit different wave amplitude and latency patterns. Especially, wave V latency and threshold are used for audiological diagnostics in children or handicapped persons who are not able to reliably report on their hearing disability. The IPL between wave I and V represent the neural transmission time between cochlea and lower brainstem. IPL is therefore used for assessing retro-cochlear pathologies between cochlea and lower brainstem.

ABR can be elicited with several stimuli as e.g. clicks, chirps or tone bursts. Broadband stimuli as click and broadband chirp are temporally concise stimuli which result in synchronized neural discharges and robust potentials. They are used to stimulate as many as possible sensory cells along the cochlea for getting highest neural activity on the nerve fibers. However, temporal specificity of the stimulus is achieved at the expense of frequency specificity. In contrast, narrowband stimuli as frequency-specific chirp or tone burst comprise limited frequency components. Sensory cells are stimulated at the site within the cochlea corresponding to the stimulus frequency components. Due to the travelling time along the basilar membrane latency varies, i.e. it increases with decreasing frequency. Typically, with decreasing frequency also the amplitude decreases and the waveform gets less sharply defined. With increasing stimulus level the frequency specificity decreases due to an increasing spread of excitation.




Figure 45: ABR stimuli: click and chirp (left: time domain; right: frequency domain)

A **click** (see *Figure 45* – blue curves) is a sound obtained by applying a DC pulse. Its abrupt onset and brief duration lead to a high synchronization of neural activity. However, earphones alter the spectrum of a DC pulse. The auditory system itself also filters the stimulus. Thus, frequency limits are always imposed on click-evoked potentials (Durrant, 1983). A **chirp** (see *Figure 45* – red curves) is a frequency-modulated sinusoidal signal with low frequencies at the beginning and high frequencies at the end of the stimulus or reversed. Due to the frequency modulation, the place-specific travelling time on the basilar membrane can be compensated resulting in an increase of the synchronization of action potentials and with that in higher ABR amplitudes (Dau *et al.*, 2000, Elberling *et al.*, 2007). The advantage of a chirp stimulus over the commonly used click stimulus is that it yields a higher synchronization of action potentials on the nerve fibers, which is especially true at stimulus levels close to threshold. Chirps may contain frequency components covering a major range of the basilar membrane (broadband chirp) or only a rather limited range (narrowband chirp). This may include

low chirp (e.g. 100 to 850 Hz), mid chirp (e.g. 850 Hz to 3 kHz), and high chirp (3 to 10 kHz). A **tone burst** is a short (about 50 to 200 μ s) sinusoidal stimulus pulse which is characterized by carrier frequency, plateau duration, rise and fall time, and the windowing function, which allows fading in and out the stimulus. The windowing of the stimulus results in additional frequency components besides the carrier frequency.

It is important to note, that due to the maturation of the neural pathway absolute latencies and IPL vary during the first year of life (Starr *et al.*, 1977). Adult latencies are reached for wave I at about 6 to 24 weeks and for wave V at about 18 months. IPL I-V decreases during the first 18 months of life. For preterm infants, latencies of all components are typically prolonged compared to term infants. The influence of age on ABRs in the elderly is a matter of controversy. While some authors reported a change in IPL, others claim that there is no delay in neural transmission time with ageing. Also, no significant correlation was observed between the score of speech discrimination tests and IPL.

Since synchronization of neural activity is low in the more apical regions of the cochlea due to the decrease of the travelling wave velocity (Mrowinski, 2009), sensitivity of ABRs is quite low at frequencies below 1000 Hz. ABRs with a broad-band stimulus (click, broadband chirp) do not allow direct frequency-specific assessment of hearing loss and exhibit nearly normal patterns in patients with low- to mid-frequency hearing loss. An ABR threshold is detected at stimulus levels corresponding to the hearing loss at medium to high frequencies. Frequency-specificity is generally improved when using ABR with narrowband stimuli (narrowband chirps, tone burst). However, low-frequency stimulation at high stimulus levels also stimulates basal sensory cells. Thus, the assessment of low-frequency functionality is only possible at low stimulus levels. Ipsilateral masking noise may be used to reduce response contributions from mid- to high-frequency cochlear regions. For getting even more frequency-specific information DPOAEs (at a hearing loss up to 50 dB HL) or ASSRs should to be used.

In patients with hearing loss wave I may be missing. In these patients, the determination of IPL is not possible and thus neural disorders cannot reliably be assessed. Correct placement of electrodes is crucial in order to yield optimal ABR patterns. Due to the fact that electromagnetic fields have impact on ABRs, proper measuring conditions are essential. Proper measuring conditions are present if ABRs are measurable down to 10 dB nHL in a young normal hearing subject.

Clinical applications of ABRs are topological diagnostics, i.e. assessment of peripheral and central sound processing on the auditory pathway up to the lower brainstem, identification of neural disorders on the auditory pathway (vestibulocochlear nerve and lower brainstem lesions: e.g. acoustic neuroma, neural disorder), non-frequency-specific determination of hearing loss within midfrequency region when using broadband stimuli (Gorga *et al.*, 1985), and rough determination of hearing loss within different frequency regions when using frequency-specific stimuli. Despite low frequency specificity ABR may provide a tool for fitting hearing aids in children.

PRACTICAL USE

Select *ABR* from the module selection screen. If more than one AEP test is licensed, *ABR* can be found in the *AEP* section. Select the preset that you would like to perform. If necessary, change the parameters (e.g. stimulus type, stimulus polarity, contralateral masking noise, stimulus level, stimulus rate, number of averages, noise stop criterion, Spread Spectrum, Auto Proceed, Auto Stop, age group for normative latency areas) and the preset name as required.

The different parameter options and possible applications are explained in the following:

Stimulus type:

- Click at low stimulus rates (≤ 20 Hz) evokes clear waves I, III, and V and can be used to obtain information on amplitude and latency. Use click for neurological issues, where the evaluation of inter-peak latency is required. Due to poor frequency specificity thresholds can be assessed only qualitatively. For hearing screening click-ABR measurements are conducted by using a fixed screening level (e.g. 35 dB nHL).
- Chirp evokes higher response amplitudes than click for wave III and V. In contrast, wave I is typically hard to identify. Use chirps for hearing threshold determination, because only wave V is evaluated.
- Narrowband (low-, mid-, high-) chirps provide latency information and are more frequency-specific than broadband stimuli (click, broadband chirp). Hence, they may be used for a more frequency-specific hearing threshold determination. However, the response amplitude is typically lower than for broadband stimuli, so that the response is harder to detect in the time domain. Please note that for low- and mid-chirps ipsilateral masking noise is presented in order to reduce response contributions from more basal cochlear regions.
- Tone bursts are more frequency-specific than narrow-band chirps but yield lower amplitudes (especially at lower levels).

Stimulus polarity:

- Alternating polarity helps reduce the stimulus artefact that is generated by the transducer itself (especially recommended for bone conduction measurements). Alternating polarity provides a broader, rounded wave V peak.
- Rarefaction and condensation provides a more peaked response and may yield higher amplitude for wave I. Latency difference between condensation and rarefaction is nearly identical in normal hearing adults. However, responses to condensation and rarefaction clicks may considerably differ in patients with cochlear hearing loss.

Rate Mode:

- If activated, the test can be performed for a fixed stimulus level at multiple stimulus rates.

Masking noise:

- Contralateral masking is recommended if there is significant asymmetry in hearing loss between ears, i.e. for differences of about 30 to 40 dB (headphones) or 50 to 60 dB (insert earphones). If a bone conductor is used, application of contralateral noise is essential.

Stimulus level:

Up to eight stimulus levels (including repetitions) can be pre-configured. Stimulus levels are given in dB nHL, i.e. relative to the hearing threshold of a collective of normal hearing subjects, which is defined as 0 dB nHL. A stimulus level can be repeated up to three times. A mute stimulus is available for comparative measurements. Measurements start at the highest level. The standard deviation of latency is typically lower and wave amplitude is higher at higher stimulus levels. At lower stimulus levels, wave I (at about 60 dB nHL) and wave III (at about 30 dB nHL) disappear. In rate mode only one stimulus level can be set.

Stimulus rate:

The higher the stimulus rate, the smaller the response amplitude (notably for I to IV and for high levels) and the longer the latency. The latency shift due to stimulus rate is compensated

on the device and hence not visible in the result graph. Higher rates improve efficiency of data collection by reducing measurement time at a fixed number of averages, but jeopardize the identification of a response, particularly in some pathological cases. Wave I and III may disappear at stimulus rates above 50 Hz. If all waves are intended to be evaluated (e.g. for neurological diagnostics) low stimulus rates should be used. For ABR threshold determination also higher stimulus rates may be used. 50 and 60 Hz are not available stimulus rates because these are typical power system frequencies. By pressing the stimulus rate box, a user-defined stimulus rate can be entered. The resulting inter-stimulus interval is displayed on the *settings* screen. In rate mode up to eight stimulus rates (including repetitions) can be selected.

Number of averages:

At fixed measurement conditions, with increasing number of averages, the noise floor decreases (number of averages increased by factor four reduces noise by half), but the measurement time increases. The time per trace and the total measurement time are shown on the *settings* screen.

Spread Spectrum:

If activated, the stimulus rate is slightly varied in order to reduce the influence of electrical interference synchronized to the stimulus rate. Also, ABR amplitude is known to decrease at a constant stimulus rate due to adaptation. Activation of this option is always recommended. Please note that the Spread Spectrum option is always active for stimulus rates exceeding 70 Hz and during binaural stimulation for decoupling the responses from the two channels.

Automated Wave V detection:

If activated, the occurrence of a statistically valid wave V is automatically detected.

Minimum Wave V amplitude:

If *Automated Wave V detection* is activated, the minimum wave V amplitude, which must be available to mark a recognized wave V as valid, can be selected.

Auto Proceed:

If activated, the recording of a trace is stopped as soon as a statistically valid wave V is detected for the given trace. The test then proceeds with the recording of the next trace. Please note that this option is only available if the *automated wave V detection* option is activated.

Auto Stop:

If activated, the test stops if for two consecutive traces no statistically valid wave V could be detected. Please note that this option is only available if the *automated wave V detection* option is activated.

Noise stop criterion:

If activated (i.e. noise stop criterion >0 nV), the recording of a trace is stopped as soon as the residual noise drops below the defined noise threshold and no response is detected. Hence, if activated this option speeds up the recording in case no response is present.


Age group (for normative latencies):

Select the appropriate age group corresponding to the age of the tested subject.

Before the test is started the subject should be instructed about the test procedure. In order to reduce muscle artefacts, the subject should be calm and fully relaxed lying comfortably on a recliner

or bed. It is also recommended that subjects keep their eyes closed during the measurement for reducing artefacts e.g. due to eye blinks. For babies, try to test the subject during sleep. In order to reduce environmental artefacts, conduct the measurement in a room with low electromagnetic radiation, i.e. in an electrically shielded metal booth or any other room without powered-on electric devices (e.g. computer, light, telephone, cell phone, power transformer) in close range to the measurement equipment. An acoustically shielded booth or a quiet room is recommended if ABR is applied for threshold determination at low sound pressure levels.

Make sure that a valid transducer (e.g. headphone, insert earphone, ear probe, bone conductor) and electrode cable are connected. Make sure that the skin is clean at the intended positions of the electrodes. If applicable, thoroughly clean the skin (e.g. using skin prepping gel) in order to remove dirt, oil, and superficial dead skin. Select appropriate electrodes and attach them on the skin of the patient. An electrolyte gel may be put on the electrode contact in order to improve the conductivity of the skin layer, which effectively increases the electrode surface area. Attach the electrode clips of the electrode cable at the correct electrode. The white and red electrode are the recording electrodes, the black electrode is the ground electrode. Do not place the ground electrode near the heart to avoid inducing electrocardiography (ECG) waves. There are several possibilities to place the electrodes. The position of the electrodes affects waveform morphology and latency. Best electrode position for yielding maximum wave amplitudes is vertex (white electrode) and ipsilateral mastoid (red electrode) as shown in Figure 46 (vertical montage).



White electrode:

Vertex or high forehead

Red electrode:

- Ipsilateral mastoid (for monaural recording)
- Nape of the neck (for binaural recording sequentially or concurrently, for bone conduction recording if bone conductor is placed on the mastoid)

Black electrode (ground):

e.g. cheek

Figure 46: ABR electrode positioning with vertical montage

Alternatively to the vertex position, the high forehead position is possible for the white electrode. However, in this case the ABR amplitude slightly decreases. Despite this fact, the forehead is preferred in practice, especially in patients where vertex electrode placement is inconvenient because of hair. Wave I amplitude may be larger in a horizontal montage (white electrode: ipsilateral mastoid; red electrode: contralateral mastoid) than in the standard vertical montage.

Select the test ear (Right+Left: simultaneous measurement of right and left ear – please note that in this case a symmetric electrode montage is recommended and that no contralateral masking noise is available). The electrode impedance measurement starts (see Figure~47). Electrode impedance is dependent on the electrode material and the surface area to which it contacts. The impedance is shown for the red and white electrode, respectively, against the black electrode ①. The ABR test can be started by pressing the play button ② if the impedances are below 6 $k\Omega$ and the difference between red and white electrode impedance is below 3 $k\Omega$. If you do not press the play button for a while, a message appears in order to remind you to start the measurement by pressing the play button. It is possible to configure an automatic start of the ABR test after adequate impedances are available in the AEP Preferences. The measurement is then automatically started if the impedances are below 4 $k\Omega$ and the impedance difference is below 2 $k\Omega$. However, in some cases these

impedances may not be obtained (e.g. in small children cleaning of the skin may not be possible as the child would wake up when rubbing the skin) and an ABR test shall be performed despite the not ideal test conditions. For these cases, a *forward* button ③ appears if the impedances are below 12 $k\Omega$ and the impedance difference is below 6 $k\Omega$. If impedances are worse no ABR test is possible. Please note that in neonates the impedance may be higher due to tenderness of the skin. If impedances are too high, check the electrode cable (see PRACTICAL USE in section *3.2.4: Overview: Auditory Evoked Potentials*) and clean the skin, use conductive gel and wait a couple of minutes until the gel is infiltrated into the skin.




Figure 47: ABR impedance measurement (left: good impedances with play button; middle: good impedances with play button and reminder to start test by pressing the play button; right: poor impedances with forward button)

When performing an ABR test with an ear probe, information about ear probe calibration ① and impedance measurement ② are shown on screen (see *Figure 48*). You can configure the ear probe calibration view mode in *AEP Preferences* to *Simple* or *Expert* mode. If the results from impedance measurement or ear probe calibration are not fully valid, a *forward* button may occur in order to proceed to the ABR test. For further information on ear probe calibration, skip criteria and differences in calibration view mode please refer to *PRACTICAL USE* in section *3.2.2: Transient Evoked Otoacoustic Emissions (TEOAEs)*.



Figure 48: ABR impedance measurement and ear probe calibration (left: good impedances and ear probe fit (simple mode); middle: good impedances and ear probe fit (binaural test, expert mode; right: bad impedance during measurement)

After successful electrode impedance measurement, you can start the test by pressing the *play* (or *forward*) button. The ABR stimulus is presented via the transducer and the response is detected via the electrode cable. During the measurement the electrode impedance is monitored. In case that impedances get too high (e.g. an electrode has fallen off), the test is automatically interrupted (see *Figure 48*). The test can only be continued if the electrode impedances get back to tolerable values.

During the ABR measurement the following information is provided on the screen (see Figure 49):




Figure 49: ABR measurement (top left: standard measurement screen; top middle: extended measurement screen; top right: level selection (for rate mode: stimulus rate selection); bottom left: Jewett marker setting; bottom middle: standard measurement screen with EEG noise graph for low and high noise conditions; bottom right: standard measurement screen for binaural measurement)

ABR traces are shown for each selected stimulus level (for rate mode: for each selected stimulus rate). In the standard measurement screen a maximum of three traces is shown at the same time. The currently tested trace is shown in red (right ear) and/or blue (left ear) ②. Already finished traces are grayed out. Normative data for wave III (red area) and V (green area) is shown if configured in the settings. When pressing the graph area an extended view of up to eight traces is shown. Pressing the graph again will switch back to the standard measurement screen with three traces. By sliding a finger horizontally over the time axis of the ABR traces graph area, the time scale can be adapted. By sliding a finger vertically over the graph area, the amplitude scale can be adapted. If the automatic wave V detection is activated, the wave V status indicator ① reflects the statistically computed probability that a wave V is present in the recorded signal. For binaural measurements, the red bar

refers to the right ear, the blue bar to the left ear. If a statistically significant wave V is detected a red or blue dot is displayed at the trace for the right or left ear, respectively. Please note that the wave V status indicator only represents the statistical detectability of a wave V signal and does not consider normal amplitude or latency.

During an ABR measurement, the pre-configured level selection can be adapted by pressing the *level setup* button ③ (for rate mode the pre-configured stimulus rates can be adapted instead). A level matrix opens and the selected levels are shown. Already tested levels are grayed out. The level setup screen can be left with the *back* ⑥ (without saving the changes) or *OK* button ⑦ (saving the changes). Also, the pre-configured number of averages can be increased by 1000 for the current trace by pressing the *increase averages* button ④. Jewett markers can be set for the currently active trace by pressing the *Jewett marker mode* button ⑤. Initially all markers are shown in gray. The active marker is shown with a circle ④. The active marker can be moved by pressing the *right/left arrow* buttons ⑱ and can be set by pressing the *set marker* button ⑤. The marker is then shown red or blue (for the right or left ear, respectively). When moving the marker, it turns gray again. You can choose which Jewett marker (V, III, and I) should be set by pressing the *Jewett marker selection* button ⑨. The Jewett marker mode can be left by pressing the *back* button ⑥.

The EEG noise bar ⑥ represents biological and external noise which is an indicator for the adequacy of measurement conditions. With good measurement conditions, the EEG noise bar should remain green. With increasing noise, the EEG noise bar gradually turns yellow and red, which indicates that biological and external artefacts may be too high, i.e. the patient may be not calm and relaxed enough (e.g. swallowing, moving, clenched jaw/teeth), or the electromagnetic radiation may be too prominent. When pressing the EEG noise bar, an EEG noise graph ② is shown. You can return to the EEG noise bar view by pressing the EEG noise graph. The progress bar ⑦ shows the progress of the measurement at the currently active trace. If the bar is full, the measurement for the current trace is finished, i.e. the number of averages that has been defined by the user is reached. Depending on the settings, a measurement may be finished before the number of averages has been reached, e.g. if Auto Proceed or Noise Stop Criterion is active. The electrode impedance ⑧ is regularly updated during the measurement and the measurement is paused if impedance gets too high. The test can be manually paused ①, continued after a pause ③, or stopped ⑨. Also, a current measurement can be restarted ⑩ or skipped ②.

After the test is finished, the ABR result is shown (see *Figure 50*). The recorded ABR traces ② are displayed for each stimulus level together with normative areas if a normative data age group has been configured. If automated wave V detection has been activated the estimated wave V is shown as a dot on the trace. Also the statistical validity of the wave V is shown in the wave V status indicator ①. After setting the Jewett markers, the wave name (I, III, and V) is displayed at the respective markers ③. The resolution of the amplitude axis can be changed by sliding a finger up (zoom in) or down (zoom out) on the result graph screen. The resolution of the time axis can be changed by sliding a finger right (zoom in) or left (zoom out) on the time axis of the result graph screen. Please note that signal delays for insert earphones and for different stimulus rates are compensated.

The latencies for wave I, III, and V can be adjusted by the user on the Jewett marker screen ④. The currently active marker is indicated by a circle. The trace can be selected with the *up/down* buttons ⑤. The active marker can be moved along the selected trace with the *right/left arrow* buttons ⑥ and can be set by pressing the *set marker* button ⑦. While moving the marker it is gray. When

setting the marker it turns red or blue (for the right or left ear, respectively). Please note that the order of the markers is always preserved with a minimum distance of 0.5 ms (e.g. when moving the Jewett I marker to the right the Jewett III marker will be moved too as soon as the Jewett I marker reaches the Jewett III marker). The latency corresponding to the active marker and the resulting amplitude is shown as a numeric value below the ABR traces graph. You can choose which Jewett marker (I, III, or V) you wish to change by pressing the *Jewett marker selection* button (§). If any marker has been changed an *OK* button (§) appears in the footer. When leaving the result view, any changes can be accepted by pressing the *OK* button (§) or discarded by pressing the *back* button (§).

By sliding a finger from right to left on the result screen you can move forward to other result views. The three screens present a level-latency-diagram (wave V: circles, wave III: triangles, wave I: squares) 9, numeric values of amplitudes, latencies and inter-peak latencies together with other test-related information (such as impedances, number of averages, and noise for each trace) 10, and general information about the used transducer and settings (as averages, stimulus rate, stimulus type, test options, normative area age group, and total time) 11. In order to ensure adequate measurement conditions noise should be <100 nV after 2000 averages. With increasing averages, noise decreases (e.g. < 70 nV after 4000 averages, < 50 nV after 8000 averages).




Figure 50: ABR result (top left: traces with estimated wave V; top middle: traces with confirmed waves I, III, and V; top right: Jewett marker setting; bottom left: level-latency diagram; bottom middle: measurement data; bottom right: general test information)

For level repetitions, there are additional options for presenting ABR traces. They can be shown sorted by time (test order) or by level and ABR traces of the same level can be displayed summarized

or can be plotted on top of each other in order to evaluate the repeatability of an ABR response for a given level. The available options are shown in *Figure 51*:




Figure 51: ABR result options for repeated levels (top left: overlay traces of same levels; top right: show summed traces for same levels; bottom left: show traces in test order; bottom right: show traces sorted by level)




Figure 52: ABR result for rate mode (left: traces with confirmed waves III and V; right: stimulus rate vs. latency diagram)

If rate mode has been selected, traces are displayed for each tested stimulus rate (instead of stimulus level as in standard mode) and a stimulus rate vs. latency diagram is available (instead of a

level-latency-diagram) (see *Figure 52*). The basic layout of result views in rate mode is similar to the result views in standard mode.

For the interpretation of ABR results, the following case examples may be considered (see *Figure 53* to *Figure 56*):

- **Normal hearing subject** (see *Figure 53*): ABRs with broadband stimuli are typically measurable down to about 10 dB nHL in a normal hearing subject. Wave V latency increases and wave V amplitude decreases with decreasing stimulus level. Wave V latencies are within normal range when an appropriate age group has been selected.




Figure 53: ABR result example for normal hearing subject (left: audiogram; middle: ABR traces; right: latency-level function)

Conductive hearing loss (see Figure 54): In a patient suffering from conductive hearing loss, the effective stimulus level decreases. As a consequence, latencies and amplitudes of all waves are changed. This includes that wave V latency typically increases out of normal latency range. The waves prior to wave V are often not detectable. The inter-peak latency is not affected. In a conductive hearing loss with a flat audiogram the level-latency function is shifted along the level axis by the amount of the conductive hearing loss. ABR threshold is reached at stimulus levels corresponding to the hearing loss at medium to high frequencies.




Figure 54: ABR result example for subject with conductive hearing loss (left: audiogram; middle: ABR traces; right: latency-level function)

- Cochlear hearing loss (see *Figure 55*): The overall effect is dependent on the severity and configuration of the hearing loss, and also on the frequency composition of the stimulus. Wave V latencies are essentially equivalent to those collected in normal hearing subjects as

long as the stimuli are at least 20 dB above the threshold at 4 kHz, the configuration of the hearing loss is not steeply sloping and of mild to moderate in severity. Inter-peak latency I-V is normal. Level-latency functions for these subjects converge on those of normal hearing subjects at high stimulus levels.

Case example 1: ABRs in a patient suffering from a moderate high frequency cochlear hearing loss (click, stimulus rate: 20 Hz): In comparison to a normal hearing subject, ABRs appear with lower amplitude and slightly increased wave V latency. The ABR threshold is reached at stimulus levels corresponding to the hearing loss at medium to high frequencies. The prolongation of latency at low stimulus levels corresponds to the propagation time of the travelling wave along the non-functioning basal part of the cochlea at these levels.

Case example 2: ABRs in a patient suffering from a severe high frequency cochlear hearing loss (click, stimulus rate: 20 Hz): ABR threshold is reached at stimulus levels corresponding to the hearing loss at high frequencies. Low-frequency thresholds cannot be assessed due to low synchronization of nerve fiber activity in the apical region of the cochlea.




Figure 55: ABR result examples for subjects with cochlear hearing loss (top: case example 1; bottom: case example 2) (left: audiogram; middle: ABR traces; right: latency-level function)

Retro-cochlear hearing loss (see *Figure 56*): ABRs in a patient with neural disorders (click, stimulus rate: 10 Hz): The ABR wave pattern is different compared to normal hearing subjects and patients with cochlear or conductive hearing loss with respect to I-V inter-peak latency (IPL). I-V IPL is prolonged (4.4 ms) due to reduced neural transmission time between cochlea and brainstem exceeding the normal I-V interval (female: 3.8 ± 0.2 ms; male: 4.0 ± 0.2 ms). If data from both ears are available, interaural wave V latency differences can be evaluated (not shown). Differences in wave V latency should not amount to more than 0.3 to 0.4 ms between ears if there is no distinct interaural asymmetry in hearing loss.



Figure 56: ABR result example for subject with retro-cochlear hearing loss

Quick ABR

An alternative to the standard ABR module is provided on Sentiero Advanced mainly for screening purposes. Select *QuickABR* from the module selection screen. If more than one AEP test is licensed, *QuickABR* can be found in the *AEP* section. If required, change the parameters. The parameter options are reduced compared to the standard ABR module and only provide a selection of stimulus level and normative latency age group. Stimulus levels are given in dB eHL, which represents an *estimated hearing level* with dB eHL = dB nHL + 10 dB. The offset was selected based on the experience that ABR can typically be recorded only down to 10 dB above hearing threshold. Hence, by introducing the additional offset, the typical ABR detection threshold for normal hearing subjects is considered to be reached at 0 dB eHL. A fixed stimulus setup (stimulus rate, stimulus polarity, number of averages) is applied. In the *QuickABR* module the statistical detection algorithm checks for wave V in the normal latency area only, which is different to the standard ABR module, which checks for a wave V in a broader time frame also including late latencies.

For preparation of the patient (instruction, skin preparation, and electrode montage) please refer to the explanations for the standard ABR module. Make sure that a valid transducer (e.g. headphone, insert earphone, ear probe) and electrode cable are connected. Select the test ear (*Right+Left*: simultaneous measurement of right and left ear). The electrode impedance measurement starts. If the test is performed with one or two ear probes also the ear probe calibration starts. For more information on impedance measurement and ear probe calibration please refer to the explanations for the standard ABR module. After successful electrode impedance measurement, you can start the test by pressing the *play* (or *forward*) button. During the *QuickABR* measurement the following information is provided on the screen (see *Figure 57*):

The recorded ABR trace is plotted for the selected stimulus level ①. The stimulus level is displayed below the ABR trace graph. Normative data for wave III (red area) and V (green area) is shown if configured in the settings. By sliding a finger vertically over the graph area, the amplitude scale can be adapted. The wave V validity indicator ② reflects the statistically computed probability that a wave V with normal latency is present in the recorded signal. For binaural measurements, the red bar refers to the right ear, the blue bar to the left ear. If a valid wave V is detected for one ear side, a green check mark is shown instead of the wave V validity indicator.

The EEG noise bar ③ represents biological and external noise which is an indicator for the adequacy of measurement conditions. With good measurement conditions, the EEG noise bar should remain

green. With increasing noise, the EEG noise bar gradually turns yellow and red, which indicates that biological and external artefacts may be too high, i.e. the patient may be not calm and relaxed enough (e.g. swallowing, moving, clenched jaw/teeth), or the electromagnetic radiation may be too prominent. The progress bar ④ shows the progress of the measurement relative to the maximum number of averages. If the bar is full, the measurement is finished. The measurement may finish earlier if a valid response is detected before reaching the maximum number of averages. The electrode impedance ⑤ is regularly updated during the measurement and the measurement is paused if impedance gets too high. The test can be manually paused ⑦, continued after a pause (pause button turns to play button – not shown in Figure 57), or stopped ⑥.




Figure 57: QuickABR measurement (left: monaural; middle/right: binaural)

After the test is finished, the *QuickABR* result is shown (see *Figure 58*). The recorded ABR trace ① is displayed for the selected stimulus level together with normative areas if a normative data age group is configured. The overall result is shown as a green check mark ② (valid response) or a red symbol ③ (invalid response). The resolution of the amplitude axis can be changed by sliding a finger up (zoom in) or down (zoom out) on the result graph screen. Please note that signal delays for insert earphones are compensated. Additional information on the applied stimulus type and level, number of averages, EEG amplitude, and impedance is shown at the bottom of the result view ④.




Figure 58: QuickABR result (left: valid response; right: no valid response)

E-ABR (module in preparation, not available for firmware 2.2)

Another version of the standard ABR module is the E-ABR module, which allows recording electrically evoked auditory brainstem responses in order to determine if the auditory nerve adequately responds to electrical stimulation. Please note that the E-ABR module is in preparation and therefore only preliminary information on the practical use of the module and on user interfaces can be presented.

The E-ABR module provides the capability to measure ABR responses at patients with a cochlear implant. The electrical stimulus is applied from an external stimulator, usually a cochlear implant fitting system (which is provided by the cochlear implant manufacturer). The electrical stimulus is transmitted directly to the cochlear implant coil where it is forwarded to the cochlear implant electrode and converted to a stimulation of the auditory nerve. The ABR responses are recorded from the patient's scalp via surface electrodes attached to the patient (as for normal ABR). The surface electrodes are attached to an electrode cable, which is connected to the Sentiero Advanced. In order to synchronize the responses correctly, the time sequence of the stimuli needs to be known on the device side. For that reason stimulus trigger impulses are transferred from the cochlear implant fitting system to the measurement device.

The E-ABR module on Sentiero Advanced allows measuring up to 15 traces in one test run and provides automated response detection and Jewett V marking. The basic ABR operation principles are the same as for the standard ABR module.

Select *E-ABR* from the module selection screen. If more than one AEP test is licensed, *E-ABR* can be found in the *AEP* section. Select the preset that you would like to perform. If necessary, change the recording parameters (e.g. automated wave V detection, deactivation of impedance check during the measurement, trigger options) and the preset name as required. Stimulus parameters must be configured at the cochlear implant fitting system. For further information on the cochlear implant fitting system please refer to the respective operation manual or contact the cochlear implant manufacturer.

Before the test is started the subject should be instructed about the test procedure. Make sure that the test conditions (patient position and relaxation, electrode montage, environmental conditions) are appropriate and a valid electrode cable is connected to the Sentiero Advanced. For further information please refer to the explanations for the standard ABR module. Please also make sure that the stimulus configuration on the cochlear implant fitting system is appropriate and an adequate transducer is connected to the cochlear implant fitting system.

Select the test ear. The electrode impedance measurement starts (see explanations for standard ABR module). The ABR test can be started by pressing the *play* button if impedances are appropriate. When starting the play button, the device checks for a valid trigger input. As soon as a trigger signal is recognized the measurement starts synchronized to the trigger signal. During the measurement the electrode impedance is monitored. In case that at least one of the impedances gets too high (e.g. an electrode has fallen off), the test is automatically interrupted. The test can only be continued if the electrode impedances get back to tolerable values.

During the E-ABR measurement the information provided on the screen (see *Figure 59*) is similar as for standard ABR, i.e. the screen shows E-ABR traces (view of all traces appears when clicking the traces graph; the traces are consecutively numbered), noise, the recording progress for the currently

measured trace, and electrode impedances. The trace starts 1 ms before stimulus onset in order to show the electrical artefact of the external stimulator. A trace is finished as soon as the progress bar is full, i.e. the pre-configured number of averages is reached. The measurement is then paused. After setting up the next test configuration on the cochlear implant fitting system, the recording can be started again on Sentiero Advanced by pressing the *play* button. If configured, a trigger pause automatically results in switching to a new trace. The recording of a trace can be skipped or restarted by pressing the *forward* or *back* button. The pre-configured number of averages can be increased by pressing the *increase averages* button. An EEG noise graph can be displayed when pressing the EEG bar. After the test is finished, an E-ABR result is shown comparably to the standard ABR result.




Figure 59: E-ABR measurement

LITERATURE

- American Speech-Language-Hearing Association (1987): Short latency auditory evoked potentials, online: http://www.asha.org/policy/RP1987-00024/.
- Bargen GA (2015): Chirp-Evoked Auditory Brainstem Response in Children: A Review, American Journal of Audiology 24, p. 573-583
- Bell SL, Allen R, Lutman ME (2002): An investigation of the use of band-limited chirp stimuli to obtain the auditory brainstem response, International Journal of Audiology 41, p. 271-278
- Cebulla M, Stürzebecher E (2013): Detectability of newborn chirp-evoked ABR in the frequency domain at different stimulus rates, International Journal of Audiology 52, p. 698–705
- Crumly W (2011): Good practices in auditory brainstem response, part 1, online: http://www.audiologyonline.com/articles/good-practices-in-auditory-brainstem-827.
- Dau T, Wegner O, Mellert V, Kollmeier B (2000): Auditory brainstem responses with optimized chirp signals compensating basilar-membrane dispersion, Journal of the Acoustical Society of America 107(3), p. 1530-1540
- Durrant JD (1983): Fundamentals of sound generation, In: Bases of auditory brainstem evoked responses (editor: Moore EJ), Grune Stratton, p. 15-49.
- Elberling C, Don M, Cebulla M, Stürzebecher E (2007): Auditory steady-state responses to chirp stimuli based on cochlear traveling wave delay, Journal of the Acoustical Society of America 122, p. 2772-2785.
- Folsom RC (1984): Frequency specificity of human auditory brainstem responses as revealed by puretone masking profiles, Journal of the Acoustical Society of America 66, p. 919-914.

- Gorga MP, Worthington DW, Reiland JK, Beauchaine KA, Goldgar DE (1985): Some comparisons between auditory brainstem response thresholds, latencies, and the pure-tone audiogram, Ear and Hearing 6, p. 105-112.
- Jewett DL, Williston JS (1971): Auditory evoked far fields averaged from the scalp of humans, Brain 94, p. 681-696.
- Lu TM, Wu FW, Chang H, Lin HC (2017): Using click-evoked auditory brainstem response thresholds in infants to estimate the corresponding pure-tone audiometry thresholds in children referred from UNHS, International Journal of Pediatric Otorhinolaryngology 95, p. 57-62
- Mauldin L, Jerger J (1979): Auditory brain stem evoked responses to bone-conducted signals, Archives of Otolaryngology 10, p. 656-661.
- Miyamoto RT (1986): Electrically evoked potentials in cochlear implant subjects, Laryngoscope 96(2), p. 178-185.
- Mrowinski D (2009): Elektrische Reaktionsaudiometrie (in German: Electrical response audiometry), In: Praxis der Audiometrie (editors: Lehnhardt E, Laszig R), Thieme.
- Thieme, Stuttgart. Lehnhardt, E. und Laszig, R., p. 203-238.
- Picton TW, Hillyard SA (1974): Human auditory evoked potentials. II. Effects of attention, Electroencephalography and Clinical Neurophysiology 36, p. 191-200.
- Picton TW, Stapells DR, Campbell KB (1981): Auditory evoked potentials from the human cochlea and brainstem, Journal of Otolaryngology 9, p. 1-41
- Starr A, Amlie RN, Martin WH, Sanders S (1977): Development of auditory function in newborn infants revealed by auditory brainstem potentials, Pediatrics 60, p. 831-839.

3.2.6 Auditory Steady-State Responses (ASSR)

Auditory Steady-State Responses (ASSRs) provide frequency-specific and quantitative information on sound processing within the brainstem and sub-cortical regions (Picton et al. 2003; Herdman and Stapells 2003). ASSRs are elicited by means of narrowband stimuli, i.e. an amplitude- (and frequency) modulated sine or a narrowband chirp, which compensates for the time delay on the basilar membrane. Compared to a modulated sine, the use of a narrowband chirp reduces frequency specificity but increases response amplitudes. Depending on the frequency components of the stimulus, different sensory cell regions within the cochlea are stimulated. Depending on the modulation frequency (i.e. for sine: amplitude modulation rate, for chirp: repetition rate – in the following only the term stimulus rate is used) onset neurons are activated in different regions along the auditory pathway. The stimulus rate is selected in such a way that the transient responses overlap at the place of generation, hence delivering a steady-state response. Typical stimulus rates are 40 and 80 Hz. 40 Hz ASSR are generated mainly in the auditory sub-cortex (Mäkelä and Hari, 1987; Pantev et al., 1996; Ross et al., 2003), whereas 80 Hz ASSR are generated in the brainstem (Herdman et al., 2002). Compared to transient AEPs, 40 Hz ASSR include both early and middle latency response components, whereas 80 Hz ASSR consist of early latency responses (Mäkelä and Hari, 1987). 40 Hz ASSR are affected by vigilance and are hence suitable for awake and alert subjects (Galambos et al., 1981), whereas 80 Hz ASSR are not affected by vigilance and are hence suitable for babies and in general for sleeping subjects (including sedated and anaesthetized subjects) (Levi et al., 1983). In newborns, middle latency responses are usually not available, so that the amplitude of 40 Hz ASSR is generally reduced. Moreover, the optimal stimulus rate also depends on the carrier frequency. 40 Hz ASSR response amplitudes increase with decreasing carrier frequency (Rodriguez et al., 1986, Picton et al., 1987), whereas 80 Hz ASSR response amplitudes are maximal at medium carrier frequencies between 1 and 2 kHz (John et al., 2001, Dimitrijevic et al., 2002). For frequencies from 6 kHz onwards, even higher stimulus rates up to 160 Hz can be used (Lodwig and Rosner, 2016). With increasing stimulus rate the background noise during recording of the response decreases. The response can be detected via far-field electrodes within the EEG as a sine signal with a frequency following the stimulus rate. The response can be statistically analyzed in the frequency domain (Stapells et al., 1987; Dobie and Wilson, 1989; Picton et al., 2001).

When using different stimulus rates for different carrier frequencies, a simultaneous measurement at multiple frequencies is possible (Lins and Picton, 1995). The suppression effect for simultaneous multiple-frequency measurements with a minimum distance of one octave between simultaneously tested frequencies is significantly lower for 80 Hz ASSR compared to 40 Hz ASSR (John *et al.*, 1998; Ross *et al.*, 2003). For 40 Hz ASSR, the amplitude may decrease by 50 %, whereas for 80 Hz ASSR, the amplitude may decrease by about 15 %. The suppression effect may vary distinctly between subjects. Moreover, the stimulus rate may be slightly changed during the test in order to reduce artefacts due to electromagnetic radiation.

In contrast to DPOAEs, ASSRs give information on the degree of the hearing loss within the entire dynamic range of hearing. Therefore, ASSRs are a suited means for frequency-specific and quantitative evaluation of hearing status especially in children and may contribute to improve hearing aid fitting. Intended use of ASSR is follow-up diagnostics after newborn hearing screening, topological diagnostics (together with tympanometry, OAE), neurological diagnostics, and frequency-specific determination of hearing thresholds.

ASSRs are especially useful for objective hearing threshold determination at low frequencies, where DPOAEs are not available. For objective determination of hearing thresholds it is recommended to measure at a first stage DPOAE Threshold at the required frequencies and to measure at a second stage ASSRs at frequencies where no DPOAE response could be detected. Average test time for eight test frequencies is about 15 minutes for normal hearing subjects and for hearing impaired subjects about 25 minutes (Rosner, 2013b).

Two ASSR workflows are available:

- ASSR Fixed allows frequency-specific measurement of ASSR at one or multiple levels. A response detection status (valid/invalid response) is delivered for each frequency/level combination.
- **ASSR Threshold** allows frequency-specific determination of hearing thresholds within a configurable level range by measurement of ASSR with adaptive level setting (bracketing algorithm). Hearing loss is estimated by means of fixed offset values. Threshold detection is stopped at a frequency if two consecutive responses at a frequency are not valid.

ASSR tests can be performed binaurally if an adequate transducer is connected (e.g. headphone, two ear probes), which may reduce test time by a factor of two. ASSR can be measured at frequencies from 250 Hz to 8 kHz with multi-frequency stimulus presentation.

PRACTICAL USE

Select *ASSR* from the module selection screen. If more than one AEP test is licensed, *ASSR* can be found in the *AEP* section. Select the preset that you would like to perform. If necessary, change the parameters (e.g. test mode, stimulus frequencies, stimulus levels (ASSR Fixed) or level range (ASSR Threshold), stimulus rate, stimulus bandwidth, averaging time) and the preset name as required.

Stimulus rate:

It is generally recommended to use 40 Hz ASSR in awake and alert subjects (e.g. children and adults) and to use 80 Hz ASSR in sleeping subjects (e.g. babies, sedated subjects, during anesthesia). The auto mode provides a performance-optimized stimulus rate paradigm that increases the stimulus rate with increasing frequency taking into account the dependency of stimulus rates on carrier frequencies and the dependency of noise floor levels on stimulus rates. The auto mode can be used for awake or sleeping subjects. Please note that stimulus rates are jittered for improving resistance to interfering noise.

Stimulus bandwidth:

With increasing stimulus bandwidth the response amplitude increases but frequency specificity decreases. The stimulus bandwidth selection impacts the frequency selection.

Masking noise

Contralateral masking is recommended if there is significant asymmetry in hearing loss between ears, i.e. differences of about 30 to 40 dB (headphones) or about 50 to 60 dB (insert earphones). If a bone conductor is used, application of contralateral noise is essential.

Averages:

At fixed measurement conditions, with increasing averaging time, the noise floor decreases and the reliability of the validity of responses and hence the quality of threshold determination increases.

Noise stop criterion:

If activated (i.e. noise stop criterion >0 nV), the recording of a trace is stopped as soon as the residual noise drops below the defined noise threshold and no response is detected. Hence, if activated this option speeds up the recording in case no response is present.

Before the test is started the subject should be instructed about the test procedure. In order to reduce muscle artefacts, the subject should be calm and fully relaxed lying comfortably on a recliner or bed. For 40 Hz ASSR subjects should stay awake, whereas for 80 Hz ASSR subjects may sleep. For babies, try to test the subject during spontaneous sleep or sedation if medically approved. In order to reduce environmental artefacts, conduct the measurement in a room with low electromagnetic radiation, i.e. in an electrically shielded metal booth or any other room without powered-on electric devices (e.g. computer, light, telephone, cell phone, power transformer) in close range to the measurement equipment. A shielded electrode cable may also help in reducing electromagnetic noise. An acoustically shielded booth or a quiet room is recommended if ASSR is applied for threshold determination at low sound pressure levels.

Make sure that a valid transducer (e.g. headphone, insert earphone, ear probe, bone conductor) and electrode cable are connected. Make sure that the skin is clean at the intended positions of the electrodes. For further instructions about skin preparation, electrode placement and impedance please refer to *PRACTICAL USE* in section 3.2.5: Auditory Brainstem Reponses (ABRs).

Select the test ear (Right+Left: simultaneous measurement of right and left ear – please note that in this case a symmetric electrode montage is recommended). The electrode impedance measurement starts. Both impedances must be <12 k Ω in order to be able to proceed with the measurement. It is recommended that impedances are <6 k Ω and that the difference between red and white electrode impedance is <3 k Ω . After successful electrode impedance measurement, you can start the test by pressing the play button. The ASSR stimuli are presented via the transducer and the response is detected via the electrode cable. During the measurement the electrode impedance is monitored. In case that at least one of the impedances gets too high (e.g. an electrode has fallen off), the test is automatically interrupted. The test can only be continued if the electrode impedances get back to tolerable values. For further information about impedance measurement see $PRACTICAL\ USE$ in section 3.2.5: Auditory Brainstem Reponses (ABRs).

During the measurement the following information is provided on the screen (see *Figure 60*):

The statistical validity of an ASSR is shown for each currently tested frequency/level combination ①. The traces are shown in red (right ear) and/or blue (left ear). For 40-Hz ASSR up to two and for 80-Hz ASSR and for automatic stimulus rate mode up to four frequency/level combinations are measured simultaneously. A trace is finished if it hits the top of the box (valid response), if it hits the gray triangular area within the box (invalid response) or if the configured maximum averaging time is reached (invalid response). The EEG noise bar ② represents biological and external noise which is an indicator for the adequacy of measurement conditions. With good measurement conditions, the EEG noise bar should remain green. With increasing noise, the EEG noise bar gradually turns yellow and red, which indicates that biological and external artefacts may be too high, i.e. the patient may be not calm and relaxed enough (e.g. swallowing, moving, clenched jaw/teeth), or the electromagnetic radiation may be too prominent. When pressing the EEG noise bar, an EEG noise graph ⑤ is shown. You can return to the EEG noise bar view by pressing the EEG noise graph. The overall progress bar

③ shows the progress of the measurement. If the bar is full, the complete test sequence is finished. The electrode impedance ④ is regularly updated during the measurement and the measurement is paused if impedance gets too high.

The test can be manually paused $\fill \cite{O}$, continued after a pause (pause button turns to play button – not shown in Figure 60), or stopped \cite{O} . Also, a current measurement can be skipped \cite{O} . If multiple measurements are running the frequency/level combination can be selected (skip \cite{O} : skip level at frequency, stop \cite{O} : skip all levels at frequency) after pressing the main forward button \cite{O} . The forward \cite{O} and stop buttons \cite{O} at the frequency/level combinations can be removed by pressing the main forward button \cite{O} again. The current measurement status for the right or left ear can be displayed when pressing the left or right graph area, respectively. The screen is colored red for the right ear and blue for the left ear. Currently tested frequency/level combinations are shown with the clock symbol \cite{O} , already tested frequency/level combinations feature a green check mark \cite{O} (valid response), a red \cite{O} (invalid response), or a gray \cite{O} (incomplete response because measurement was stopped for this frequency/level combination).




Figure 60: ASSR measurement (top: monaural 40 Hz ASSR measurement; top left: standard measurement screen; top right: measurement screen with EEG noise graph; bottom: binaural 80 Hz ASSR measurement; bottom left: standard measurement screen; bottom middle: skip/stop selection; bottom right: measurement status and result matrix)

The statistical response algorithm is based on the magnitude squared coherence test (Dobie and Wilson, 1989) and determines whether a response is present or absent. The response detection algorithm compares the mean square (= power) of the averaged waveform to the sum of mean

squares of the single sweeps and calculates a score. A valid response is detected if the detection score of 4 (which is equivalent to a significance level of 5%) is achieved after a minimum of 45 seconds has been averaged. This means the power of the grand average must be four times higher than the power added up of the single sweeps.

After the test is finished, the ASSR result is shown (see *Figure 61* for general ASSR result views and *Figure 62* for audiogram result view specific for ASSR Threshold). For all ASSR measurements, three screens are available. The response status matrix ① shows the response status for each tested frequency/level combination, i.e. if the response is valid (green check mark), invalid (red Ø) or undefined/not finished (gray Ø). When pressing any of the response status symbols, the related detailed response information ② is shown including information on frequency, level, stimulus rate and response status at the top of the screen. Below, the statistical validity trace is shown in red (right ear) or blue (left ear), respectively. In *Figure 61* you can see three typical response examples (bottom left: valid response — the validity trace hits the top of the box before the maximum measurement time is reached; bottom middle: invalid response — the validity trace hits the gray triangle within the box; bottom right: invalid response — the timeout is reached without getting a valid response). At the bottom, the response amplitude, noise amplitude, and electrode impedances are displayed. You can return to the response status matrix by pressing the screen or the *back* button. The general information screen ③ shows the applied transducer, the test configuration (test mode, stimulus bandwidth, stimulus rate), and the total measurement time.




Figure 61: ASSR result (top left: response status matrix; top right: general test information; bottom: detailed response information for valid and invalid responses)

In order to ensure adequate measurement conditions noise should be below 20 nV after 360 s averaging time. With increasing averaging time, noise decreases and vice versa, i.e. if a valid ASSR response is detected quickly, a higher noise is expected to be shown.

For ASSR Threshold an additional ASSR audiogram is available (see *Figure 62*). The audiogram shows estimated ASSR thresholds, which are derived by subtracting an empirically determined offset from the lowest stimulus level at which a valid ASSR response has been detected (Rosner, 2013a). Thresholds are shown in red (right ear) or blue color (left ear). The width of the threshold symbol corresponds to the configured stimulus bandwidth. An *arrow up* symbol at the threshold symbol means that the threshold is equal or better to the shown threshold. An *arrow down* symbol means that the threshold exceeds the level at which this symbol is shown. By sliding a finger from right to left on the audiogram screen, the response matrix is shown (see *Figure 61*).




Figure 62: ASSR Threshold result (top left: ½ octave stimulus bandwidth; top right: octave stimulus bandwidth; bottom: two octaves stimulus bandwidth).

The following case examples (see *Figure 63*) show that ASSRs are capable of appropriately estimating behavioral thresholds. The relationship between hearing thresholds estimated from ASSR measurements and behavioral hearing threshold has been investigated by various study groups (e.g. Dimitrijevic *et al.*, 2002, Ahn *et al.*, 2007).




Figure 63: ASSR case examples (left top/bottom: behavioral pure-tone thresholds; right top: ASSR thresholds; right bottom: ASSR response status matrix)

LITERATURE

- Ahn JH, Lee HS, Kim YJ, Yoon TH, Chung JW (2007): Comparing pure-tone audiometry and auditory steady state response for the measurement of hearing loss, Otolaryngology Head and Neck Surgery 136, p. 966-971
- Dimitrijevic A, John MS, van Roon P, Purcell DW, Adamonis J, Ostroff J, Nedzelski JM, Picton TW (2002): Estimating the audiogram using multiple auditory steady-state responses, Journal of the American Academy of Audiology 13, p. 205-224.
- Dobie RA, Wilson MJ (1989): Analysis of auditory evoked potentials by magnitude squared coherence, Ear and Hearing 10, p. 2-13.
- Elberling C, Don M, Cebulla M, Stürzebecher E (2007): Auditory steady-state responses to chirp stimuli based on cochlear traveling wave delay. Journal of the Acoustical Society of America 122, p. 2772-2785
- Galambos R, Makeig S, Talmachoff PJ (1981): A 40-Hz auditory potential recorded from the human scalp, Proceedings of the National Academy of Sciences of the United States of America 78, p. 2643-2647.
- Herdmann AT, Lins O, van Roon P, Stapells DR, Scherg M, Picton TW (2002): Intracerebral sources of human auditory steady-state responses, Brain Topography 15, p. 69-86.
- Herdman AT, Stapells DR (2003): Auditory stead state response thresholds of adults with sensorineural hearing impairment, International Journal of Audiology 42, p. 237-248.
- John MS, Lins OG, Boucher BL, Picton TW (1998): Multiple auditory steady-state responses (MASTER): stimulus and recording parameters, Audiology 37, p. 59-82.

- John MS, Dimitrijevic A, van Roon P, Picton TW (2001): Multiple auditory steady-state responses to AM and FM stimuli, Audiology and Neurootology 6, p. 12-27.
- Levi EC, Folsom RC, Dobie RA (1993): Amplitude-modulation following response (AMFR): effects of modulation rate, carrier frequency, age, and state, Hearing Research 68, p. 42-52.
- Lins OG, Picton TW (1995): Auditory steady-state responses to multiple simultaneous stimuli, Electroencephalography and Clinical Neurophysiology 96, p. 420-432.
- Lodwig A (2015): Multi-Stimulus suppression in DPOAE vs. ASSR, IERASG presentation.
- Lodwig A (2017): Multi-Rate Spread Spectrum ASSR, IERASG presentation.
- Lodwig A, Rosner T (2016): Frequenzspezifische ASSR mit multiplen, nicht-konstanten Reizraten (in German: frequency-specific ASSR with multiple, non-constant stimulus rates), AGERA presentation.
- Mäkelä JP, Hari R (1987): Evidence for cortical origin of the 40 Hz auditory evoked response in man, Electroencephalography and Clinical Neurophysiology 66, p. 539-546.
- Pantev C, Roberts LE, Elbert T, Ross B, Wienbruch C (1996): Tonotopic organization of the sources of human auditory steady-state responses, Hearing Research 101, p. 62-74.
- Picton TW, Skinner CR, Champagne SC, Kellett AJ, Maiste AC (1987): Potentials evoked by the sinusoidal modulation of the amplitude or frequency of a tone, Journal of the Acoustical Society of America 82, p. 165-178.
- Picton TW, Dimitrijevic A, John MS, van Roon P (2001): The use of phase in the detection of auditory steady-state responses, Clinical Neurophysiology 112, p. 1698-1711.
- Picton TW, John MS, Dimitrijevic A, Purcell D (2003): Human auditory steady-state responses, International Journal of Audiology 42, p. 177-219.
- Rodriguez R, Picton T, Linden D, Hamel G, Laframboise G (1986): Human auditory steady state responses: effects of intensity and frequency, Ear and Hearing 7, p. 300-313.
- Rosner T (2011): Zeitoptimierte Hörschwellenbestimmung mit atoakustischen Emissionen und auditorischen Potentialen (in German: Time-optimised threshold determination with otoacoustic emissions and auditory potentials), dissertation, Technische Universität München.
- Rosner T (2013a): Chirp-evozierte ASSR mit einem Handgerät (in German: Chirp evoked ASSR on a handheld device), AGERA annual meeting.
- Rosner T (2013b): Chirp evoked ASSR on a handheld device, IERASG annual meeting.
- Ross B, Draganova R, Picton TW, Pantev C (2003): Frequency specificity of 40-Hz auditory steady-state responses, Hearing Research 186, p. 57-68.
- Stapells DR, Makeig S, Galambos R (1987): Auditory steady-state responses: threshold prediction using phase coherence, Electroencephalography and Clinical Neurophysiology 67, p. 260-270.
- Van Maanen A, Stapells DR (2005): Comparison of multiple auditory steady-state responses (80 versus 40 Hz) and slow cortical potentials for threshold estimation in hearing-impaired adults, International Journal of Audiology 44, p. 613-624.

3.2.7 Electrocochleography (ECochG)

Electrocochleography (ECochG) records neural activity from the cochlea and vestibulocochlear nerve by means of a near-field electrode positioned in the ear canal or resting on the surface of the eardrum (extratympanic recording) or with a needle electrode that is perforating the eardrum to rest on the cochlear promontory (transtympanic recording). Due to the invasive nature of transtympanic recording, this measurement technique has only limited applications and extratympanic recording is usually favored. A typical electrode montage includes an ear canal electrode deeply inserted into the ear canal (or a tympanic membrane electrode resting on the eardrum) (recording electrode), a surface electrode attached to the contralateral earlobe or mastoid (reference electrode), and a surface electrode attached to the cheek (ground electrode). Impedances may be relatively high for ear canal electrodes which may also result in a poor balancing of electrode impedances and hence may yield less efficient common mode rejection and noise suppression. Moreover, the amplitude and the signal-to-noise ratio (SNR) of the recorded potentials depend on the insertion depth of the ear canal electrode, i.e. the deeper the insertion the higher the amplitudes and the lower the SNR. It is also possible to use a surface electrode attached to the ipsilateral earlobe or mastoid as recording electrode. However, this electrode montage usually suffers from highly reduced amplitudes and is therefore not recommended. In general, when using an extratympanic recording electrode, extended signal averaging is needed in order to yield adequate responses.

The following potentials can be recorded and can be roughly allocated to the following generation places: (1) cochlear microphonic (CM): outer hair cells, (2) summating potential (SP): basilar membrane, and (3) action potential (AP): vestibulocochlear nerve. The action potential including its main component N1 is commonly the most prominent peak in the recorded response and is usually dominated by activity of basal nerve fibers. N1 virtually corresponds to wave I from an ABR measurement. The summating potential and the action potential can be enhanced by summing up rarefaction and condensation responses. The cochlear microphonic resembles the waveform of the stimulus especially at low to medium stimulus levels and therefore can be easily confused with stimulus artefacts. Moreover, the cochlear microphonic can be covered by the stimulus artefact when using extratympanic electrodes. The cochlear microphonic is suppressed by stimulation with alternating polarity. The visibility of cochlear microphonic is improved by subtracting condensation from rarefaction responses.

As stimulus a click or tone burst is typically used. The response latency amounts to about 1 to 3 ms. Action potential latencies increase with decreasing stimulus level reflecting the decreasing contribution of basal cochlear regions. In contrast, cochlear microphonic and summating potential latencies are rather independent of stimulus level. The amplitudes of all components increase with increasing stimulus level. The summating potential exhibits a distinct saturation at high stimulus levels. Especially the action potential amplitude also decreases with increasing stimulus rate. Masking is usually not applied for ECochG measurements as the recorded potentials are near-field potentials, which are largest close to the recording electrode.

Potential clinical applications include the determination of hearing thresholds, identification and monitoring of Menière's disease or endolymphatic hydrops, and intraoperative monitoring. However, please note that hearing threshold determination by means of ECochG (especially when recorded with an extratympanic electrode) is not as reliable as hearing threshold determination by means of ABR.

PRACTICAL USE (module in preparation, not available for firmware 2.2)

Please note that the ECochG module is in preparation, so that only preliminary information on the practical use of the module and on user interfaces can be presented. The module is not available in Firmware 2.2.

Select *ECochG* from the module selection screen. If more than one AEP test is licensed, *ECochG* can be found in the *AEP* section. Select the preset that you would like to perform. If necessary, change the parameters (e.g. stimulus level, stimulus rate) and the preset name as required.

Before the test is started the subject should be instructed about the test procedure. In order to reduce muscle artefacts, the subject should be calm and fully relaxed lying comfortably on a recliner or bed. For babies, try to test the subject during spontaneous sleep or sedation if medically approved. In order to reduce environmental artefacts, conduct the measurement in a room with low electromagnetic radiation, i.e. in an electrically shielded metal booth or any other room without powered-on electric devices (e.g. computer, light, telephone, cell phone, power transformer) in close range to the measurement equipment. A shielded electrode cable may also help in reducing electromagnetic noise. An acoustically shielded booth or a quiet room is recommended if ECochG is applied for threshold determination at low sound pressure levels.

Make sure that a valid transducer (e.g. headphone, insert earphone, ear probe) and electrode cable are connected and that both cables are separated from each other in order to reduce crosstalk from the transducer. Make sure that the skin is clean at the intended positions of the surface electrodes (reference, ground). For further instructions about skin preparation, electrode placement and impedance please refer to *PRACTICAL USE* in section 3.2.5: Auditory Brainstem Reponses (ABRs). Select a suitable recording electrode and position the recording electrode in the ear canal. For further information on electrode montage and aspects regarding electrode selection please refer to the initial paragraph of the ECochG section.

Select the test ear. The electrode impedance measurement starts. Both impedances must be suitable for performing a measurement. For further information about impedance measurement see *PRACTICAL USE* in section 3.2.5: Auditory Brainstem Reponses (ABRs). After successful electrode impedance measurement, you can start the test by pressing the play button. The ECochG stimulus is presented via the transducer and the response is detected via the electrode cable. During the measurement the electrode impedance is monitored. In case that at least one of the impedances gets too high (e.g. an electrode has fallen off), the test is automatically interrupted. The test can only be continued if the electrode impedances get back to tolerable values.

During the measurement the following information is provided on the screen (see Figure 64):

ECochG traces are shown for the configured stimulus level and rate in three graphs: responses from rarefaction and condensation stimuli shown together in one plot ①, rarefaction and condensation responses added to yield a single trace ②, and rarefaction and condensation responses subtracted to yield a single trace ③. The tested trace is shown in red (right ear) or blue (left ear).

During an ECochG measurement, the EEG noise bar ④ represents biological and external noise which is an indicator for the adequacy of measurement conditions. With good measurement conditions, the EEG noise bar should remain green. With increasing noise, the EEG noise bar gradually turns yellow and red, which indicates that biological and external artefacts may be too

high, i.e. the patient may be not calm and relaxed enough (e.g. swallowing, moving, clenched jaw/teeth), or the electromagnetic radiation may be too prominent. When pressing the EEG noise bar, an EEG noise graph is shown. You can return to the EEG noise bar view by pressing the EEG noise graph. The progress bar 5 shows the progress of the measurement. If the bar is full, the measurement is finished, i.e. the number of averages that has been defined by the user is reached. The electrode impedance 6 is regularly updated during the measurement and the measurement is paused if impedance gets too high. The test can be manually paused 7 or stopped 8.




Figure 64: ECochG measurement (left) and result (middle, right)

After the test is finished, the ECochG result is shown (see *Figure 64*). The recorded ECochG traces are displayed again in three graphs (rarefaction and condensation, rarefaction plus condensation, rarefaction minus condensation) as explained above. Below the graph additional test information is displayed as e.g. the stimulus level, noise, number of averages, and stimulus rate. By swiping over the time axis you can adjust the resolution of the time axis. By sliding a finger from right to left on the result screen you can move forward to the general information result screen, which presents information about the used transducer and settings (averages, stimulus rate, and other test options). In order to ensure adequate measurement conditions noise should be <100 nV after 2000 averages. With increasing averages, noise decreases.

LITERATURE

- American Speech-Language-Hearing Association (1987): Short latency auditory evoked potentials, online: http://www.asha.org/policy/RP1987-00024/.
- Eggermont JJ (1974): Basic Principles For Electrocochleography, Acta Oto-Laryngologica, 77: Sup. 316, p. 7-16.
- Ferraro JA (2010): Electrocochleography: A Review of Recording Approaches, Clinical Applications, and New Findings in Adults and Children, Journal of the American Academy of Audiology 21, p. 145–152.
- Sohmer H, Feinmesser M (1973): Routine Use Of Electrocochleography (Cochlear Audiometry) On Human Subjects, Audiology, 12:3, p. 167-173.

3.2.8 Overview: Immittance Audiometry

The function of the middle ear is to minimize the loss of acoustic energy that occurs when sound is transferred from air in the outer ear canal (low density) to fluid in the inner ear (high density). Without the specific middle ear features, approximately 99.9 % of the sound energy would be reflected at the fluid of the inner ear due to the different densities. The middle ear helps to improve the energy transmission to the inner ear by increasing sound pressure and force (see *Figure 65*). The increase of sound pressure is simply due to the fact that the tympanic membrane area is seventeen times larger than the area of the footplate of the stapes which is the connecting link between middle and inner ear at the oval window. The increase of sound pressure becomes clear when looking at the physical equation, which defines pressure as force divided by area (p = F / A). So, with reduced area and same force the pressure increases. The increase of force is due to the different length of the malleus and the incus providing a lever action of the ossicular chain (malleus, incus, and stapes). Both mechanisms yield an impedance matching which allows for a transmission of up to 60 % of the sound energy to the inner ear at about 1 to 3 kHz.




Figure 65: Schematic drawing of the middle ear. A_1 and A_2 are the areas of the stapes footplate and tympanic membrane, I_1 und I_2 are the lever arm lengths of malleus and incus

The middle ear is able to increase the impedance if necessary for providing protection against loud sounds. In case of a sound exceeding about 80 dB HL, the stapedius muscle is activated resulting in an increased stiffness of the middle ear. As a consequence, the energy transmitted to the inner ear decreases. The stapedius reflex is active both on the ipsilateral and contralateral ear.

The term immittance comprises impedance Z (unit: acoustic ohm) or admittance Y (= 1/impedance, unit: acoustic mho) and their complex components. Regarding the middle ear, the impedance defines the opposition of the middle ear to the flow of acoustic energy, whereas the admittance defines the amount of energy that flows into the middle ear. In middle ear diagnostics, typically the admittance is evaluated. Admittance is determined by compliance (= 1/stiffness, spring load), mass, and friction or resistance. Mathematically, the admittance Y is a complex value consisting of conductance G (real part) and susceptance B (imaginary part), i.e. Y = G + jB. Friction influences conductance, whereas compliance and mass influence susceptance. Conductance (friction) is independent of frequency, whereas susceptance (compliance, mass) is dependent on frequency with compliant susceptance

being inversely proportional to frequency and mass susceptance being directly proportional to frequency. With increasing frequency, the total susceptance progresses from positive values (stiffness controlled) towards 0 mmho (resonance) to negative values (mass controlled). The resonance frequency is directly proportional to the stiffness of the middle ear, i.e. with increasing stiffness the resonance frequency increases (e.g. with otosclerosis), and inversely proportional to the mass of the middle ear, i.e. with increasing mass the resonance frequency decreases.

The typical measurement setup consists of an ear probe which is sealed in the ear canal and contains a loudspeaker that produces the stimulus signal, a microphone that monitors the response (i.e. the portion of the stimulus reflected from the tympanic membrane), and a pneumatic port for applying air pressure via an air tube. The device to which the ear probe is connected provides the pump for generating varying air pressure. The impedance calibration is typically done with several test cavities of different volumes. The test cavities have a defined volume of air so that their admittance is known. This allows calibrating the device for the typical range of acoustic admittances.

Immittance audiometry determines the reflectance of sound at the ear drum and is hence able to provide information on middle ear dysfunction caused by different middle ear pathologies including otosclerosis (increase of stiffness), otitis media (increase of mass), and Eustachian tube dysfunction (decrease of static pressure in the tympanic cavity). Different pathologies may result in similar tympanograms (e.g. flat tympanogram may occur due to malleus fixation or middle ear effusion), and the same pathology in different subjects or with different characteristics may not result in identical tympanogram shapes (e.g. tympanogram in middle ear effusion is dependent on the amount of fluid in the tympanic cavity) so that additional diagnostics may be necessary. Tympanometry and acoustic reflex tests cannot be reasonably performed if the ear drum is perforated (e.g. if a tympanostomy tube is in place). For more information on tympanometry, acoustic reflex tests, and Eustachian tube function tests please refer to sections 3.2.9: Tympanometry, 3.2.10: Acoustic Reflex Test, and 3.2.11: Eustachian Tube Function Tests, respectively.

PRACTICAL USE

It is recommended to check the functionality of the pneumatic system (i.e. air tube and pump unit) regularly. This can be done via the pump unit test (see Figure 66), which checks correct pressure generation and retention. Connect the tympanometry ear probe to the device and make sure that also the air tube is properly attached to the device. Place the probe tip in an appropriate test cavity and make sure that the probe tip is deeply inserted into the test cavity so that there is no leakage between probe tip and test cavity. First, pressure generation is tested. The static air pressure starts at 0 daPa and continuously increases until +200 daPa, i.e. a continuously increasing curve should be displayed. Second, pressure retention is tested. The initial static air pressure of +200 daPa shall be kept without major reduction of air pressure during the test time, i.e. a rather flat line should be displayed. If the pneumatic system functions appropriately, the test passes (i.e. a green icon is shown). Please note that if the second subtest (pressure retention) fails it may be due to a leaky connection, which primarily may occur between the ear probe and the test cavity. In this case try another test cavity or try to seal the ear probe by detaching the probe tip and placing a finger on top of the output channels. Alternatively, conduct the test for the pump unit only as follows: detach the tympanometry ear probe from the device and seal the air pressure socket at the device e.g. with a finger or by screwing the cap on the pressure outlet socket. In this setup only the pump unit of the device is tested.




Figure 66: Pump unit test (left: test with slightly leaky probe tip seal; middle: test with closed pressure outlet socket; right: failed test with leaky seal when probe tip is not inserted into the test cavity)

It is recommended to readjust the admittance calibration on the device regularly via the tympanometry calibration procedure (see *Figure 67*), which sets the correct admittance reference values for the three test cavities of known volume (0.5, 2.0, and 5.0 ml).




Figure 67: Admittance calibration at three different volumes

For this test, use the provided test cavity box. Place the correct probe tip (PT-S) in the first test cavity (e.g. 0.5 ml) and make sure that the probe tip is deeply inserted into the test cavity so that there is no leakage between probe tip and test cavity. Wait until the volume value settles. If the displayed

value differs from the nominal value press the *re-calibrate* button. Press the *volume* (0.5 ml) button to start the calibration for the 0.5 ml volume. The respective admittance is adjusted accordingly. Proceed with the two other volumes as indicated on the screen. During calibration for a volume a progress bar is shown. If all three test volumes have been calibrated a message box will appear indicating either successful or failed calibration. If calibration failed please repeat the procedure and make sure that you use the correct volumes as indicated on the screen. Please remember that the volume of the cavity next to the silver screw represents the 0.5 ml volume. After successful calibration you can verify the admittance calibration by re-entering the tympanometry calibration. Check that the displayed value corresponds to the respective nominal value.

If there is any deviation from the behavior described above or if you suspect any dysfunction, please retry with another tympanometry ear probe and/or contact your distributor.

LITERATURE

- American Speech-Language-Hearing Association (1988): Tympanometry, online: http://www.asha.org/policy/RP1988-00027.htm
- Shanks J, Shohet J (2009): Tympanometry in clinical practice, In: Handbook of clinical audiology (editors: Katz J, Medwetsky L, Burkard R, Hood L), Lippincott, Williams & Wilkins, p. 157-188.

3.2.9 Tympanometry

METHODOLOGY

Sound conduction through the middle ear (from tympanic membrane via ossicles to footplate of stapes) can be tested by measuring the acoustic admittance (= 1/impedance) of the tympanic membrane while varying the static pressure in the outer ear canal using an air pump. The function of the tympanic membrane is best if pressure is equal on both sides, i.e. in the tympanic cavity and in the outer ear canal. In this state, the vibration of the tympanic membrane is maximal and the largest part of a sound signal is transmitted to the inner ear and only a small part is reflected. In case of differing pressure on both sides of the tympanic membrane, more energy is reflected. This situation may occur e.g. due to a ventilation disorder of the tympanic cavity, fluid in the tympanic cavity or perforation of the tympanic membrane. In general, it must be taken into account that impedance tests only refer to the acoustic energy that flows into the middle ear but not how much energy is transmitted through the system. There are specific pathologies (e.g. tympanosclerotic plaques) that may result in an abnormal tympanogram but are not associated with a distinct hearing loss. In contrast, otosclerosis may produce a relatively normal tympanogram with a marked hearing loss.

Due to the measurement position, which is remote to the tympanic membrane, the admittance measurement is influenced by the middle ear *and* outer ear canal properties and hence also by the insertion depth of the ear probe. In order to derive a measure of middle ear admittance alone, the ear canal admittance must be removed from the overall admittance. Under sufficient tension of the tympanic membrane by a high positive or negative pressure, the middle ear admittance is driven towards zero, i.e. the admittance at these conditions provides an estimate for ear canal admittance. The ear canal admittance can be subtracted from the overall admittance due to the assumption that the impedances of the ear canal and middle ear are configured in parallel (Djupesland and Zwislocki, 1972), which means that the admittances are summed up $(Y_{total} = Y_{ear\ canal} + Y_{middle\ ear})$. However, this assumption is only an approximation, so that the ear canal admittance derived from tympanometric measurements is an estimate. In general, this procedure is considered adequate for clinical use.

Standard clinical tympanometry is usually performed with a low frequency stimulus at 226 Hz (or 220 Hz). At low frequencies, the normal middle ear system is stiffness-controlled and susceptance (stiffness element) contributes more to overall admittance than conductance (frictional element). Other typical probe tone frequencies include 678 Hz (or 630, 660 Hz), 800 Hz, and 1000 Hz.

Typically static air pressure is varied from +300 daPa to -300 daPa. The direction of pressure change (i.e. from positive to negative pressure or vice versa) may influence static admittance (Wilson *et al.*, 1984). At higher frequency probe tones (e.g. 678 Hz) notched tympanograms are more frequent with increasing pressure change (Wilson *et al.*, 1984). Also, the rate of pressure change can have an effect on tympanograms. Single-peaked tympanograms typically increase in amplitude with increasing rates of pressure change, but also the incidence of multiple-peaked tympanograms increases (Creten and van Camp, 1974). Moreover, the incidence of notches increases with successive runs of tympanometry measurements maybe due to the viscoelasticity of the tympanic membrane.

The result of the admittance measurement is a graph called a tympanogram which plots middle ear admittance as a function of static air pressure in the outer ear canal (see *Figure 68*). Different middle ear pathologies exhibit different tympanogram shapes. The following rough description refers to low frequency (226 Hz) tympanograms. In case of normal middle ear function the tympanogram shape

corresponds to a Gaussian bell curve with its maximum being around zero static pressure (black solid line), i.e. maximum energy is transferred into the middle ear at atmospheric pressure without any static pressure offset. If there is a Eustachian tube dysfunction the peak of the Gaussian bell curve is shifted in the direction of negative pressure values (red dashed line). This is due to the fact that the tympanic membrane moves best in its normal position, i.e. when the static pressure in the ear canal and the static pressure in the tympanic cavity are the same. If the static pressure in the tympanic cavity is negative then the static pressure in the outer ear canal has to be negative with the same value. As a result, the peak of the Gaussian bell curve is present exactly at the static pressure which is obtained in the tympanic cavity. In case of middle ear effusion, middle ear mass is increased. In this case middle ear movement is considerably reduced resulting in a lower compliance (light blue dotted line), which is nearly independent of static pressure. Also, in case of otosclerosis middle ear movement is reduced. As a consequence, the peak of the Gaussian bell curve is small, however located within the zero static pressure range (gray solid curve). An activated stapedius muscle yields reduced compliance as well.




Figure 68: Schematic overview of tympanometry (top) with tympanogram examples (bottom: black solid line: normal; gray solid line: otosclerosis; red dashed line: Eustachian tube dysfunction; light blue dotted line: effusion)

Typical parameters from a tympanogram are the tympanometric shape, ear canal volume, static admittance, tympanometric peak pressure, and tympanometric width. Jerger and Northern (1980) introduced three types of tympanogram shapes, which refer to low frequency (226 Hz) tympanograms. Type A represents a normal tympanogram with a pronounced peak around 0 daPa, Type B shows a flat tympanogram without pronounced peak, and Type C refers to tympanograms with the peak shifted to negative static pressure. For low probe tone frequencies commonly a single-

peaked tympanogram occurs. However, in neonates and for higher probe tone frequencies tympanograms often exhibit multiple peaks and notches (dependent on which immittance component is measured).

Ear canal volume is typically estimated as the admittance at the negative or positive maximum pressure (e.g. at +200 daPa). For low frequency probe tones only a small error occurs due to a phase difference between the admittance vector of middle ear and ear canal. At higher frequencies, this error becomes more prominent. For a 226 Hz probe tone, the ear canal volume is commonly given in ml (which is similar to mmho at this frequency and at specific environmental conditions). For higher frequency probe tones, the ear canal volume is given in mmho. The tympanogram can be normalized by subtracting the ear canal volume from the curve yielding the compensated static acoustic admittance Y_{tm} , which is an estimate of the acoustic admittance at the lateral surface of the tympanic membrane. It is typically higher when maximum negative (rather than maximum positive) pressure is used to estimate ear canal admittance (Margolis and Smith, 1977). The tympanometric peak pressure is the pressure at which the maximum admittance occurs. The tympanometric width is the pressure difference at one half of the compensated static acoustic admittance. It quantifies the relative sharpness of the peak.

PRACTICAL USE

Select *Tymp* from the module selection screen. The *Tymp* module can be found in the *Immittance* section. Select the preset that you would like to perform. If necessary, change the parameters (e.g. minimum/maximum pressure, speed of pressure change, probe tone frequency, auto start, auto stop, reflex options) and the preset name as required. The different parameter options and possible applications are explained in the following:

Pressure setup:

Configures pressure range and pump speed. Pressure range is up to -600 to +400 daPa for class 1, and -300 to +300 for class 2. Pump unit speed can be set in increments of 50 daPa/s to values of 50 to 200 daPa/s. The "afap" option changes pressure as fast as possible with 600 daPa/s at maximum reducing the speed to 200 daPa/s as soon as a steep rise of the tympanometric curve is reached.

Probe tone frequency:

The default value is 226 Hz. Class 1 offers 678, 800, and 1000 Hz probe tones in addition. 226 Hz is commonly used as the tympanogram shows the distinct patterns which allow a classification of tympanograms (type A, B, C). 1000 Hz is frequently used for small children. Please note that the multi-frequency (MF) option allows testing of all four frequencies simultaneously without influencing test results or increasing test time.

Auto Start:

The tympanogram recording is started as soon as an appropriate ear probe fit is detected. If the option is disabled, the test must be started manually by pressing the *play* button.

Auto Stop:

The tympanogram recording can be stopped as soon as the tympanometric peak is detected, which means that the pressure will not be decreased any further than necessary.

Manual pressure control:

The tympanogram can be recorded manually, i.e. the pressure can be set by the examiner.

The tympanometry test can optionally be performed with consecutive automated acoustic reflex testing. If you wish to perform an acoustic reflex test immediately after tympanometry select reflex options in the tympanometry settings accordingly (see 3.2.10: Acoustic Reflex Test).

Make sure that a valid transducer (tympanometry ear probe EP-TY) is connected to the device and check that the air tube connector is tightly screwed on the pressure outlet socket of the device (see *Figure 69*). You may use a clip to fix the ear probe cable to the subject's clothes. Select an ear tip with appropriate size matching the probe tip size and the subject's ear canal size. Make sure that the ear probe is inserted without any leakage between ear probe and ear canal. Foam tips are usually not suitable for performing tympanometry because they are not airtight.



Figure 69: Tympanometry setup (left: attaching the tympanometry ear probe EP-TY to a Sentiero Desktop device; top right: EP-TY status light, bottom right: fixation clip)

After selecting the test ear, the measurement is ready to be started. The ear probe status light indicates the current measurement condition:

- Steady light: Ready for testing please place the ear probe in the ear
- Slowly blinking (heartbeat): Measurement in progress
- Fast blinking: Leakage, unable to generate the required pressure in the ear canal

The device continuously monitors the probe fit, i.e. the acoustical volume seen at the ear probe (see *Figure 70*).




Figure 70: Monitoring of probe fit for tympanometry (left: invalid probe fit when starting the test; middle: valid probe fit when starting the test; right: invalid probe fit during test)

The tympanometry measurement starts automatically (if configured) once the ear probe is appropriately positioned in the ear canal and the acoustical volume has stabilized. The ear probe status light indicator will change from steady to heartbeat mode while the admittance curve is recorded. In case of leakage during the test, a message box appears asking the user to commence or cancel the test. The ear probe status light will change to fast blinking, and instructions will appear on the device display. The device will periodically retry to continue the measurement automatically.

During the tympanometry measurement the following information is provided on the screen (see *Figure 71*):




Figure 71: Tympanometry measurement (left: before measurement; middle: during measurement with cartoon mode enabled; right: manual measurement)

Before a test is started manually (if auto start option is disabled), the estimated ear canal volume (ECV) is continuously displayed as bar on the left side of the tympanogram graph ①. The numeric ECV value ④ is additionally shown below the tympanogram graph. For class 1 tympanometers, the tympanogram view mode (admittance Y, conductance G, and susceptance B) can be switched by pressing the YBG button ②. As soon as the tympanometry measurement starts (either automatically if the auto start option is enabled or after pressing the play button ⑤), the admittance as a function of the applied static pressure offset range is plotted as tympanogram ③. If the cartoon mode is enabled, a plane moves along the tympanogram trace. The cartoon mode is especially meant for focusing a child's interest and hence improving noise conditions during a test. If a peak is detected, the peak amplitude and static pressure offset is displayed below the ECV estimate ④. A gray box in the middle of the tympanogram graph represents the normal area ⑨. If the tympanogram peak is located within this area, the tympanogram can be expected to be normal. An ongoing measurement can be stopped by pressing the stop button ⑧.

After a tympanometry measurement is finished, there are several options: In case of artefacts or unclear results, the measurement can be restarted by pressing the *play* button (5), i.e. the current trace is discarded and replaced by a newly recorded trace. For class 1 tympanometers, up to three tympanometry traces can be recorded in a single measurement session. An additional tympanogram, which will be displayed on top of the existing one(s), can be added by pressing the *plus* button (6). This may be useful if tympanograms are to be compared under different conditions, e.g. before and after performing a Valsalva maneuver. The other ear can be tested with the same parameter settings by pressing the *ear* button (7). This is the quickest way of recording tympanograms for both ears. Since the tympanometry measurement starts automatically once the probe fit is stable (if auto start

is activated), remove the ear probe from the current test ear before pressing the *ear* button, to prevent the same ear being measured twice (once as left and once as right ear).

A tympanogram can also be recorded manually if configured in the parameter settings. The static pressure offset can be set by pressing the *left* or *right arrow* buttons or by using the slider ①.

Please note that if you have configured the tympanometry to be followed by an automatic acoustic reflex measurement, the acoustic reflex test starts after a valid tympanogram has been recorded. For more information on the acoustic reflex test please refer to section 3.2.10: Acoustic Reflex Test.

After the test is finished, the tympanometry result is shown (see Figure 72).




Figure 72: Tympanometry result (left: Admittance with subtracted ear canal volume; middle, right: admittance and conductance for 1000 Hz probe tone)




Figure 73: Tympanometry 3D result view (not available in firmware 2.2)

The tympanogram is shown as red (right ear) or blue trace (left ear). The resolution of the ml/mmho axis can be adjusted by swiping a finger up or down the screen. A vertical bar at the left side of the tympanogram graph represents the ear canal volume (ECV) ①. The numeric value of the ECV is shown below the tympanogram graph together with the tympanometric width (TW, for 226 Hz probe tones only), and the peak amplitude and static pressure offset ④. The normative area is shown as gray box ②. For class 1 tympanometers, admittance, susceptance and conductance graphs can be toggled via the *YBG* button ③. For 226 Hz probe tones, the tympanogram can be shown normalized

by subtracting the estimated ear canal volume from the trace yielding the compensated static acoustic admittance Y_{tm} (= admittance of the tympanic membrane alone). The probe tone frequency and view mode are indicated in the upper right corner. For multi-frequency measurements you can chose to plot the tympanogram for each probe tone frequency by pressing the respective small admittance graph 6. When pressing the 3D tympanogram graph 5 (not available in firmware 2.2), a 3D plot calculated from the traces at the four probe tone frequencies is shown (see *Figure 73*). You can rotate the 3D plot by sliding a finger horizontally or vertically over the screen. By pressing the *magnifying glass* buttons you can zoom in or out of the graph.

LITERATURE

- American Speech-Language-Hearing Association (1988): Tympanometry, online: http://www.asha.org/policy/RP1988-00027.htm
- Creten W, van Camp K (1974): Transient and quasistatic tympanometry, Scandinavian Audiology 3, p.39-42.
- Djupesland G, Zwislocki J (1972): Sound pressure distribution the outer ear, Scandinavian Audiology 1, p. 197-203.
- Jerger J, Northern J (1980): Clinical Impedance Audiometry, Thieme.
- Kandzia F (2015): An optimized protocol for acoustic impedance measurements: "Simultaneous multi component multi frequency tympanometry", DGA poster presentation.
- Margolis R, Smith P (1977): Tympanometric asymmetry, Journal of Speech and Hearing Research 20, p. 437-446.
- Wilson R, Shanks J, Kaplan S (1984): Tympanometric changes at 226 Hz and 678 Hz across ten trials and for two directions of ear canal pressure change, Journal of Speech and Hearing Research 27, p. 257–266.

3.2.10 Acoustic Reflex Test

METHODOLOGY

The acoustic reflex (or stapedius reflex, attenuation reflex, auditory reflex) is an involuntary contraction of the stapedius muscle in response to high-level sound stimuli (above about 80 dB HL). The stapedius muscle stiffens the ossicular chain by pulling the stapes away from the oval window of the cochlea. As a consequence the transmission of vibrational energy to the cochlea is decreased and the impedance is increased. The pathway involved in the acoustic reflex is complex and can involve the ossicular chain itself, the cochlea, the auditory nerve and the brainstem. The stapedius reflex is considered to be a protective procedure. In cases where this protective procedure of the middle ear muscles is disabled (dysfunction of the efferent hearing system, otosclerosis) cochlear micromechanics (especially hair cell stereocilia) may be more easily damaged.

The acoustic reflex can be elicited ipsilaterally or contralaterally with short tone pulses. When elicited ipsilaterally, impedance measurement and acoustic reflex stimulation are done in the same ear. When elicited contralaterally, impedance measurement and acoustic reflex stimulation are done in opposite ears. Commonly, the threshold at which the stapedius muscle begins to change middle ear impedance is determined. The test is conducted at maximum compliance, i.e. at the static pressure offset, which yields a peak during tympanometry. Therefore, it is recommended to conduct tympanometry before acoustic reflex testing. The result is a graph which plots middle ear admittance as a function of time for a defined stimulus level. The duration of the stimulus may vary from 2 s (determination of acoustic reflex) to about 10 s (determination of decay of acoustic reflex over time). If the stimulus level is high enough to elicit an acoustic reflex, a distinct change in admittance can be observed.

Acoustic reflex measurements are done to identify otosclerosis, cochlear recruitment and neural disorders. Otosclerosis means fixation of the footplate within the oval window of the inner ear due to sclerotic lesions. In case of a fixation of the stapes footplate no change in impedance is measurable when delivering a loud sound to the ear. Recruitment is characterized by normal cochlear hearing at high and limited cochlear hearing at low sound levels. In hearing-impaired ears (e.g. revealed by pure-tone audiometry) having normal acoustic reflex thresholds at about 80 dB HL cochlear recruitment is most likely. The acoustic reflex threshold increases in case of middle ear dysfunction.

PRACTICAL USE

Automatic acoustic reflex

Before performing acoustic reflex testing, a tympanometry test should be performed in order to determine the tympanometric peak pressure at which reflex testing should be conducted. This sequence can be run automatically by selecting or setting up an appropriate preset in the *Tymp* module in the *Immittance* modules section. The following options referring to automatic acoustic reflex testing are available in the tympanometry settings:

Acoustic reflex mode:

Select "always" if you would like to start automated acoustic reflex testing irrespective of the tympanometry result or select "if peak within norm (226 Hz only)" if you would like to start

acoustic reflex testing only if the tympanogram peak is within the normal area. If you do not require acoustic reflex testing select "never".

Acoustic reflex frequencies and start/stop level:

Select the frequencies at which you would like to perform acoustic reflex testing. If you wish to perform a screening test at a specific level, set the start and stop level to the same value. If you wish to determine an acoustic reflex threshold, select the start and stop level for the level range to be checked. The maximum level during automatic mode is limited to 100 dB HL.

Acoustic reflex validity criterion and repeat to confirm:

Select the criterion for a change in admittance which is considered a valid acoustic reflex. If a valid acoustic reflex is determined you can chose to repeat the test at the respective stimulus configuration for confirmation.

Acoustic reflex pressure setup:

Select the static pressure offset at which the acoustic reflex testing shall be performed. You can either use the pressure at the peak of the previously measured tympanogram or you can specify a fixed static pressure offset. If a stable airtight seal cannot be achieved even with an ear tip of suitable size, the measurement can be performed at ambient pressure (i.e. 0 daPa). Although it is recommended to perform acoustic reflex testing at the peak tympanometric pressure, the test at ambient pressure may provide useful information, especially if the tympanometric peak pressure is close to zero anyways, i.e. when the middle ear function is normal.

Make sure that a valid transducer (tympanometry ear probe EP-TY) is connected. The acoustic reflex measurement is quite sensitive to movement of the probe. This is because the change in admittance which is caused by the acoustic reflex is quite small. Movements of the ear probe can influence the insertion depth and with it, the admittance that is seen by the probe. These changes are in the same order as that of the reflexes. It is therefore recommended not to hold the ear probe during acoustic reflex testing. You may use a clip to fix the ear probe cable to the subject's clothes. Select an ear tip with appropriate size matching the probe tip size and the subject's ear canal size. Make sure that the ear probe is inserted without any leakage between ear probe and ear canal. Foam tips are not suitable for performing tympanometry or acoustic reflex testing because they are not airtight. An adequate ear tip is necessary that allows a stable probe fit.

Contralateral acoustic reflexes can be measured if a suitable transducer for contralateral stimulation is connected in addition to the tympanometry ear probe. The transducer should be connected prior to selecting the acoustic reflex module in order for it to be recognized by the device. Suitable contralateral transducers include any headphones (HP-xx), insert ear phones (IP-xx), or the EP-VIP ear probe. Please note that there are monaural transducers which are specifically meant to be applied for contralateral acoustic reflex testing.

In automatic mode, the static pressure offset for performing acoustic reflex testing is taken over from the latest tympanometry result of the same ear (or depending on the parameter configuration a fixed static pressure offset is used). Acoustic reflex testing will start automatically without further user interaction required. Acoustic reflexes can be determined at various sine frequencies between 0.5 and 4 kHz and with noise (broadband, high-pass, low-pass). Please note that the deflection of the acoustic reflex traces (positive or negative) can be configured in *Tymp Preferences*.

During the automatic acoustic reflex measurement the following information is provided on the screen (see *Figure 74*):




Figure 74: Automatic acoustic reflex measurement (left: threshold determination with enabled cartoon mode; middle: manual mode after test is finished; right: current result overview)

The currently measured acoustic reflex trace is shown in the main graph area (1). Previously measured reflex traces for the same stimulus frequency are displayed in the miniature graph areas (2) below the main graph area. The reflex growth at the stimulus frequency under test can be deduced from the traces in the miniature graph areas if multiple levels are tested. A black marker indicates the miniature graph area in which the currently active trace will be plotted. If the cartoon mode is enabled, a flying object (e.g. a plane) is moving along the trace. The type of flying object changes from trace to trace. The cartoon mode is especially meant for focusing a child's interest and hence improving noise conditions during a test. Controls are grayed out and disabled while the automatic measurement is ongoing. The automatic measurement at a stimulus frequency is finished if the highest configured level is reached or a valid acoustic reflex is determined. The automated test can also be interrupted by the user at any time by pressing the stop button while the test is running (5). As soon as the automatic test is stopped by the user or finished at all selected stimulus frequencies the controls are enabled again and the user can replace or add additional traces in manual mode. Frequency, level and pressure controls will then be available for manually performing measurements, e.g. if a certain recording has to be repeated because of artefacts or unclear results. For further information on the practical use of the manual mode please refer to the section on manual acoustic reflex below.

An overview that shows all current measurements with the respective status (green check mark in case of a valid reflex) ④ can be displayed by pressing the *magnifying glass* button ⑥. If you wish to trigger the test sequence including tympanometry and automatic acoustic reflex testing on the other ear, please press the *ear* button ⑩. This is the quickest way of recording tympanograms and acoustic reflexes for both ears. Since the tympanometry measurement starts automatically once the probe fit is stable (if auto start is activated), remove the ear probe from the current test ear before pressing the *ear* button, to prevent the same ear being measured twice (once as left and once as right ear).

The acoustic reflex result view (see *Figure 75*) shows an overview of all ipsilateral or contralateral acoustic reflex traces. By selecting any of the miniature graphs the respective detailed view is

opened. If the trace is valid a green check mark is shown together with the trace. The change in admittance due to the activation of the acoustic reflex is shown at a valid trace. Below the graph the stimulus time frame (stimulus turned on from 1 to 3 s after the start of the measurement) is plotted.




Figure 75: Acoustic reflex result (left: result overview; right: detailed information)

Manual acoustic reflex

If licensed as class 1 tympanometer, *Manual Reflex* is available in the *Immittance* section as a separate measurement mode, so that there is no need to go through a sequence of tympanometry and automatic reflex threshold in advance. For all devices, the manual mode is available after finishing or interrupting an automated acoustic reflex test. Please note that this section focuses on the specific options available for manual testing. The general workflow and user interface is similar to automatic acoustic reflex testing explained above.

During the manual acoustic reflex measurement the following information is provided on the screen (see *Figure 76*).




Figure 76: Manual acoustic reflex measurement (left: measurement; right: current result overview)

As for automatic acoustic reflex testing, the current acoustic reflex trace ① and previously recorded traces for the same stimulus type ② are shown. Valid traces are marked with a green check mark. The manual mode allows full control of measurement parameters and measurement order. Measurement parameters can be setup for each individual reflex measurement. The parameters are

the stimulus frequency or noise (including additional high-pass or low-pass noise) ⑤, the stimulus level ③ (stimulus levels exceeding 100 dB HL are available depending on frequency), and the presentation ear side (ipsilateral or contralateral if contralateral transducer is connected) ④. The stimulus level step size can be scaled down to 1 dB by activating the respective checkbox ①. The static pressure offset can be adjusted by pressing the *static pressure offset* button ⑥. Once all parameters are set as requested, press the *play* button ⑧ in order to perform the test. If you wish to continue testing on the other ear press the *ear* ⑩ button. The acoustic reflex testing can be finished by pressing the *stop* button ⑦.

The device offers four storage slots per stimulus frequency and presentation modality (ipsilateral/contralateral). The slots are depicted as miniature traces and are annotated with the presentation level. A small black marker at the graph corner indicates the storage slot of the next measurement to be performed. Clicking on any of the slots' graph will make the corresponding curve appear in the main graph and will select that slot to be overwritten with the next measurement. An overview of all recorded traces can be shown by pressing the *magnifying glass* button (9). The overview (12) shows miniature graphs for all stimulus frequencies, whereas in the main view, only the reflex growth of the currently selected stimulus frequency is visible. Please note that stimulus levels are shown for each stimulus frequency in the sequence chosen by the examiner.

Acoustic reflex decay

If licensed as class 1 tympanometer, acoustic reflex decay testing is possible. An acoustic reflex decay test examines if the change in admittance due to the acoustic reflex can be maintained during continuous ipsilateral or contralateral stimulation. The acoustic reflex decay is often determined with contralateral stimulation in order to reduce the impact of potential stimulus artefacts on the ipsilateral recording. An acoustic reflex decay test is considered positive if the amplitude of the acoustic reflex decreases to 50 % within a given time period (usually around 10 seconds). Lower frequencies (500 Hz or 1000 Hz) are most commonly used for this test because also normal ears may exhibit an acoustic reflex decay at higher frequencies.

The operation of the acoustic reflex decay module is principally the same as for the manual acoustic reflex module, as the test only differs with regard to the stimulus duration and evaluation of recorded traces (see *Figure 77*).




Figure 77: Acoustic reflex decay measurement (left: during a measurement; middle: after a measurement with negative decay; right: after a measurement with positive decay)

Acoustic reflex decay testing is usually performed at stimulus levels 10 dB above the acoustic reflex threshold. Therefore, an acoustic reflex threshold determination, either manually or automatically, should be performed in advance. The static pressure offset should be set to the tympanometric peak pressure, as with normal acoustic reflex testing. Please be aware that intense stimulus levels may be uncomfortable or even painful to the patient. Good clinical judgement should be used before applying stimuli at high intensities, especially during acoustic reflex decay testing, where the actual duration of the stimulus presentation amounts to 14 seconds.

If an acoustic reflex decay test is positive, i.e. the trace reaches half the amplitude of the initial acoustic reflex amplitude before the end of the measurement, the time at which this amplitude is reached is shown as *half life* at the trace. The gray line at the trace represents 50 % amplitude of the acoustic reflex.

The acoustic reflex decay result view (see *Figure 78*) is similar to the acoustic reflex module and shows an overview of all ipsilateral or contralateral acoustic reflex decay traces. By selecting any of the miniature graphs the respective detailed view is opened. Below the graph the stimulus time frame (stimulus turned on from 1 to 15 s after the start of the measurement) is plotted.




Figure 78: Acoustic reflex decay result (left: result overview; right: detailed information)

Please note that admittance changes can be monitored using the patulous Eustachian tube subtest (see 3.2.11: Eustachian Tube Function Tests). Admittance changes caused by acoustic reflexes that were elicited by external stimulation, like free field sound sources, or hearing aid / cochlear implant based stimulation, can be recorded. It is up to the user to interpret the recorded curves, since the device has no knowledge about the external stimulation.

LITERATURE

- American Speech-Language-Hearing Association (1988): Tympanometry, online: http://www.asha.org/policy/RP1988-00027.htm
- Chermak GD, Dengerink JE, Dengerink HA (1983): Test-retest reliability of the acoustic stapedius reflex, Audiology 22, p. 136-143
- Deutsch LJ (1972): The Threshold of the Stapedius Reflex for Pure-Tone and Noise Stimuli, Acta Otolaryngology 74

3.2.11 Eustachian Tube Function Tests

METHODOLOGY

In addition to performing tympanometry and acoustic reflex testing, workflows are available to assess the Eustachian tube function (ETF) on a class 1 tympanometer. ETF test are used to check the operability of the Eustachian tube under different conditions. The Eustachian tube provides pressure equalization between the middle ear and the nasopharynx and with that to the external air pressure. The pressure equalization can be triggered by swallowing or yawning. Pressure differences cause a temporary conductive hearing loss due to reduced flexibility of the tympanic membrane and the ossicles. Certain middle ear diseases as e.g. otitis media may affect Eustachian tube function.

PRACTICAL USE

Select *ETF* and the desired subtest from the *Immittance* section of the measurement module selection menu. All ETF tests require the tympanometry ear probe EP-TY. The handling of the tympanometry ear probe is described in section *3.2.9: Tympanometry*. In the following the practical use of the different ETF subtests is explained.

Non-Perforated Eardrum

This Eustachian tube function test can be used if the tympanic membrane is not perforated. During the test three tympanograms are recorded. In between the recordings, the patient is instructed to swallow with nose and mouth closed (Toynbee maneuver) and to breathe out of the nose with nose and mouth closed (Valsalva maneuver). Given a normal Eustachian tube function, the tympanogram peaks are expected to shift between the recordings as a result of the middle ear pressure change caused by performing the Toynbee and Valsalva maneuver. Please note that the difference in static pressure offset between the three tympanogram peaks also depends on the ability of the patient to perform the maneuvers. Information on how to instruct the patient will be displayed on screen during the progress of the measurement (see *Figure 79*).




Figure 79: Non-perforated eardrum Eustachian tube function measurement (left: first tympanometry measurement without any intervention; middle: message box for second tympanometry measurement with Toynbee maneuver; right: message box for third tympanometry measurement with Valsalva maneuver)

The result of the non-perforated eardrum ETF test is plotted showing all three tympanograms (see *Figure 80*), i.e. the standard tympanogram and the two tympanograms with the Toynbee and Valsalva maneuver, respectively.




Figure 80: Non-perforated eardrum Eustachian tube function result

Perforated Eardrum

This Eustachian tube function test can be used if the tympanic membrane is perforated. This test pressurizes the ear canal and the middle ear if the eardrum is perforated. To test the Eustachian tube function, the patient should be instructed to swallow multiple times during the test. If the Eustachian tube opens, pressure will drop (see *Figure 81*). Special care must be taken to ensure an airtight fit of the ear probe during testing in order to make sure that the observed pressure drop is due to the Eustachian tube function and not due to leakage.




Figure 81: Perforated ear drum Eustachian tube function result

Patulous Eustachian Tube

This Eustachian tube function test can be used for patients with a patulous Eustachian tube, i.e. the Eustachian tube stays permanently or intermittently open. The test performs a high resolution admittance measurement, similar to an acoustic reflex measurement, but without stimulus. The admittance change is recorded and plotted for 20 seconds. It can be used to monitor admittance changes caused by Eustachian tube opening during breathing, chewing, or swallowing. If the

Eustachian tube is permanently open, a higher fluctuation of admittance is expected compared to normal middle ear status. It is common that the rhythms of heartbeat or breathing can be identified in the admittance curve in these cases.




Figure 82: Patulous Eustachian tube result with the heartbeat rhythm of a subject visible in the admittance curve

LITERATURE

- Smith ME, Tysome JR (2015): Tests of Eustachian tube function: a review, Clinical Otolaryngology 40, p. 300-311

Made in Germany

PATH MEDICAL GmbH Landsberger Straße 65 82110 Germering Germany

Tel.: +49 89 800 76 502 Fax: +49 89 800 76 503 Internet: www.pathme.de

