Path Medical Sentiero Quick Screen Instructions

To Screen – 3 button pushes!

- 1. Press "On" button to turn on
- 2. Press the top yellow "DPOAE Quick Screen" on the touch pad
- 3. Press the right or left options and the test will start



Other Quick functions:

- To retest from the results screen simply touch the blue "retest" option located below the results
- To screen the other ear press the green arrow (faces left) at the bottom of the screen. This will return you to the "right" or "left " ear options. Touch the option for the ear you want to screen and the test will begin

Power

- 2 hours to charge it takes 2 hours to fully charge the re-chargeable battery
- 6-8 hours of testing the screening unit will screen for 6-8 hours of testing before it needs re-charged
- You can test with the unit plugged into a power outlet while it re-charges

PE Tubes -

 You can test children with known PE tubes and should expect a pass if the ears are clear and PE tubes are functioning.

Maintenance

- Annual calibration Complete annually to ensure accuracy and to obtain software updates etc.
- Probe tips should be replaced when they become clogged. Do not attempt to clean them - they are disposable and must be replaced when they become clogged.
- Cleaning:
 - The device must be turned off and removed from all connected 1. components to be cleaned
 - May be wiped clean with a damp cloth using a mild antiseptic 2. solution (e.g. cetylcide). Do not clean with alcohol or other solvent.
 - 3. Do not put excessive pressure on the clear display window or allow any utensil to puncture the display window or keypad.
 - Do not allow any fluid to enter the device or immerse the 4.
 - instrument in fluids or attempt to sterilize the instrument or any or its accessories.



