

November 20, 2025

**Dear customer,**

PATH MEDICAL is pleased to announce the release\* of new nanoT ymp software (1.0.0.956) and device firmware (build 953).

**NANO TYMP FEATURE UPDATES AND FIXES**

The following new features/fixes will be available for nanoT ymp with the update:

**• Improvements to GDT workflow and connectivity**

- Data is automatically exported if the application is closed. There is no more need to manually trigger the export before closing the application.
- As an alternative to PDF data generation, GDT can be configured to export XML files instead.
- A missing version field in the GDT request file is now tolerated by the nanoT ymp software.

**• Noah integration improvements and fixes**

- "File did not validate" issue: Under certain conditions, reflex measurements which were interrupted during test prevented data storage to Noah. This behavior has been corrected.
- nanoT ymp serial number and last calibration date are now stored within the Noah database.

- **Tympanometry and Reflex testing**

- The nanoTymp probe and device LED blinking pattern has been simplified. It now provides a clear indication of when a measurement is ongoing and when it is finished.
- Numerical values of tympanogram peak compliance and peak pressure are always provided, even if the tympanogram is considered flat according to the built-in criteria.

- **Service and functional test**

- The pump unit test can be launched directly from the settings page without the need to navigate through the Windows start menu (applies to Windows version of the nanoTymp software only)

*\*Please note that regulatory requirements in your territory may result in delays to this update, as notification and clearance may be necessary.*