

April 20, 2023

Dear customer,

PATH MEDICAL is pleased to announce the release of nanoTymp software (1.0.0.804) and device firmware (build 801).

NANO TYMP FEATURE UPDATES AND FIXES

The following new features/fixes are available for nanoTymp

- **Tympanometry**
 - Flat tympanograms are correctly classified as type B (previously empty)
 - Scale range of 1000 Hz tympanometry adjusted
- **Reflex**
 - Recorded reflex traces can be deleted to start over
- **Miscellaneous**
 - Report date when printing from NOAH is shown correctly
 - Pump Unit Test application now included (Windows platform)
 - Ability to enter license key for 1000 Hz update
 - French translation added

For further information about PATH MEDICAL products please refer to <http://www.pathme.de> or contact info@pathme.de.
