

August 21, 2023

Dear customer,

PATH MEDICAL is pleased to announce the release of new Mira PC software (2.4.5).

The update of Mira 2.4.5 is described below. Please pay particular attention to the important information described in chapter 2 ("Important Notices").

1. PC SOFTWARE

- **Mira:**
 - Add Chinese user interface translation.

2. IMPORTANT NOTICES

Please read the following important information carefully and make sure to inform your customers about this information.



PATH MEDICAL STRONGLY RECOMMENDS TO USE THE LATEST MIRA AND THE LATEST PATH SERVICE TOOL INCLUDING THE LATEST SPEECH SAMPLES AND FIRMWARE AND TO UPDATE ALL DEVICES WITH THE LATEST FIRMWARE BY USING MIRA OR THE PATH SERVICE TOOL UPDATE FUNCTION. PLEASE NOTE THAT SOME DEVICE TYPES ARE NOT RECOGNIZED BY OLDER MIRA AND PATH SERVICE TOOL VERSIONS.



WHEN UPDATING DEVICE FIRMWARE, THE TEST MODULE PARAMETERS ON THE DEVICE MAY BE RESET TO FACTORY PRESETS UNDER SPECIAL CIRCUMSTANCES (E.G. WHEN UPDATING DEVICES WITH FIRMWARE 1.6.4485 OR LOWER), I.E. INDIVIDUALLY SET TEST MODULE PARAMETERS MAY BE OVERWRITTEN. THEREFORE, PLEASE MAKE SURE TO KEEP TRACK OF INDIVIDUAL PARAMETERS BEFORE UPDATING THE FIRMWARE.

3. GENERAL INFORMATION

For service requests (i.e., technical or handling problems) regarding Senti and Sentiero hardware, firmware or software and related accessories and service tools please contact PATH MEDICAL service at service@pathme.de.

For further information about PATH MEDICAL products please refer to <http://www.pathme.de> or contact info@pathme.de.
